CONSORT

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

K120

January 1997

Table of contents

Introduction	1
Manufacturer	1
Warranty	1
Servicing	1
Specifications	2
Inputs	2
Error Codes	2
Recorder output	2
Maintenance	3
Calibration	3
Conductivity measurement	3
Temperature measurement	3
Platinizing the cell	3

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This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

Introduction

The normal operating temperature should be between 4° and 40°C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

Connect the instrument only to an earthed power line. The required power source is indicated on the label at the back of the instrument. Do not cut and splice the power cord. When removing the power cord from the wall outlet, be sure to unplug by holding the plug attachment and not by pulling the cord. Do not hold the plug by wet hand.

CONSORT nv
Parklaan 36
B2300 Turnhout
Belgium

Tel (++32)(14)41 12 79
Fax (++32)(14)42 91 79

Manufacturer

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

Warranty

CONSORT is not liable for consequential damages arising out of the use or handling of its products.

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Servicing

Manual
Cables
Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

Specifications

Conductivity	0...100 mS/cm
Temperature	0...100°C
Resolution	0.5°C, 0.1 µS/cm
Inputs	BNC input for conductivity cell BANANA input for a Pt1000 temperature probe BNC input for platinising a conductivity cell
Temp. Comp.	automatic with Pt1000 or manual
Analogue output	0...2 V, ca 8 kΩ
Display	4 digit/character, 13 mm LED display, shows readings, °C & messages
Ambient temp.	4...40°C
Rel. Humidity	0...90 % (non-condensing!)
Power supply	210-250 VAC (* 105-125 V~), 50/60 Hz, max. 6 VA
Cabinet	rugged plastic/metal cabinet
Dimensions	300 x 180 x 120 mm
Weight	2 kg

Each time the instrument is switched on, it will check briefly its internal calibrations while showing all functions of the display and the software-version. An error message appears when something is wrong.

The measuring cell should be connected to the coaxial INPUT-1. Automatic temperature compensation and temperature measurements are possible by plugging a Pt1000 thermocompensator into the °C terminals. You can also use a combination conductivity electrode with built in Pt1000. Its banana plug should be inserted in the upper °C terminal. Without compensator, the manual temperature compensation is automatically switched on.

Inputs

- [Or. EC]** = Overrange Conductivity (e.g. short circuited electrode cable).
[Or. °C] = Overrange °C (measured degrees out of range).
[Err.MEM] = Fatal memory error (all factory settings are cleared. Have your instrument serviced!).

Error codes

A recorder can be connected to the red (+) and the black (-) terminals. The 0 to 2 V output range corresponds with 0 to 1000 digits. Use only laboratory recorders with a high input impedance!

Recorder output

A conductivity cell is active and stable only after wetting! For this purpose it must be immersed for **at least one hour** in distilled water. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

A polluted cell may be cleaned with a soft detergent or diluted nitric acid. Greasy substances may be removed with acetone or alcohol (never do this with plastic electrodes!).

The cell constant depends on the physical dimensions of the platinum plates. They are coated with a spongy black platinum deposit. This increases greatly the effective surface and reduces polarising effects. Since the coating of platinum black is quite resistant to chemical action, the cell constant remains fairly constant during the electrode life time.

1. Select the **EC** mode.
2. After rinsing the cell several times with a 0.01 M KCl standard solution (1413 $\mu\text{S}/\text{cm}$), immerse it in that standard.
3. Put the $^{\circ}\text{C}$ adjuster to the temperature of the standard (not necessary when using a Pt1000 thermocompensator).
4. Turn the cm^{-1} adjuster until the instrument reads the conductivity of the standard solution at 25°C (1.41 mS/cm).

1. Select the **EC** mode.
2. Rinse the electrodes with the sample and immerse them in it.
3. Read the display (remember: compensated 25° !). The autoranging feature gives you the largest possible number of digits.
4. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

- * Stirring the solution during the measurements promotes the homogeneity and is therefore always recommended.
- * When starting a titration in the EC mode, switch to Titr. to lock the actual range and avoid cross-over errors due to differences in the measuring frequency.

Provided a Pt1000 thermocompensator is used, read temperature by selecting the $^{\circ}\text{C}$ mode.

1. Clean the cell thoroughly by immersion in a warm solution of potassium dichromate in concentrated sulfuric acid.
2. Connect the cell to the coaxial input-2 (cathode).
3. Connect another platinum electrode (anode) to the banana connector POL.
4. After rinsing the electrodes with distilled water, immerse them in a platinising solution.
5. After ca. 15 min. the cell is replatinised. Rinse the cell with distilled water.

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