
Consort

**MANUAL
HANDLEIDING
MODE D'EMPLOI
ANLEITUNG**

**C1010
C1020**

March 2013

C15010

Input pH/mV: pH, mV

Input EC: Conductivity, Dissolved oxygen

Input °C: Temperature (Pt1000)

C1020

Input pH/mV: pH, mV

Input EC: Conductivity, Dissolved oxygen

Input °C: Temperature (Pt1000)

Index

General informations.....	1
Installation of the instrument.....	2
pH/mV/temperature measurement.....	2
Conductivity/TDS/Salinity measurement.....	3
Dissolved oxygen measurement.....	3
Maintenance of electrodes.....	4
Good Laboratory Practice.....	5
Data Storage in the internal memory.....	5
Algemene informatie.....	7
Instelling van het toestel.....	8
pH/mV/temperatuur-meting.....	8
Geleidbaarheid/TDS/Saliniteitsmeting.....	9
Opgeloste-zuurstof-meting.....	9
Onderhoud van de elektroden.....	10
Goede Laboratorium Praktijk.....	11
Data-opslag in het intern geheugen.....	11
Informations générales.....	13
Installation de l'appareil.....	14
Mesure de pH/mV/Température.....	14
Mesure de Conductivité/TDS/Salinité.....	15
Mesure d'Oxygène dissous.....	15
Entretien des électrodes.....	16
Bonnes Pratiques de Laboratoire.....	17
Collection de données dans la mémoire interne.....	17
Allgemeine Informationen.....	19
Installation des Gerätes.....	20
pH/mV/Temperatur-Messung.....	20
Leitfähigkeit/TDS/Salinitätsmessung.....	21
Gelöster-Sauerstoff-Messung.....	21
Wartung der Elektroden.....	22
Gute Laborpraxis.....	23
Datenspeicherung im internen Speicher.....	23
Specifications.....	25
NIST pH buffers.....	26
Standard pH buffers.....	26
Conductivity standards.....	26
Redox standards.....	26

General informations

Introduction

This instrument is manufactured with the latest technology and needs no particular maintenance. **Consort** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 0° and 40° C. Never use the instrument in a room with high humidity (>95 %) or at very low temperatures (condensation water!).

Manufacturer

Consort bvba
Parklaan 36
B2300 Turnhout
Belgium

Tel (+32)(14)41 12 79
Fax (+32)(14)42 91 79
E-mail: info@consort.be

Warranty

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **Consort** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **Consort** before returning any product for warranty repair on a freight prepaid basis!

Consort is not liable for consequential damages arising out of the use or handling of its products.

Servicing

In the event of this instrument being returned for servicing, the owner is requested **NOT** to send the following items unless they are suspect:

Manual

Cables

Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **Consort** dealer.

AC adaptor

Connect the mains adapter to the DC socket. Do not hold the adapter by wet hand.

Keyboard

MODE = Selects all modes or escapes from error traps, calibration procedures, etc. by returning to the original mode.

CAL = Starts or proceeds a calibration or a function.

↑↓ = Button for entering a value or for selecting a function.

ON/OFF = Switches the instrument on or off.

Inputs

The measuring electrodes should be connected to the corresponding coaxial pH/mV or EC connectors. Automatic temperature compensation and temperature measurements are possible by plugging a Pt1000 temperature probe into the °C terminals. Without Pt1000, the manual temperature compensation is automatically switched on.

- *Never immerse a pH electrode and a dissolved oxygen electrode together in the same solution!*
- *A blinking decimal point warns you for unstable measurements. Wait to read the display! Stirring the solution during the measurements promotes the homogeneity and is obligatory!*
- *The instrument will refuse automatic calibration when the electrode is unstable. Insufficient stirring or a worn electrode may be the cause.*

Set-up

1. Select [SETTINGS] by pressing **MODE**.
2. Follow the instructions on the screen to adjust language, contrast and automatic power-off timer or to erase the memory. Press **CAL** to confirm or **MODE** to exit.

pH/mV/temperature measurement

pH measurement

1. Select the desired mode (pH) by pressing **MODE**. The display will immediately show the measured value according to the previous calibration. Should you want to recalibrate, press **CAL**.
2. The instrument permits to choose between the buffers in memory (1.68, 2.00, 4.00, 4.01, 6.87, 7.00, 9.18, 9.21, 10.01, 12.00, 12.45). Select the proper values and press **CAL**. The unused buffers should be switched off.
3. Rinse the electrodes with distilled water and immerse them in the first buffer solution. Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
4. After rinsing the electrodes with distilled water, immerse them in the samples and read the display.
5. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.

mV measurement

1. Select the desired mode (mV) by pressing **MODE**. The display will immediately show the measured value according to the previous calibration. Should you want to recalibrate, press **CAL**.
 2. Immerse the electrodes in a standard solution of known potential. Select [CALIBRATE] and press **CAL**. Adjust to the proper value and press **CAL**.
- *Select [RESET] and press CAL to reset the calibration.*

Temperature measurement

1. Select the desired mode (°C) by pressing **MODE**. Should you want to recalibrate, press **CAL**.
 2. Immerse the Pt1000 in a solution of known temperature. Select [CALIBRATE] and press **CAL**. Adjust to the proper value and press **CAL**.
- *Select [RESET] and press CAL to reset the calibration.*
 - *Without Pt1000, press CAL, adjust the manual temperature compensation and proceed by pressing CAL.*

Conductivity/TDS/Salinity measurement

Conductivity measurement

1. Select the desired mode (S/cm) by pressing **MODE**. The display will immediately show the measured value according to the previous calibration. Should you want to recalibrate, press **CAL**.
2. The instrument permits to choose between the standards in memory (0.01, 0.1, 1 M KCl). Select the proper values and press **CAL**. The unused standards should be switched off.
3. After rinsing the electrode several times with the first standard solution, immerse it in that standard. The solution temperature is not so critical but should lie between 0°C and 40°C. When no Pt1000 is used, do not forget to compensate manually first! Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
4. Rinse the electrode several times with the sample, immerse it in that solution and read the display.
5. Rinse the electrode always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

TDS measurement (C1020 only)

1. Select the desired mode (mg/l)(TDS) by pressing **MODE**. The display will immediately show the measured value.

Salinity measurement (C1020 only)

1. Select the desired mode (SAL) by pressing **MODE**. The display will immediately show the measured value.

Dissolved oxygen measurement

Dissolved oxygen measurement

1. Select [O₂] by pressing **MODE**. The display will immediately show the measured value according to the previous calibration. Should you want to recalibrate, press **CAL**.
2. The display shows the air pressure correction and the salinity correction. Select the proper values and press **CAL**. Leave salinity correction to zero unless you are going to measure in heavily salted solutions such as e.g. sea-water (35 g/l). Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
3. The electrode exposed to the air reaches an equilibrium corresponding to the partial pressure of oxygen and thus to saturation in water at the given temperature. The instrument shows the measured saturation, current, temperature, and will calibrate automatically when readings are stable.
4. After rinsing the electrodes with distilled water, immerse them in the samples and read the display. Stirring the solution during the measurements promotes the homogeneity and is **obligatory!** The advection rate must be at least 10 cm/s.
5. Rinse the electrode always after use and store it in distilled water.

Interferences

All substances which can diffuse through the membrane and for which 800 mV potential suffices for polarographic reduction, will be reduced in the electrode. This will give a corresponding current contribution, if they are present. Interference can be caused by ions entering the electrode through porous or mechanically damaged membranes and by diffusion of other reactive gases apart from oxygen, e.g. CO₂, Cl₂, SO₂, and H₂S. These substances react in undesired manner with the electrode. Acidic or basic gases change the pH value of the electrolyte solution and thus disturb the reading, particularly when measuring small oxygen concentrations. High salt concentrations in the sample solution can falsify readings too.

pH electrodes

A pH electrode is active and stable only after wetting! For this purpose it must be immersed for at least ten hours in a 3...4 M KCl solution. During short interruptions (e.g. storage) the electrode should be immersed in a 3...4 M KCl solution. In doing this it is always kept ready for use. When the interruption is longer than a month, refill the closing cap with 3...4 M KCl and plug it on the electrode tip in order to protect the glass bulb. Before use, ensure that the reference part of the electrode is topped up with a 3...4 M KCl solution.

Avoid a low pressure inside the electrode! Therefore always remove the closure from the refilling aperture during the measurements as well as during the calibration. This allows the saltbridge solution to flow through the ceramic liquid junction and prevents contamination of the electrolyte. For the same reason, the inside level should always be higher than the outside level of the measuring solution. Close the refilling aperture again when storing the electrode.

A polluted electrode may be cleaned with a soft detergent or 0.1 M HCl. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

If the electrode is polluted by proteinaceous materials (such as blood), it should stand in a cleaning solution overnight and then be cleaned with distilled water before use. The pH electrode wears away by being used. If the electrode tends to respond slower and calibration becomes difficult, even after cleaning, it should be replaced by a new one.

Metal electrodes

Metal electrodes (Pt, Ag, Au) are always ready for use. During short interruptions they are immersed in distilled water. **They should be cleaned regularly:**

- Silver electrodes are immersed in a concentrated ammonia solution during one hour.
- Platinum or gold electrodes are immersed in concentrated nitric acid during one hour.

Conductivity electrodes

A conductivity cell is active and stable only after wetting! For this purpose it must be immersed for at least one hour in distilled water. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

A polluted cell may be cleaned with a soft detergent or diluted nitric acid. Greasy substances may be removed with acetone or alcohol (never do this with plastic electrodes!).

Dissolved oxygen electrodes

A dissolved oxygen electrode is active and stable only when polarised! The electrode body has almost unlimited life and requires no maintenance. After prolonged use of the electrode, it may become deactivated. An indication is that the electrode no longer responds correctly to calibration. In this case:

1. Carefully screw off the measuring head.
 2. Cautiously remove the precipitated deposits with the aid of filter paper. **Do not use grinding paper or a glass fibre!**
 3. Rinse several times with distilled water and shake carefully off the water drops.
 4. Take the measuring head, fill it with new electrolyte solution and very slowly screw it onto the electrode while holding it vertically.
 5. **Let the electrode rest for a few hours!** The electrode is now ready for use again.
- Replace the membrane only when damaged, **NOT** when calibration is no longer possible!

Good Laboratory Practice

1. Select the desired range by pressing **MODE**.
2. Press **CAL**, select GLP and press **CAL** to display a complete calibration report.

Data Storage in the internal memory

Storing data (C1020 only)

1. Press **↓** to store the displayed value in memory. The display shows a next storage address.
1. Select the desired address and press **CAL** to enter the measured value in memory.

Recalling data (C1020 only)

1. Press **↑** to recall the stored data. The display shows a list of stored values.
2. Select the desired address to read the stored values.
3. Press **MODE** to return to the normal measurements.

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Specifications

Specifications	C10x0	
pH	Range	0...14 pH
	Resolution	0.01 pH
	Accuracy	0.2% ± 1 digit
	Calibration	1...3 points
	Buffers	11 pre-programmed
	Temperature compensation	0...100 °C
	ISO-pH	6...8 pH
Slope		80...120%
mV	Range	±1000 mV
	Resolution	1 mV
	Accuracy	0.2% ± 1 digit
	Calibration	1 point
rH ₂	Range	
	Resolution	
CONDUCTIVITY	Range (cc dependent)	0...100 mS/cm (C1010) 0...1000 mS/cm (C1020)
	Resolution (cc dependent)	0.1 µS/cm (C1010) 0.01 µS/cm (C1020)
	Accuracy	1% f.s. of range
	Calibration	1 point
	Standards	3 pre-programmed
	Cell constant (cc)	1 cm ⁻¹ ±30% (C1010) 0.1/1/10 cm ⁻¹ ±30% (C1020)
	Temperature compensation	0...100 °C
	Reference temperature	20° or 25 °C
	Temperature coefficient	natural waters (EN27888)
SALINITY (C1020 only)	Range	0.0...70.0
	Reference temperature	15 °C
TDS (C1020 only)	Range	0...100 g/l
	Resolution	0.1 mg/l
DISSOLVED OXYGEN	Range	0...20 mg/l (0...200%)
	Resolution	0.01 mg/l (0.1%)
	Accuracy	1% ± 1 digit
	Calibration	1 point
	Temperature compensation	0...10 °C
	Salinity compensation	0...40
	Air pressure compensation	800...1200 hPa
TEMPERATURE	Range	0...100 °C
	Resolution	0.1 °C
	Accuracy	0.5 °C
	Calibration	1 point
INPUTS	pH/mV	BNC, 10 ¹² Ω
	Conductivity/Dissolved oxygen	BNC
	Temperature	2 banana, for Pt1000
STORAGE MEMORY (C1020 only)	Data sets	300
DISPLAY	LCD	128x64 pixels
	White backlight	✓
AMBIENT CONDITIONS	Temperature	0...40 °C
	Humidity	0...95%, non condensing
POWER SUPPLY	Mains	100...240 VAC, 50/60 Hz
	Low voltage	9...15 VDC
DIMENSIONS	WxDxH	13x18x10 cm
WEIGHT	Meter	600 g

NIST pH buffers

Temp.	1.68	4.01	6.87	9.18	12.45
0° C	1.666	4.003	6.984	9.464	13.423
5° C	1.668	3.999	6.951	9.395	13.207
10° C	1.670	3.998	6.923	9.332	13.003
15° C	1.672	3.999	6.900	9.276	12.810
20° C	1.675	4.002	6.881	9.225	12.627
25° C	1.679	4.008	6.865	9.180	12.454
30° C	1.683	4.015	6.853	9.139	12.289
40° C	1.694	4.035	6.838	9.068	11.984
50° C	1.707	4.060	6.833	9.011	11.705
60° C	1.723	4.091	6.836	8.962	11.449
70° C	1.743	4.126	6.845	8.921	
80° C	1.766	4.164	6.859	8.885	

Standard pH buffers

Temp.	2.00	4.00	7.00	9.21	10.01	12.00
5° C	1.85	3.99	7.08	9.45	10.24	12.56
10° C	1.89	3.99	7.06	9.38	10.18	12.35
20° C	1.96	3.99	7.02	9.26	10.06	12.14
25° C	2.00	4.00	7.00	9.21	10.01	12.00
30° C	2.03	4.01	6.99	9.16	9.87	11.89
40° C	2.08	4.03	6.98	9.06	9.89	11.71
50° C	2.12	4.06	6.97	8.99	9.83	11.52
60° C	2.14	4.09	6.98	8.93	9.79	11.33

Conductivity standards

Temp.	0.01 M	0.1 M	1 M	(mS/cm)
0° C	0.776	7.15	65.4	
5° C	0.896	8.22	74.1	
10° C	1.020	9.33	83.2	
15° C	1.147	10.48	92.5	
20° C	1.278	11.67	102.1	
25° C	1.413	12.88	111.8	
30° C	1.552	14.12	121.7	
35° C	1.694	15.38	131.7	
40° C	1.839	16.67	142.1	

Redox standards

Temp.	124	358	(mV)
0° C	156	388	
5° C	148	381	
10° C	141	376	
15° C	135	370	
20° C	130	364	
25° C	124	358	
30° C	122	352	
35° C	116	346	
40° C	99	340	

DECLARATION OF CONFORMITY

We declare under our sole responsibility that the product

**Multi-channel analysers
content of the type numbers**

C1010, C1020

to which this declaration relates is in conformity
with the following standards

LOW VOLTAGE DIRECTIVE 2006/95/EC

EN61010-1

EMC DIRECTIVE 2004/108/EC

EN61326-1

Consort bvba

Parklaan 36, B-2300 Turnhout, Belgium

Tel: (+32) (0)14 41 12 79

Fax: (+32) (0)14 42 91 79

Sales: sales@consort.be

The electro-magnetic susceptibility has been chosen at a level that gains proper operation in residential areas, on business and light industrial premises and on small-scale enterprises, inside as well as outside of buildings. All places of operation are characterised by their connection to the public low voltage power supply system.

WARRANTY CERTIFICATE

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory.

Accessories and breakable items such as electrodes are not warranted unless proven to be defective before shipment.

The original purchase order numbers, Consort invoice numbers and serial numbers of the products must be provided.

CONSORT will repair all defective equipment returned to it during the warranty period without charge (CIF Turnhout prepaid by sender), provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures described in the instruction manual and when not having been subject to accident, alteration, misuse or abuse.

If the products have been used with or have come into contact with fluids, an MSDS (material safety data sheet) must be supplied prior to issuing a return authorisation.

A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

Consort bvba

Parklaan 36 • B-2300 Turnhout • Belgium

Tel: (+32) (0)14 41 12 79

Fax: (+32) (0)14 42 91 79

Sales: sales@consort.be

Support: support@consort.be

Information: info@consort.be