
Consort

MANUAL

HANDLEIDING

MODE D'EMPLOI

D230

June 2011

This instrument is manufactured with the latest technology and needs no particular maintenance. **Consort** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 0° and 40°C. Never use the instrument in a room with high humidity (>95 %) or at very low temperatures (condensation water!).

Connect the instrument only to an earthed power line. The required power source is indicated on the label at the back of the instrument. Do not cut and splice the power cord. When removing the power cord from the wall outlet, be sure to unplug by holding the plug attachment and not by pulling the cord. Do not hold the plug by wet hand.

Manufacturer

Consort bvba

Parklaan 36
B2300 Turnhout
Belgium

Tel (+32)(14)41 12 79
Fax (+32)(14)42 91 79
E-mail: info@Consort.be

Warranty

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **Consort** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **Consort** before returning any product for warranty repair on a freight prepaid basis!

Consort is not liable for consequential damages arising out of the use or handling of its products.

Servicing

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

- Manual
- Cables
- Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **Consort** dealer.

Connections

1. The electrodes should be connected to the coaxial BNC connectors.
2. Automatic temperature compensation and temperature measurements are possible by plugging a Pt1000 thermocompensator into the corresponding banana terminals. Without compensator, the manual temperature compensation is automatically switched on.
3. Connect the first rack to the serial port of the computer using the RS232 cable (DP9) supplied with the system.
4. Connect additional racks with each other and to the first rack using the special cables D298.
5. Install the software on the computer and read the help file for further instructions.

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Specifications		D230
pH	<i>Range</i>	0...14 pH
	<i>Resolution</i>	0.001 pH
	<i>Accuracy</i>	0.1% ± 1 digit
	<i>Calibration</i>	1...2 points
	<i>Buffers</i>	9 pre-programmed 2 user specified
	<i>Temperature Compensation</i>	-30...+130 °C
	<i>ISO-pH</i>	6...8 pH
	<i>Slope</i>	80...120%
mV	<i>Range</i>	±2000 mV
	<i>Resolution</i>	0.1 mV
	<i>Accuracy</i>	0.1% ± 1 digit
	<i>Calibration</i>	1 point
ION	<i>Range</i>	0.01 ng/l...100 g/l
	<i>Resolution</i>	3 digits
	<i>Accuracy</i>	0.5% ± 1 digit
	<i>Calibration</i>	2 points + blank
DISS. OXYGEN	<i>Range</i>	0...60 mg/l (0...600%)
	<i>Resolution</i>	0.01 mg/l (0.1%)
	<i>Accuracy</i>	1% ± 1 digit
	<i>Calibration</i>	1 point
	<i>Temperature Compensation</i>	0...50 °C
CONDUCTIVITY	<i>Range</i>	0...2000 mS/cm
	<i>Resolution</i>	0.001 µS/cm
	<i>Accuracy</i>	0.5% f.s. of range
	<i>Calibration</i>	1 point
	<i>Standards</i>	0.01/0.1/1 M KCl 3 user specified
	<i>Cell Constant</i>	0.01/0.1/1/10 cm ⁻¹ ± 50%
	<i>Temperature Compensation</i>	-30...+130 °C
	<i>Reference Temperature</i>	20° or 25 °C
	<i>Temperature Coefficient</i>	natural waters (EN27888)
	<i>Capacitive Compensation</i>	✓
	TDS	<i>Range</i>
<i>Resolution</i>		0.01 mg/l
TEMPERATURE	<i>Range</i>	-30...+130 °C
	<i>Resolution</i>	0.1 °C
	<i>Accuracy</i>	0.3 °C
	<i>Calibration</i>	±10 °C
CHANNELS	<i>pH/mV/Ion/Oxygen</i>	4...28
	<i>Conductivity</i>	4...28
	<i>Temperature</i>	4...28
INPUTS	<i>pH/mV/Ion/Oxygen</i>	BNC, 10 ¹² Ω
	<i>Conductivity</i>	BNC
	<i>Temperature</i>	2 banana, for Pt1000
CALIBRATION	<i>GLP</i>	✓
DATA-LOGGING	<i>Values/channel</i>	unlimited + °C/date/time
	<i>Interval</i>	4 s ... 24 h
SOFTWARE	<i>Languages</i>	EN, NL, FR
SECURITY	<i>Password Protection</i>	✓
POWER SUPPLY	<i>Mains</i>	210...250 VAC
DIMENSIONS	<i>WxDxH</i>	48x24x13 cm
WEIGHT	<i>Meter</i>	12 kg

WARRANTY CERTIFICATE

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory.

CONSORT will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse.

A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

DECLARATION OF CONFORMITY

We declare under our sole responsibility that the product

**Multi-parameter analyser
content of the type numbers
D230, D291, D292**

to which this declaration relates is in conformity
with the following standards

**EN61010
LOW VOLTAGE DIRECTIVE 73/23/EEG
EN50081-1
EN50082-1
EN60555-2
EMC DIRECTIVE 89/336/EEG**

Turnhout, June 20, 2002

A handwritten signature in black ink, appearing to be 'J. G. G.', is located below the date and above the company name.

on behalf of Consort bvba