CONSORT

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

C231
C531

October 1997

Table of contents

Introduction	1
Manufacturer	1
Warranty	1
Servicing	1
Specifications	2
Keyboard	2
Error Codes	2
AC Adapter	2
Inputs	3
DATA ACQUISITION	3
Logging data	3
Recalling data	3
pH-MEASUREMENT	3
Maintenance	3
pH-measurement	4
mV-MEASUREMENT	5
Maintenance	5
mV-measurement	5
CONDUCTIVITY MEASUREMENT	5
Maintenance	5
Conductivity measurement	5
TROUBLE SHOOTING	6
Simple test	6
Calibrate	6
Reset	6

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This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

Introduction

The normal operating temperature should be between 4° and 40°C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

If it is equipped with rechargeable batteries, the instrument should be connected permanently to the mains in order to keep the batteries in good condition. Do not store the instrument for longer than 1 month without recharging the batteries. The batteries can be replaced (C531 only) by opening the bottom compartment of the cabinet.

CONSORT nv

Parklaan 36
B2300 Turnhout
Belgium

Tel (++32)(14)41 12 79
Fax (++32)(14)42 91 79

Manufacturer

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

Warranty

CONSORT is not liable for consequential damages arising out of the use or handling of its products.

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Servicing

Manual
Cables
Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

Specifications

Ranges	pH 0...14 pH mV ±1000 mV Conductivity 0...100 mS/cm °C 0..100°C
Resolution	0.01 pH, 0.1 µS/cm, 1 mV, 0.5°C
Inputs	BNC input for pH/mV electrodes (separate inputs on C231) BANANA (C531) or BNC (C231) input for conductivity cell BANANA input for Pt1000 temperature probe
Temp. Comp.	automatic with Pt1000 or manual
Storage memory	<i>C531 only</i> : 100 values
Display	4 ½ digit/character, 12 mm LCD display, shows readings, °C & messages
Keys	6 tactile membrane keys
Ambient temp.	4...40°C
Rel. Humidity	0...90 % (non-condensing!)
Power supply	<i>C231</i> : built-in NiMH batteries <i>C531</i> : 4 NiMH batteries, size AA AC adapter for 230 V: model A4070 AC adapter for 115 V: model A4071
Cabinet	IP65 cabinet
Dimensions	195 x 101 x 44 mm
Weight	400 g

Each time the instrument is switched on, it will check briefly its internal calibrations while showing all functions of the display and the software-version. An error message appears when something is wrong.

MODE	= Selects all modes or escapes from error traps, calibration procedures, etc..by returning to the original mode.
CAL	= Starts or proceeds a calibration or a function.
▲/▼	= Button for entering a value or for selecting a function.
ON/OFF	= Switches the instrument on or off.

[Or.]	= Overrange.
[CAL]	= Calibration failed (e.g. defective electrode or used buffers).
[MEM]	= Fatal memory error (all factory settings are cleared. Have your instrument serviced!).

Remove the protective silicone stopper from the REF socket and connect the jack of a suitable AC adapter (model **A4070** for 230 V~ or model **A4071** for 115 V~) to the socket for recharging the NiMH batteries. **Avoid to use an AC adapter without NiMH batteries being inserted in the instrument!** For field work, remove the adapter and reinsert the protective silicone stopper into the socket.

Keyboard

Error codes

AC Adapter

The pH/mV -electrode should be connected to the coaxial INPUT connector. The redox-electrode (C231) should be connected to the coaxial mV-input. If separate electrodes are used, connect the reference electrode to the REF. terminal. The conductivity cell should be connected to the AUX terminals (C531) or EC-input (C231). Automatic temperature compensation and temperature measurements are possible by plugging a Pt1000 thermocompensator into the °C terminals. Without compensator, the manual temperature compensation is automatically switched on.

DATA-ACQUISITION

Up to 100 values can be stored in the non-volatile memory of the instrument.

1. Press ∇ to store the displayed value in memory. The display shows briefly [STO] after which a next storage address (1...100) and a blinking [L] appear.
2. Select the desired address with \blacktriangle or ∇ and press **CAL** to enter the measured value in memory.

Logging data

(C531 only)

1. Press \blacktriangle to recall the stored data. The display shows briefly [rCL] after which the next storage address (1...100) and a blinking [L] appears.
2. Select the desired address with \blacktriangle or ∇ and press **CAL** to read the stored value. A blinking [r] appears in the display.
3. Proceed with point 1) to display the value stored in the next address or press **MODE** to return to the measurements.

Recalling data

(C531 only)

pH-MEASUREMENT

A pH electrode is active and stable only after wetting! For this purpose it must be immersed for **at least ten hours** in a 3...4 M KCl solution. During short interruptions (e.g. storage) the electrode should be immersed in a 3...4 M KCl solution. In doing this it is always kept ready for use. When the interruption is longer than a month, refill the closing cap with 3...4 M KCl and plug it on the electrode tip in order to protect the glass bulb. Before use, ensure that the reference part of the electrode is topped up with a 3...4 M KCl solution.

Avoid a low pressure inside the electrode! Therefore always remove the closure from the refilling aperture during the measurements as well as during the standardisation. This allows the saltbridge solution to flow through the ceramic liquid junction and prevents contamination of the electrolyte. For the same reason, the inside level should always be higher than the outside level of the measuring solution. Close the refilling aperture again when storing the electrode.

A polluted electrode may be cleaned with a soft detergent or 0.1 M HCl. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

If the electrode is polluted by proteinaceous materials (such as blood), it should stand in a pepsin solution overnight and then be cleaned before use. The pH electrode wears away by being used. If the electrode tends to respond slower and calibration becomes difficult, even after cleaning, it should be replaced by a new one.

Maintenance

A pH measurement is always relative which makes it necessary to use standards. Each time the pH meter is switched-on or electrodes are changed, the instrument should be standardised with 2 buffers differing from each other by at least 2 units, in order to compensate for:

Zero point:

Each electrode has its own specific pH value, corresponding to the zero potential point (ISO pH). This specific pH value can differ from its theoretical zero point (e.g. 7 pH) by a few tenths of a pH.

Slope:

A good electrode follows the Nernst equation within 90 to 100 %. If the slope lies outside these limits, the electrode becomes slow and may not be optimal for most applications.

1. Select the pH range by pressing **MODE**. The display will immediately show the measured pH according to the previous standardisation. Should you want to restandardise, press **CAL**.
2. Rinse the electrodes with distilled water and immerse them in the first buffer solution.
3. The display shows one of the 9 buffers in memory, e.g. [4.01] while upper indicator is blinking. Select the proper buffer with ▲ or ▼ and press **CAL**. The instrument shows the measured buffer and will standardise automatically when readings are stable ([c] stops blinking).
4. After rinsing the electrodes with distilled water, immerse them in the second buffer solution.
5. The display shows one of the 9 buffers in memory, e.g. [9.18], while the second indicator is blinking. Select the proper buffer with ▲ or ▼ and press **CAL**. The instrument shows the measured buffer and will standardise automatically when readings are stable ([c] stops blinking).
6. After rinsing the electrodes with distilled water, immerse them in the samples and read the display.
7. Rinse the electrodes always with distilled water after use and store them in a 3.5 M KCl solution.

- * Never immerse the pH electrode and the conductivity cell at the same time in the solutions!
- * See page 25 for the buffers in memory.
- * Stirring the solution during the measurements promotes the homogeneity and is therefore always recommended.
- * Read temperature by pressing **MODE**. Without Pt1000, adjust the manual temperature compensation with ▲ or ▼ and proceed by pressing **MODE**.
- * The instrument will refuse automatic standardisation when the electrode is unstable. Insufficient stirring or a worn electrode may be the cause.

mV-MEASUREMENT

Metal electrodes (Pt, Ag, Au): Metal electrodes are always ready for use. During short interruptions they are immersed in distilled water. **They should be cleaned regularly:**

- * Silver electrodes are immersed in a concentrated ammonia solution during one hour.
- * Platinum or gold electrodes are immersed in concentrated nitric acid during one hour.

1. Select the mV range by pressing **MODE**.
2. After rinsing the electrodes with distilled water, immerse them in the solution to be measured and read the display.
3. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.

Maintenance

**mV
measurement**

CONDUCTIVITY MEASUREMENT

A conductivity cell is active and stable only after wetting! For this purpose it must be immersed for **at least one hour** in distilled water. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

A polluted cell may be cleaned with a soft detergent or diluted nitric acid. Greasy substances may be removed with acetone or alcohol (never do this with plastic electrodes!).

Maintenance

The cell constant depends on the physical dimensions of the platinum plates. They are coated with a spongy black platinum deposit. This increases greatly the effective surface and reduces polarising effects. Since the coating of platinum black is quite resistant to chemical action, the cell constant remains fairly constant during the electrode life time.

**Conductivity
measurement**

1. Select the conductivity-range by pressing **MODE**.
2. After rinsing the cell several times with a 0.01 M KCl standard solution (1413 $\mu\text{S}/\text{cm}$), immerse it in that standard. The solution temperature is not so critical but should lie between 0°C and 30°C. When no Pt1000 is used, do not forget to compensate manually first! Press **CAL**.
3. The display shows the reference temperature [r.20] or [r.25]. Select the desired value with \blacktriangle or \blacktriangledown . Press **CAL** to continue.
4. The display will show the computed cell constant e.g. [1.045c] and will calibrate automatically when readings are stable ([c] stops blinking).
5. Rinse the cell several times with the sample, immerse it in that solution and read the display.
6. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

- * Never immerse the pH electrode and the conductivity cell at the same time in the solutions!
- * Stirring the solution during the measurements promotes the homogeneity and is therefore always recommended.
- * Read temperature by pressing **MODE**. Without Pt1000, adjust the manual temperature compensation with ▲ or ▼ and proceed by pressing **MODE**.

TROUBLE SHOOTING

Try this procedure if you suspect the electronic performance of the instrument.

Simple test

1. Short-circuit the pH/mV input(s). Leave the AUX (C531) or EC (C231) input open.
2. Switch the instrument ON.
3. **C231 only**: Select the pH range by pressing **MODE**. The display should show $[7.00] \pm 1$ pH. Run the Calibrate procedure if erroneous values are read.
4. Select the mV ranges by pressing **MODE**. The display should show $[0] \pm 1$ mV. Run the **Calibrate** procedure if erroneous values are read.
5. Select the conductivity range by pressing **MODE**. The display should show $[0.0] \pm 0.1$ $\mu\text{S/cm}$. Run the **Calibrate** procedure if erroneous values are read.

This procedure permits to calibrate automatically all amplifiers of the instrument.

Calibrate

1. Short-circuit the pH/mV input(s). Leave the AUX (C531) or EC (C231) input open.
2. Switch the instrument ON while holding **CAL** pressed.
3. The display shows [SELF] for a while and returns to the measurements.

This procedure resets the instrument to the original factory settings. All personal memorised data (including access code) will be erased!

Reset

1. Switch the instrument ON while holding **MODE** pressed.
2. The display shows [rSEt] for a while and returns to the measurements.

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