CONSORT

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

C732

August 2002

QUALITYMASTERS

CERTIFICATE OF REGISTRATION

This is to certify that:

Consort NV

Parklaan 36
B-2300 Turnhout

Has been assessed by QualityMasters in respect of its Quality Management System and found to comply with:

NEN-EN-ISO 9001 : 2000

Approval is hereby granted for registration providing the rules and conditions relating to certification are observed at all times.

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Lack of fulfillment of conditions as set forth in the certification regulations may render this certificate invalid. The use of the accreditation mark indicates accreditation in respect to the activities covered by the accreditation number EU9706010F

Q

M

REGISTRATION

CERTIFICATE OF

Authorized signatory

A. LEUWATY
QualityMasters B.V.

This certificate remains the property of QualityMasters B.V.

Instruments manufactured by CONSORT

Electrophoresis power supplies

•

Ion analysers

•

pH meters

•

Conductometers

•

Dissolved oxygen meters

•

Thermometers

•

Industrial controllers

CONSORT nv
Parklaan 36
B2300 Turnhout
Belgium

Find out more at www.consort.be

Tel (+32)(14)41 12 79

Fax (+32)(14)42 91 79

E-mail: info@consort.be

DECLARATION OF CONFORMITY

We declare under our sole responsibility that the product

**Multi-parameter analyser
content of the type numbers**

C732

to which this declaration relates is in conformity
with the following standards

EN61010

LOW VOLTAGE DIRECTIVE 73/23/EEG

EN50081-1

EN50082-1

EN60555-2

EMC DIRECTIVE 89/336/EEG

Turnhout, November 10, 1995

A handwritten signature in black ink, appearing to be 'J. G. G.', is located below the date.

on behalf of **CONSORT nv**

Table of contents

Introduction	5
Manufacturer	5
Warranty	5
Servicing.....	5
Specifications.....	6
Keyboard.....	6
Error codes.....	7
Inputs	7
Recorder output.....	7
Digital output.....	7
Reference	8
Date	8
Time.....	8
Temperature	8
Recorder.....	8
Specific features	9
Alternating measurements	9
pH measurement.....	9
mV measurement	10
rH ₂ measurement.....	10
Conductivity measurement.....	11
µW measurement	11
Data-logging	12
Computer.....	13

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This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 4° and 40° C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

Connect the instrument only to an earthed power line. The required power source is indicated on the label at the back of the instrument. Do not cut and splice the power cord. When removing the power cord from the wall outlet, be sure to unplug by holding the plug attachment and not by pulling the cord. Do not hold the plug by wet hand.

Manufacturer

CONSORT nv

Parklaan 36
B2300 Turnhout
Belgium

Tel (+32)(14)41 12 79
Fax (+32)(14)42 91 79
E-mail: info@consort.be

Warranty

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

CONSORT is not liable for consequential damages arising out of the use or handling of its products.

Servicing

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Manual

Cables

Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

Ranges	pH	0...14 pH
	mV	±2000 mV
	rH₂	0...42 rH ₂
	μW	0...400000 μW
	Conductivity	0...20/200/2000 mS/cm (cc = 0.1/1/10)
	Resistivity	0...2/20/200 MΩ.cm (cc = 0.1/1/10)
	Salinity	0...10 %
	°C	0...100°C
Resolution	0.01 pH, 1 mV, 0.001 μS/cm, 1 Ω.cm, 0.01 ppm, 0.1°C	
Inputs	BNC input for a pH electrode	
	BNC input for a mV electrode	
	BNC input for a 2-pole conductivity cell	
	DIN-8 input for a 4-pole conductivity cell	
	BANANA input for a Pt1000 temperature probe	
Temp. Comp.	automatic with Pt1000 or manual	
Digital output	programmable RS232, 150...4800 b/s, for bi-directional communication with a computer or printer	
Analogue output	0...2 V, ca 8 kΩ, programmable 8 bit D/A	
Display	8+4 digit/character, 13 mm LED display	
Keys	8 tactile membrane keys	
Ambient temp.	4...40°C	
Rel. Humidity	0...90 % (non-condensing!)	
Power supply	210-250 V~ (* 105-125 V~), 50/60 Hz, max. 4 VA	
Cabinet	rugged ABS/metal cabinet	
Dimensions	300 x 180 x 120 mm	
Weight	2.4 kg	

Keyboard

MODE	=	Selects all modes or escapes from error traps, calibration procedures, etc..by returning to the original mode.
CAL	=	Starts or proceeds a calibration or a function.
▲/▼	=	Button for entering a value or for selecting a function.
SET	=	Installs the meter according your application.
HOLD	=	Holds display when measuring.
RES	=	Changes resolution/units of the display.
PRINT	=	Displayed value is printed through the RS232 output.
ON/OFF	=	Switches the instrument on or off.

[Or. pH]	=	Overrange pH.
[Or. mV]	=	Overrange mV.
[Or. μ S]	=	Overrange μ S.
[Or. MS]	=	Overrange mS.
[Or. OhM]	=	Overrange Ω .
[Or. °C]	=	Overrange °C.
[no dAtA]	=	Data-logging memory empty.
[not.CAL]	=	Calibration procedure interrupted.
[Err.SLP]	=	Slope not within 80...120 %.
[Err.pHo]	=	ISO pH not within 6...8 pH.
[Err.cc]	=	Cell constant not within 70...130% of the nominal value.
[Err.CAL]	=	Calibration error.
[Err.COd]	=	Invalid code.
[Err.COM]	=	General communication error.
[Err.MEM]	=	Fatal memory error. Have your instrument serviced!
[ILLEGAL]	=	An unauthorised person has failed to repair the unit. Warranty is elapsed from now on!

Inputs

The pH electrode should be connected to the coaxial pH connector. The mV electrode should be connected to the coaxial mV connector. If separate electrodes are used, connect the reference electrodes to the REF. terminals. Connect a 2-pole conductivity cell to the coaxial EC input or connect a 4-pole conductivity cell to the DIN input. Never connect both types at the same time to the instrument! Automatic temperature compensation and temperature measurements are possible by plugging a Pt1000 thermocompensator into the °C terminals. Without compensator, the manual temperature compensation is automatically switched on.

Recorder output

A recorder can be connected to the red (+) and black (-) REC terminals. Use only laboratory recorders with a high input impedance!

Digital output

A standard RS232 output terminal (DP9) is provided for interfacing the instrument with a printer or computer (using the optional AK2310 cable). Data is sent in the ASCII code at a BAUD rate of 150...4800 bps (8 bit, no parity, 1 stopbit).

Serial port pinout specifications:

- pin 1 : connected to pin 4 and pin 6
- pin 2 : TxD, transmit data
- pin 3 : RxD, receive data
- pin 4 : connected to pin 1 and pin 6
- pin 5 : Gnd, signal ground
- pin 6 : connected to pin 1 and pin 4
- pin 7 : connected to pin 8
- pin 8 : connected to pin 7
- pin 9 : not connected

1. Press **SET** until the display shows [SEt rEF] and press **CAL** to enter, with \blacktriangle or \blacktriangledown , the reference type of the redox electrode [CALOMEL] or [AgCL]. Press **CAL** to continue or **SET** to return.
2. The display shows e.g. [rEF=25°C]. Select the temperature to which all future conductivity measurements will be referred with \blacktriangle or \blacktriangledown and press **CAL** to continue.

Date

1. Press **SET** until the display shows [SEt dAtE] and press **CAL** to continue.
2. The display shows [dAtE=on] or [dAtE=OFF]. Select the desired mode with \blacktriangle or \blacktriangledown and press **CAL** to continue.
3. Program the date with \blacktriangle or \blacktriangledown (respectively year, month and day) by pressing **CAL** each time.

Time

1. Press **SET** until the display shows [SEt tiME] and press **CAL** to continue.
2. The display shows [tiME=on] or [tiME=OFF]. Select the desired mode with \blacktriangle or \blacktriangledown and press **CAL** to continue.
3. Program the time with \blacktriangle or \blacktriangledown (respectively hour, minutes and seconds) by pressing **CAL** each time.

Temperature

1. Press **SET** until the display shows [SEt °C] and press **CAL** to continue.
2. The display will show the measured temperature e.g. [°C=23.7][Aut]. If necessary, calibrate your probe with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to return. Press **RES** to reset the temperature calibration to the original factory calibration.
3. When no Pt1000 temperature probe is connected to the instrument the display will show the manual temperature e.g. [°C=23.7][MAn]. Shift to the proper temperature with \blacktriangle or \blacktriangledown and press **CAL** to continue.

Recorder

1. Press **SET** until the display shows [SEt rEC] and press **CAL** to continue.
2. The display shows for all channels e.g. [rEc=PH], [rEc=MV], [rEc=EC], [rEc=°C]...etc. Select the desired mode with \blacktriangle or \blacktriangledown . Press **CAL** to continue.
3. This step is **skipped** unless [rEc=EC] has been chosen. Press **CAL** and the display shows the conductivity measuring range & units, e.g. [rEc 0.01MS][Unit]. Select the desired range with \blacktriangle or \blacktriangledown and press **CAL** to continue.
4. The display shows e.g. [rEc=6.00][LoW]. Select the desired minimum level corresponding to a 0 V recorder output with \blacktriangle or \blacktriangledown and press **CAL** to continue.
5. The display shows e.g. [rEc=8.00][HiGH]. Select the desired maximum level corresponding to a 2 V recorder output with \blacktriangle or \blacktriangledown and press **CAL** to continue.

1. Press **SET** until the display shows [SEtUP] and press **CAL** to continue.
2. The display lets you choose with \blacktriangle or \blacktriangledown between [bEEP on] or [bEEP OFF]. Decide whether or not the beeper should sound each time the keyboard is touched. Press **CAL** to continue.
3. The display shows e.g. [br=2400]. Pre-set the desired baud-rate (150...4800 b/s) with \blacktriangle or \blacktriangledown and press **CAL** to continue.
4. The display shows e.g. [rS=120]. Select the desired time interval (0...9999 s) between the transmitted data (RS232), with \blacktriangle or \blacktriangledown . Pre-set to zero if no automatic printing is required. Press **CAL** to continue.
5. The display shows e.g. [Id.no=5]. Select the desired identification number (0...999) for your meter with \blacktriangle or \blacktriangledown and press **CAL** to continue. The instrument can be identified by e.g. a computer when a specific number is allocated to it. Pre-set to zero to ignore this number.
6. The display shows [Code on], [Code OFF] or [Code CAL]. Decide, with \blacktriangle or \blacktriangledown , whether all future access to the instrument should be by entering a code (on) or not (off). Press **CAL** to continue. Select [Code CAL] to enter your personal code to avoid undesired access to the instrument. The display shows [E=-----] while [=] is blinking. Enter your secret sequence of 5 keys to proceed. Press **CAL** to continue.

Alternating display

1. Press **SET** until the display shows [ALT] and press **MODE** to start the alternating mode. The display alternates between the selected modes with an interval of 4 s. Pressing **PRINT** in this mode will print all values of the selected channels.
2. Press **MODE** again to stop this mode and continue with the actual channel.

pH measurement

1. Select the pH range by pressing **MODE**. The display will immediately show the measured pH according to the previous calibration. Should you want to recalibrate, press **CAL**.
2. Rinse the electrodes with distilled water and immerse them in the first buffer solution.
3. The display shows one of the 9 buffers in memory (b1) or the manual input (M1) of user specified buffers, e.g. [b1=4.01] while [b1] or [M1] is blinking. Select the proper buffer with \blacktriangle or \blacktriangledown and press **CAL**. If you wish to change the manual buffer, first press **RES** to reset to the proper value with \blacktriangle or \blacktriangledown and then press **CAL**.
4. The instrument shows the measured buffer and will calibrate automatically when readings are stable ([=] stops blinking). If the buffer has a different temperature, compensate the indicated value manually with \blacktriangle or \blacktriangledown (not necessary when using a Pt1000).
5. Go on in the same way with the next buffer (b2) by pressing **CAL** each time or press **MODE** to start immediately the measurements.
6. After rinsing the electrodes with distilled water, immerse them in the samples and read the pH-value from the display.
7. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.

Maintenance:

A **pH electrode is active and stable only after wetting!** For this purpose it must be immersed for **at least ten hours** in a 3...4 M KCl solution. During short interruptions (e.g. storage) the electrode should be immersed in a 3...4 M KCl solution. In doing this it is always kept ready for use. When the interruption is longer than a month, refill the closing cap with 3...4 M KCl and plug it on the electrode tip in order to protect the glass bulb. Before use, ensure that the reference part of the electrode is topped up with a 3...4 M KCl solution. **Avoid a low pressure inside the electrode!** Therefore always remove the closure from the refilling aperture during the measurements as well as during the standardisation. This allows the saltbridge solution to flow through the ceramic liquid junction and prevents contamination of the electrolyte. For the same reason, the inside level should always be higher than the outside level of the measuring solution. Close the refilling aperture again when storing the electrode. A polluted electrode may be cleaned with a soft detergent or 0.1 M HCl. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**). If the electrode is polluted by proteinaceous materials (such as blood), it should stand in a pepsin solution overnight and then be cleaned before use. The pH electrode wears away by being used. If the electrode tends to respond slower and calibration becomes difficult, even after cleaning, it should be replaced by a new one.

mV measurement

1. Select the mV range by pressing **MODE**. The display will immediately show the measured mV according to the previous calibration. Should you want to recalibrate, press **CAL**.
 2. Immerse the electrodes in a standard solution of known potential.
 3. The display shows e.g. [MV=204]. Calibrate to the proper value with \blacktriangle or \blacktriangledown and press **CAL** again or press **RES** to reset the absolute readings.
 4. After rinsing the electrodes with distilled water, immerse them in the solution to be measured and read the display.
 5. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.
- * Press RES to switch between mV readings referred to a standard hydrogen electrode [MV.H] or not [MV].

rH₂ measurement

1. Select the rH₂ [rH] range by pressing **MODE**.
2. After calibrating and rinsing the pH and redox electrodes with distilled water, immerse them in the solution to be measured and read the display.

1. Select the conductivity range by pressing **MODE**. The display will immediately show the measured mV according to the previous calibration. Should you want to recalibrate, press **CAL**.
 2. The display will show the cell constant (e.g. [cc=1.073]). Press **CAL** to start the automatic determination procedure.
 3. The display shows one of the KCl standard solutions in memory [St=0.01M][cc] (1413 $\mu\text{S}/\text{cm}$), [St=0.1M][cc] (12.88 mS/cm), [St=1M][cc] (111.8 mS/cm) or the manual input [St=MAAn]. Select the proper standard with \blacktriangle or \blacktriangledown .
 4. This step is **skipped** unless the manual input has been chosen. The display shows e.g. [MS=10.00][cc] while the units are blinking. Choose with \blacktriangle or \blacktriangledown between 1.000 (units), 10.00 (tens) or 100.0 (hundreds) the desired calibration range and press **CAL** again. The display shows e.g. [MS=1.234][cc] while [=] is blinking. Select the desired standard value (at real temperature!) with \blacktriangle or \blacktriangledown .
 5. After rinsing the cell several times with the standard solution, immerse it the standard solution. The solution temperature is not so critical but it should lie between 0°C and 30°C. If the standard has a different temperature, compensate the indicated value manually with \blacktriangle or \blacktriangledown (not necessary when using a Pt1000). Press **CAL** to continue.
 6. The display will show the computed cell constant e.g. [cc=1.045] and will calibrate automatically when readings are stable ([=] stops blinking). Press **MODE** to start the measurements.
 7. Rinse the electrodes with the sample and immerse them in it.
 8. Read the display (remember: compensated to the reference temperature!). The autoranging feature gives you the largest possible number of digits + their unit.
 9. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).
- * Press RES to switch between ppm (salinity), S/cm (conductivity) or $\Omega\cdot\text{cm}$ (resistivity: E6 = M Ω , E3 = k Ω , E0 = Ω).

Maintenance:

A conductivity cell is active and stable only after wetting! For this purpose it must be immersed for **at least one hour** in distilled water. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition). A polluted cell may be cleaned with a soft detergent or diluted nitric acid. Greasy substances may be removed with acetone or alcohol (never do this with plastic electrodes!).

1. Select the μW range by pressing **MODE**.
2. After calibrating and rinsing the pH, redox and conductivity electrodes with distilled water, immerse them in the solution to be measured and read the power (E-6 = μW) (E-3 = mW) (E-0 = W).

Pre-set and start the data-logging:

1. Press **SET** until the display shows [SEt LoG] and press **CAL** to continue.
2. Select [Log dAtA] with \blacktriangle or \blacktriangledown and press **CAL** to continue.
3. Choose with \blacktriangle or \blacktriangledown between [LOG PH], [LOG MV], [LOG MS],... etc. Select the desired mode and press **CAL**.
4. The display shows e.g. [dt=100]. Select the desired time interval (0...9999 s) between the data-logging, with \blacktriangle or \blacktriangledown and press **CAL**.
5. This step is **skipped** unless manual data-logging is chosen. Pre-set the time-interval (dt) to zero and press **CAL** to continue (point 6 is skipped).
6. The display shows e.g. [n=40]. Select the desired number of values (1...500) to be data-logged with \blacktriangle or \blacktriangledown and press **CAL**.
7. The display shows [StArt ?]. Press **CAL** to start the data-logging according the previous settings (a blinking [L] appears) or press **MODE** to return. After the last value has been data-logged the display will show [End LOG]. When manual data-logging has been selected, press **PRINT** to put a next measurement into memory. In the meantime the display shows the logging-number e.g. [LOG 7]. Press **MODE** to return.

Continue with a next series of measurements, without erasing the previous data-logging:

1. Press **SET** until the display shows [SEt LoG] and press **CAL** to continue.
2. Select [Cnt LoG] with \blacktriangle or \blacktriangledown provided the memory has enough capacity left! Press **CAL** to continue to continue with a next series of measurements, without erasing the previous data-logging.
3. Press **SET** to return.

View the stored values on the display:

1. Press **SET** until the display shows [SEt LoG] and press **CAL** to continue.
2. Select [LSt dAtA] with \blacktriangle or \blacktriangledown and press **CAL**.
3. The display shows the first data-logged value in memory (a blinking [n] appears). Press \blacktriangle or \blacktriangledown to view all other values.
4. While pressing **HOLD** the corresponding date is shown.
5. While pressing **RES** the corresponding time shown.
6. Press **PRINT** to send the recalled data to a printer.
7. Press **SET** to return.

Send the stored values to a printer or computer:

1. Press **SET** until the display shows [SEt LoG] and press **CAL** to continue.
2. Select [Prt dAtA] with \blacktriangle or \blacktriangledown and press **CAL** to print the data.
3. The display shows [PrintinG] while all data-logged values, including corresponding dates and times, are sent to the RS232 output in the following format, compatible with EXCEL™, LOTUS 123™,....The columns are separated by a TAB character:

Identification No: 037

#	VALUE	UNIT	C	CH	H	D
0001	7.34	pH	23.7	2	13:19:57	29/01/95
0002	7.21	pH	23.7	2	13:20:57	29/01/95
0003	7.08	pH	23.8	2	13:21:57	29/01/95
0004	6.71	pH	23.8	2	13:22:57	29/01/95

Interfacing with a computer:

Send a single character to the instrument and it will execute the corresponding command as follows:

- 1 = press **MODE**.
- 2 = press **CAL**.
- 3 = press **▲**.
- 4 = press **▼**.
- 5 = press **SET**.
- 6 = press **HOLD**.
- 7 = press **RES**.
- 8 = press **PRINT**.
- ? = send display to computer.
- + = keyboard on (manual operation remains possible).
- = keyboard off (manual operation disabled).
- V = start the input of a value (only valid in certain routines).

Follow exactly the same measuring or calibration procedures, as described in this manual, to include any desired command in your computer data-acquisition program. This example continuously prints the measured values on the computer screen:

10	ON ERROR GOTO 70	'makes sure no error occurs
20	OPEN "COM1:2400,N,8,2" AS #1	'prepares computer input
30	PRINT #1,"8"	'PRINT button pressed
40	INPUT #1,V\$	'measured values transmitted
50	PRINT V\$	'received values printed
60	GOTO 30	'repeat procedure
70	RUN	

Input of a value:

The direct input of values ("V" command) is only possible when the instrument is running a routine where also a manual input of these values could be made. If not, it will ignore any "V" command from the computer!

Examples of allowed routines are:

- * temperature input (e.g. manual temperature compensation).
- * control parameter input such as Low, High, ...
- * time settings.
- * input of manual buffer values.

Examples of NOT allowed routines are:

- * during measurements.
- * when values have to be selected rather than changed (e.g. buffer values in memory).
- * while an electrode is being calibrated.

A 16-bit value (2's complement) should be transmitted in the following sequence:

1st character	=	"V" (start a value input)
2nd character	=	highest byte of the value in ASCII
3rd character	=	lowest byte of the value in ASCII
4th character	=	checksum of 2nd and 3rd character in ASCII
5th character	=	LF, linefeed (ASCII-10)

When a correct checksum has been received, the instrument will send a confirmation character "!" to the computer. If not, a question mark "?" is sent. Both are eventually preceded by the identification number. This example makes the computer to enter the value "1000" on instrument #7, in the sophisticated mode:

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
Data acquisition software to transfer data to a personal computer:

- For Windows 95 or higher.
- Allows to log, view and print the incoming values in real time or to retrieve data-logged values from the memory of the instrument.
- All channels of multi-channel instruments can be processed simultaneously.
- Values are shown in separate windows, tables or graphs.
- A graph window also allows to set alarm limits and to view several statistic values.

To receive your free copy of this Data Information System, make sure to fill in all lines, sign and send or fax this page to Consort nv, Parklaan 36, B-2300 Turnhout, Belgium.

Fax.: +32/14/429179

Please PRINT in capitals!

1. purchase date:

2. purchased from:

2. indicate model:

4. serial no:

5. organisation:

6. name:

7. division:

8. address:

9. postal code:

10. city:

11. telephone:

12. telefax:

13. e-mail:

14. I understand that I may use this software only for evaluation and that any other use requires payment of the registration fee.

Signature:

We would like to know your opinion, provide us with your remarks and fax a copy of this page to Consort: +32/14/429179.

1. purchase date:

2. purchased from:

2. indicate model:

4. serial no:

5. organisation:

6. name:

7. division:

8. address:

9. postal code:

10. city:

11. telephone:

12. telefax:

13. e-mail:

- Are you generally satisfied with the performance of the instrument ?
- Are you satisfied with the service of the company you purchased from ?
- What is the degree of satisfaction with:
 - the design of the instrument ?
 - the operating instructions ?
 - the supplied accessories ?
 - the delivery ?
 - the price/quality ratio ?
- What are your suggestions:

<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....