# Consort

---

MANUAL

HANDLEIDING

MODE D'EMPLOI

**D130**

---

October 2006

**Table of contents**

Introduction .....	3
Manufacturer .....	3
Warranty .....	3
Servicing.....	3
Keyboard .....	3
Inputs .....	3
Connections.....	3
Batteries .....	4
Language .....	4
Channels.....	4
Date and Time.....	4
Password .....	4
Data-logging .....	4
Measurement .....	5
Information .....	5
Maintenance of pH electrodes.....	6
Maintenance of metal electrodes .....	6
Maintenance of oxygen electrodes.....	6
Maintenance of ion selective electrodes.....	6

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This instrument is manufactured with the latest technology and needs no particular maintenance. **Consort** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 0° and 40°C. Never use the instrument in a room with high humidity (>95 %) or at very low temperatures (condensation water!).

---

**Manufacturer****Consort nv**

Parklaan 36  
B2300 Turnhout  
Belgium

Tel (+32)(14)41 12 79  
Fax (+32)(14)42 91 79  
E-mail: info@Consort.be

---

**Warranty**

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **Consort** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **Consort** before returning any product for warranty repair on a freight prepaid basis!

**Consort** is not liable for consequential damages arising out of the use or handling of its products.

---

**Servicing**

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Manual  
Cables  
Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **Consort** dealer.

---

**Keyboard**

**MODE** = Selects all modes or escapes from error traps, calibration procedures, etc..by returning to the original mode.  
**↑↓** = Button for adjusting a value or for selecting a function.  
**CAL** = Starts or proceeds a calibration or a function.  
**ON/OFF** = Switches the instrument on or off.

---

**Inputs**

The measuring electrodes should be connected to the BNC input connectors (channel-1 to -8 from left to right). Automatic temperature compensation and temperature measurements are possible by plugging a Pt1000 temperature probe into each of the corresponding black banana terminals. You can also use a combination pH electrode with built-in Pt1000. Its banana plugs should be inserted in the black banana terminals. Without Pt1000, the manual temperature compensation is automatically switched on.

---

**Connections**

1. Put the CD with software in the computer and open "ReadMe.txt" for further instructions to install the system.
2. Connect the instrument to the computer using the USB cable and switch it on.
3. Proceed with the set-up as indicated in the help file.

## THE NEXT SECTIONS SHOULD ONLY BE USED WHEN THE INSTRUMENT IS OPERATED INDEPENDENTLY AND EQUIPPED WITH BATTERIES !

---

### *Batteries*

The batteries can be replaced by opening the cabinet. Lift the two covers, left and right, next to the frontpanel to reach the screws. Only use NiMH batteries type C! The instrument will automatically switch off when no button has been pressed during the last 15 minutes, except when it is calibrating or data-logging.

---

### *Language*

1. Press **MODE**.
2. Select [LANGUAGE] and press **CAL**.
3. Select the desired language and press **CAL**.

---

### *Channels*

1. Press **MODE**.
2. Select [CHANNELS] and press **CAL**.
3. Select the desired channel and press **CAL** to enable (✓) or disable (-) it.

---

### *Date and Time*

1. Press **MODE**.
2. Select [CLOCK] and press **CAL**.
3. Select [date] and press **CAL** to adjust the date.
4. Select [time] and press **CAL** to adjust the time.

---

### *Password*

1. Press **MODE**.
2. Select [PASSWORD] and press **CAL**.
3. A private code can be programmed to avoid undesired programming or calibration of the instrument. Press **CAL** and follow the instructions on the screen to enter your secret sequence of 5 keys.

---

### *Data-logging*

1. Press **MODE**.
2. Select [DATA-LOGGER] and press **CAL**.

#### **Set-up:**

1. Select [PARAMETERS] and press **CAL**.
2. Select [INTERVAL] and press **CAL**.
3. Select the desired time interval between the data-logging and press **CAL**.
4. Select [STOP] and press **CAL**.
5. Select [Continuous] (erases automatically the first values when the memory is full to free memory for the next ones) or [Memory limit] (data-logging stops when the memory is full) and press **CAL**.

#### **Start/Stop data-logger:**

1. Select [START] and press **CAL**. When the data-logging is still active, press [STOP] to stop it.

#### **Erase table:**

1. Select [ERASE TABLE] and press **CAL**. Follow the instructions on the screen to erase the stored values.

#### **Information:**

1. Select [INFO] to see whether the data-logger is active, the number of recorded values, when the next values will be recorded...

- *The instrument cannot be switched off when the data-logger is active!*

1. Press **MODE** until mV values are shown. The display will always show the raw signals of the electrodes in mV. Select the desired channel with the arrows.
2. Should you want to recalibrate, press **CAL**.
3. The display shows all the standards in memory:  
pH: 1.68, 4.00, 4.01, 6.87, 7.00, 9.18, 9.21, 10.00, 12.45, Manual.  
Ion: 1, 10, 100 ng/l, 1, 10, 100 µg/l, 1, 10, 100 mg/l, 1 g/l.  
O<sub>2</sub>: 100% saturation.  
mV: adjust potential.  
°C: adjust temperature.  
Select the proper value suitable for to the connected electrode and press **CAL**.
4. **pH:**  
rinse the electrodes with distilled water and immerse them in the first standard solution. Press **CAL** and wait until the measurement is stable (mV stops blinking). Press **CAL** to confirm. Repeat eventually the same steps with a second standard solution.  
**mV or °C:**  
adjust the value with the arrows and press **CAL**.  
**Ion:**  
rinse the electrodes with distilled water and immerse them in the first standard solution. Press **CAL** and wait until the measurement is stable (mV stops blinking). Press **CAL** to confirm. Repeat the same steps with a second standard solution. Eventually perform a blank correction after the second calibration.  
**O<sub>2</sub>:**  
hold the electrode in the air and wait until the measurement is stable (mV stops blinking). Press **CAL** to confirm.
5. After rinsing the electrodes with distilled water, immerse them in the samples.
6. Rinse the electrodes always with distilled water after use and store them according to the maintenance instructions.

1. Press **MODE**.
2. Select **[INFO]** to see general information about the instrument.

**A pH electrode is active and stable only after wetting!** For this purpose it must be immersed for **at least ten hours** in a 3...4 M KCl solution. During short interruptions (e.g. storage) the electrode should be immersed in a 3...4 M KCl solution. In doing this it is always kept ready for use. When the interruption is longer than a month, refill the closing cap with 3...4 M KCl and plug it on the electrode tip in order to protect the glass bulb. Before use, ensure that the reference part of the electrode is topped up with a 3...4 M KCl solution.

**Avoid a low pressure inside the electrode!** Therefore always remove the closure from the refilling aperture during the measurements as well as during the calibration. This allows the saltbridge solution to flow through the ceramic liquid junction and prevents contamination of the electrolyte. For the same reason, the inside level should always be higher than the outside level of the measuring solution. Close the refilling aperture again when storing the electrode.

A polluted electrode may be cleaned with a soft detergent or 0.1 M HCl. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

If the electrode is polluted by proteinaceous materials (such as blood), it should stand in a cleaning solution overnight and then be cleaned with distilled water before use. The pH electrode wears away by being used. If the electrode tends to respond slower and calibration becomes difficult, even after cleaning, it should be replaced by a new one.

---

*Maintenance of metal electrodes*

**Metal electrodes (Pt, Ag, Au):** Metal electrodes are always ready for use. During short interruptions they are immersed in distilled water. **They should be cleaned regularly:**

- Silver electrodes are immersed in a concentrated ammonia solution during one hour.
- Platinum or gold electrodes are immersed in concentrated nitric acid during one hour.

---

*Maintenance of oxygen electrodes*

**A dissolved oxygen electrode is active and stable only when polarised!** The electrode body has almost unlimited life and requires no maintenance. After prolonged use of the electrode, it may become deactivated. An indication is that the electrode no longer responds correctly to calibration. In this case:

1. Carefully screw off the measuring head.
  2. Cautiously remove the precipitated deposits with the aid of filter paper. **Do not use grinding paper or a glass fibre!**
  3. Rinse several times with distilled water and shake carefully off the water drops.
  4. Take the measuring head, fill it with new electrolyte solution and very slowly screw it onto the electrode while holding it vertically.
  5. **Let the electrode rest for a few hours!** The electrode is now ready for use again.
- *Replace the membrane only when damaged, not when calibration is no longer possible!*

---

*Maintenance of ion selective electrodes*

We refer to the manual supplied with the ion specific electrodes.

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Specifications		D130
pH	<i>Range</i>	0...14 pH
	<i>Resolution</i>	0.001 pH
	<i>Accuracy</i>	0.1% ± 1 digit
	<i>Calibration</i>	1...2 points
	<i>Buffers</i>	9 pre-programmed 2 user specified
	<i>Temperature Compensation</i>	-30...+130 °C
	<i>ISO-pH</i>	6...8 pH
	<i>Slope</i>	80...120%
mV	<i>Range</i>	±2000 mV
	<i>Resolution</i>	0.1 mV
	<i>Accuracy</i>	0.1% ± 1 digit
	<i>Calibration</i>	1 point
ION	<i>Range</i>	0.01 ng/l...100 g/l
	<i>Resolution</i>	3 digits
	<i>Accuracy</i>	0.5% ± 1 digit
	<i>Calibration</i>	2 points + blank
DISS. OXYGEN	<i>Range</i>	0...60 mg/l (0...600%)
	<i>Resolution</i>	0.01 mg/l (0.1%)
	<i>Accuracy</i>	1% ± 1 digit
	<i>Calibration</i>	1 point
	<i>Temperature Compensation</i>	0...50 °C
TEMPERATURE	<i>Range</i>	-30...+130 °C
	<i>Resolution</i>	0.1 °C
	<i>Accuracy</i>	0.3 °C
	<i>Calibration</i>	±10 °C
CHANNELS	<i>pH/mV/Ion/Oxygen</i>	8
	<i>Temperature</i>	8
INPUTS	<i>pH/mV/Ion/Oxygen</i>	BNC, 10 <sup>12</sup> Ω
	<i>Temperature</i>	2 banana, for Pt1000
CALIBRATION	<i>GLP</i>	✓
DISPLAY	<i>Set-up/Calibration</i>	122x32 pixels
	<i>Languages</i>	EN, NL, FR
DATA-LOGGING	<i>Values (computer)</i>	unlimited + °C/date/time
	<i>Values (internal)</i>	32000 + °C/date/time
	<i>Interval</i>	4 s ... 24 h
SECURITY	<i>Password Protection</i>	✓
POWER SUPPLY	<i>Computer</i>	via USB connection
	<i>Batteries</i>	4 x1.2 V, C, NiMH
DIMENSIONS	<i>WxDxH</i>	24x9x16 cm
WEIGHT	<i>Meter</i>	1.2 kg

## WARRANTY CERTIFICATE

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory.

**CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse.

A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

## DECLARATION OF CONFORMITY

We declare under our sole responsibility that the product

**Multi-parameter analyser**  
**content of the type numbers**

**D130**

to which this declaration relates is in conformity

with the following standards

**EN61010**

**LOW VOLTAGE DIRECTIVE 73/23/EEG**

**EN50081-1**

**EN50082-1**

**EN60555-2**

**EMC DIRECTIVE 89/336/EEG**

Turnhout, October 14, 2004

A handwritten signature in black ink, appearing to be 'J. G. G.', is located below the date of the declaration.

on behalf of **CONSORT** nv