CONSORT

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

K410

January 1997

Table of contents

Introduction	1
Manufacturer	1
Warranty	1
Servicing	1
Error Codes	2
Maintenance	2
Conductivity measurement	2

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Introduction

This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 4° and 40°C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

The batteries can be replaced by opening the bottom compartment of the cabinet. Avoid to use rechargeable batteries!

CONSORT nv
Parklaan 36
B2300 Turnhout
Belgium

Tel (++32)(14)41 12 79
Fax (++32)(14)42 91 79

Manufacturer

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

Warranty

CONSORT is not liable for consequential damages arising out of the use or handling of its products.

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Servicing

Manual
Cables
Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

- [Or] = Overrange (e.g. input short circuited).
[Err.cc] = Cell constant not within 0.7...1.3 cm⁻¹.
[CAL] = Calibration error (check electrodes and standard solutions).
[Err.MEM] = Fatal memory error (all factory settings are cleared. Have your instrument serviced!).

Error codes

A conductivity cell is active and stable only after wetting! For this purpose it must be immersed for **at least one hour** in distilled water. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

Maintenance

A polluted cell may be cleaned with a soft detergent or diluted nitric acid. Greasy substances may be removed with acetone or alcohol (never do this with plastic electrodes!).

The cell constant depends on the physical dimensions of the platinum plates. They are coated with a spongy black platinum deposit. This increases greatly the effective surface and reduces polarising effects. Since the coating of platinum black is quite resistant to chemical action, the cell constant remains fairly constant during the electrode life time.

Conductivity measurement

1. Adjust the manual temperature compensation by pressing **MODE**, e.g. [21.5C]. Then press **CAL**, select the desired temperature with **MODE**, and press **CAL** again.
2. Select the conductivity range by pressing **MODE** and press **CAL**.
3. After rinsing the cell several times with a 0.01 M KCl standard solution (1413 µS/cm), immerse it in that standard.
4. The display will show the computed cell constant e.g. [1.04c] and will calibrate automatically when readings are stable ([c] stops blinking).
5. Rinse the cell with the sample and immerse it in that solution.
6. Read the display (remember: compensated to 25°C!).

- * Stirring the solution during the measurements promotes the homogeneity and is therefore always recommended.
- * Adjust the manual temperature compensation by pressing **MODE**. Then press **CAL**, select the desired temperature with **MODE**, and press **CAL** again.
- * The instrument will refuse automatic calibration when the cell is unstable. Insufficient stirring or a worn cell may be the cause.

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