CONSORT

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

K511

October 1997

Table of contents

Introduction	1
Manufacturer	1
Warranty	1
Servicing	1
Specifications	2
Keyboard	2
Error Codes	2
AC Adapter	3
Inputs	3
Recorder output	3
DATA AQUISITION	3
Logging data	3
Recalling data	3
CONDUCTIVITY MEASUREMENT	3
Maintenance	3
Conductivity measurement	4
TROUBLE SHOOTING	4
Simple test	4
Calibrate	4
Reset	4

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This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

Introduction

The normal operating temperature should be between 4° and 40°C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

If it is equipped with rechargeable batteries, the instrument should be connected permanently to the mains in order to keep the batteries in good condition. Do not store the instrument for longer than 1 month without recharging the batteries. The batteries can be replaced by opening the bottom compartment of the cabinet.

CONSORT nv
Parklaan 36
B2300 Turnhout
Belgium

Tel (++32)(14)41 12 79
Fax (++32)(14)42 91 79

Manufacturer

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

Warranty

CONSORT is not liable for consequential damages arising out of the use or handling of its products.

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Servicing

Manual
Cables
Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

Specifications

Conductivity	0...100 mS/cm
Temperature	0...100°C
Resolution	0.5°C, 0.1 µS/cm
Inputs	BNC input for conductivity cell BANANA input for a Pt1000 temperature probe
Temp. Comp.	automatic with Pt1000 or manual
Storage memory	200 values
Analogue output	raw electrode signal
Display	4 ½ digit/character, 12 mm LCD display, shows readings, °C & messages
Keys	6 tactile membrane keys
Ambient temp.	4...40°C
Rel. Humidity	0...90 % (non-condensing!)
Power supply	4 NiMH batteries, size AA AC adapter for 230 V: model A4070 AC adapter for 115 V: model A4071
Cabinet	IP65 cabinet
Dimensions	195 x 101 x 44 mm
Weight	400 g

Each time the instrument is switched on, it will check briefly its internal calibrations while showing all functions of the display and the software-version. An error message appears when something is wrong.

MODE	= Selects all modes or escapes from error traps, calibration procedures, etc..by returning to the original mode.
CAL	= Starts or proceeds a calibration or a function.
▲/▼	= Button for entering a value or for selecting a function.
ON/OFF	= Switches the instrument on or off.

Keyboard

[Or]	= Overrange (e.g. input short circuited).
[Err.cc]	= Cell constant not within 0.7...1.3 cm ⁻¹ .
[CAL]	= Calibration error (check electrodes and standard solutions).
[Err.MEM]	= Fatal memory error (all factory settings are cleared. Have your instrument serviced!).

Error codes

Remove the protective silicone stopper from the REF socket and connect the jack of a suitable AC adapter (model **A4070** for 230 V~ or model **A4071** for 115 V~) to the socket for recharging the NiMH batteries. **Avoid to use an AC adapter without NiMH batteries being inserted in the instrument!** For field work, remove the adapter and reinsert the protective silicone stopper into the socket.

AC Adapter

The measuring cell should be connected to the coaxial INPUT connector. Automatic temperature compensation and temperature measurements are possible by plugging a Pt1000 thermocompensator into the °C terminals. You can also use a combination conductivity electrode with built in Pt1000. Its banana plug should be inserted in the upper °C terminal. Without compensator, the manual temperature compensation is automatically switched on.

Inputs

A recorder can be connected to the AUX terminals. The 20 mV output range corresponds with ca. 1000 digits. Use only laboratory recorders with a high input impedance!

Recorder output

DATA ACQUISITION

Up to 200 values can be stored in the non-volatile memory of the instrument.

1. Press ∇ to store the displayed value in memory. The display shows briefly [STO] after which a next storage address (1...200) and a blinking [L] appear.
2. Select the desired address with \blacktriangle or ∇ and press **CAL** to enter the measured value in memory.

Logging data

1. Press \blacktriangle to recall the stored data. The display shows briefly [rCL] after which the next storage address (1...200) and a blinking [L] appears.
2. Select the desired address with \blacktriangle or ∇ and press **CAL** to read the stored value. A blinking [r] appears in the display.
3. Proceed with point 1) to display the value stored in the next address or press **MODE** to return to the measurements.

Recalling data

CONDUCTIVITY MEASUREMENT

A conductivity cell is active and stable only after wetting! For this purpose it must be immersed for **at least one hour** in distilled water. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

Maintenance

A polluted cell may be cleaned with a soft detergent or diluted nitric acid. Greasy substances may be removed with acetone or alcohol (never do this with plastic electrodes!).

Conductivity measurement

The cell constant depends on the physical dimensions of the platinum plates. They are coated with a spongy black platinum deposit. This increases greatly the effective surface and reduces polarising effects. Since the coating of platinum black is quite resistant to chemical action, the cell constant remains fairly constant during the electrode life time.

1. Select the conductivity-range by pressing **MODE**.
2. After rinsing the cell several times with a 0.01 M KCl standard solution (1413 $\mu\text{S/cm}$), immerse it in that standard. The solution temperature is not so critical but should lie between 0°C and 30°C. When no Pt1000 is used, do not forget to compensate manually first!
3. Press **CAL**. The display shows the reference temperature [r.20] or [r.25]. Select the desired value with \blacktriangle or \blacktriangledown . Press **CAL** to continue.
4. The display will show the computed cell constant e.g. [1.045c] and will calibrate automatically when readings are stable ([c] stops blinking).
5. Rinse the cell several times with the sample, immerse it in that solution and read the display.
6. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

* Stirring the solution during the measurements promotes the homogeneity and is therefore always recommended.

* Read temperature by pressing **MODE**. Without Pt1000, adjust the manual temperature compensation with \blacktriangle or \blacktriangledown and proceed by pressing **MODE**.

TROUBLE SHOOTING

Try this procedure if you suspect the electronic performance of the instrument.

1. Leave the conductivity input open.
2. Switch the instrument ON.
3. Select the conductivity range by pressing **MODE**. The display should show [0.0] \pm 0.1 $\mu\text{S/cm}$. Run the **Calibrate** procedure if erroneous values are read.

This procedure permits to calibrate automatically all amplifiers of the instrument.

1. Leave the conductivity input open.
2. Switch the instrument ON while holding **CAL** pressed.
3. The display shows [SELF CAL] for a while and returns to the measurements.

This procedure resets the instrument to the original factory settings.

1. Switch the instrument ON while holding **MODE** pressed.
2. The display shows [rESET] for a while and returns to the measurements.

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