CONSORT

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

K610

January 1997

Table of contents

Introduction	1
Manufacturer	1
Warranty	1
Servicing	1
Specifications	2
Keyboard	2
Inputs	2
Error Codes	3
AC Adapter	3
CONDUCTIVITY MEASUREMENT	3
Maintenance	3
Calibration	3
Conductivity measurement	4
TEMPERATURE MEASUREMENT	4
Measurement	4
DATA COMMUNICATION	4
Recorder output	4
Digital output	4
Printer connection	5
TROUBLE SHOOTING	5
Simple test	5
Calibrate	5
Reset	5

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This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

Introduction

The normal operating temperature should be between 4° and 40°C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

If it is equipped with rechargeable batteries, the instrument should be connected permanently to the mains in order to keep the batteries in good condition. Do not store the instrument for longer than 1 month without recharging the batteries. The batteries can be replaced by opening the bottom compartment of the cabinet.

CONSORT nv
Parklaan 36
B2300 Turnhout
Belgium

Tel (++32)(14)41 12 79
Fax (++32)(14)42 91 79

Manufacturer

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

Warranty

CONSORT is not liable for consequential damages arising out of the use or handling of its products.

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Servicing

Manual
Cables
Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

Specifications

Conductivity	0...200 mS/cm
Resistivity	0...20 MΩ.cm
Temperature	0...100°C
Resolution	0.1°C, 0.1 μS/cm, 10 Ω.cm
Inputs	BNC input for conductivity cell BANANA input for a Pt1000 temperature probe
Temp. Comp.	automatic with Pt1000 or manual
Digital output	programmable RS232, 2400 b/s
Analogue output	0...2 V, ca 8 kΩ, 8 bit D/A
Display	8 digit/character, 12 mm LCD display, shows readings, °C & messages
Keys	10 tactile membrane keys
Ambient temp.	4...40°C
Rel. Humidity	0...90 % (non-condensing!)
Power supply	210-250 V~ (* 105-125 V~), 50/60 Hz, max. 2 VA, & 4 NiCd batteries, size AA
Cabinet	IP65 cabinet
Dimensions	252 x 121 x 50 mm
Weight	600 g

Each time the instrument is switched on, it will check briefly its internal calibrations while showing all functions of the display and the software-version. An error message appears when something is wrong.

MODE	= Selects all modes or escapes from error traps, calibration procedures, etc..by returning to the original mode.
CAL	= Starts or proceeds a calibration or a function.
↖/↗	= Button for entering a value or for selecting a function.
°C	= Reads temperature, when a Pt1000 is present, or manual temperature compensation.
HOLD	= Holds display when measuring.
RES	= Changes resolution/units of the display.
PRINT	= Displayed value is printed through the RS232 output.
ON/OFF	= Switches the instrument on or off.

Keyboard

The measuring cell should be connected to the coaxial INPUT connector. Automatic temperature compensation and temperature measurements are possible by plugging a Pt1000 thermocompensator into the °C terminals. You can also use a combination conductivity electrode with built in Pt1000. Its banana plug should be inserted in the upper °C terminal. Without compensator, the manual temperature compensation is automatically switched on.

Inputs

- [Or. EC] = Overrange Conductivity (e.g. short circuited electrode cable).
- [Or. OhM] = Overrange Resistivity (e.g. interrupted electrode cable).
- [Or. °C] = Overrange °C (measured degrees out of range).
- [not.CAL] = Calibration procedure interrupted (the previous calibration remains valid!).
- [Err.cc] = Cell constant not within 70...130% of the nominal value.
- [Err.CAL] = Calibration error (check electrodes and standard solutions).
- [Err.Stb] = Electrode not stable.
- [Err.Prt] = Communication error with the attached printer.
- [Err.bAt] = Battery charger defective (disconnect the AC-adapter immediately from the mains, remove the batteries and have your instrument serviced!).
- [Err.MEM] = Fatal memory error (all factory settings are cleared. Have your instrument serviced!).
- [ILLEGAL] = An unauthorised person has failed to repair the unit.
Warranty is elapsed from now on!

Error codes

Remove the protective silicone stopper from the REF socket and connect the jack of a suitable AC adapter (model **A4070** for 230 V~ or model **A4071** for 115 V~) to the socket for recharging the NiCd batteries. **Avoid to use an AC adapter without NiCd batteries being inserted in the instrument!** For field work, remove the adapter and reinsert the protective silicone stopper into the socket.

AC Adapter

CONDUCTIVITY MEASUREMENT

A conductivity cell is active and stable only after wetting! For this purpose it must be immersed for **at least one hour** in distilled water. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

Maintenance

A polluted cell may be cleaned with a soft detergent or diluted nitric acid. Greasy substances may be removed with acetone or alcohol (never do this with plastic electrodes!).

The cell constant depends on the physical dimensions of the platinum plates. They are coated with a spongy black platinum deposit. This increases greatly the effective surface and reduces polarising effects. Since the coating of platinum black is quite resistant to chemical action, the cell constant remains fairly constant during the electrode life time.

Calibration

1. Select the conductivity range by pressing **MODE** (the display shows [μS] or [MS]) and press **CAL**.
2. The display will show the cell constant (e.g. [cc=1.073]). Press **CAL** to start the automatic determination procedure.
3. After rinsing the cell several times with the standard solution, immerse it the standard solution. The solution temperature is not so critical but it should lie between 0°C and 30°C. If the standard has a different temperature press **°C** here (not necessary when using a Pt1000). Press **CAL** to continue.
4. The display will show the computed cell constant e.g. [cc1=1.045] and will calibrate automatically when readings are stable ([=] stops blinking).

1. Select the conductivity range by pressing **MODE** (the display shows [μ S] or [MS]).
2. Rinse the electrodes with the sample and immerse them in it.
3. Read the display (remember: compensated to the reference temperature!). The autoranging feature gives you the largest possible number of digits + their unit.
4. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

Conductivity measurement

- * Disconnect the AC-adapter to increase accuracy in the lower conductivity range!
- * See page 21 for the standards in memory.
- * A blinking decimal point warns you for unstable measurements. Wait to read the display!
- * Press **RES** to switch between ppm (salinity), S/cm (conductivity) or Ω .cm (resistivity).
- * Stirring the solution during the measurements promotes the homogeneity and is therefore always recommended.

TEMPERATURE MEASUREMENT

1. Read temperature by pressing $^{\circ}\text{C}$.
2. Without Pt1000, adjust the manual temperature compensation with \blacktriangle or \blacktriangledown and proceed by pressing **MODE**.

Measurement

DATA COMMUNICATION

A recorder can be connected to the red (+) and the black (-) terminals. The 0 to 2 V output range corresponds with 0 to 2000 digits. Use only laboratory recorders with a high input impedance!

Recorder output

A standard RS232 output terminal (DP9) is provided for interfacing the instrument with a printer or computer. Data is sent in the ASCII code at a BAUD rate of 2400 b/s (8 bit, no parity, 1 start and 2 stopbits, automatic line feed ON, data transfer protocol Xon/Xoff).

Digital output

Serial port pinout specifications:

- pin 1 : RLSD, received line signal detector
- pin 2 : TxD, transmit data
- pin 3 : RxD, receive data
- pin 4 : not connected
- pin 5 : Gnd, signal ground
- pin 6 : DSR, data set ready
- pin 7 : not connected
- pin 8 : CTS, clear to send
- pin 9 : not connected

Interfacing with a printer:

Printer connection

Press **PRINT** each time you want to print the measured values. The columns are separated by a TAB character, e.g.:

Log#	VALUE	UNIT	C	H	D
001	16.49E-03	S/cm	23.7	13:19:57	29/01/95
002	15.05E-03	S/cm	23.7	13:20:57	29/01/95
003	13.39E-03	S/cm	23.8	13:21:57	29/01/95
004	1.92E-03	S/cm [^]	23.8	13:22:57	29/01/95
005	0.44E-03	S/cm [^]	23.9	13:23:57	29/01/95

TROUBLE SHOOTING

Try this procedure if you suspect the electronic performance of the instrument.

Simple test

1. Leave the conductivity input open.
2. Switch the instrument ON.
3. Select the conductivity range by pressing **MODE**. The display should show $[0.0\mu\text{S}] \pm 0.1 \mu\text{S/cm}$. Run the **Calibrate** procedure if erroneous values are read.

This procedure permits to calibrate automatically all amplifiers of the instrument.

Calibrate

1. Leave the conductivity input open.
2. Switch the instrument ON while holding **CAL** pressed.
3. The display shows [SELF CAL] for a while and returns to the measurements.

This procedure resets the instrument to the original factory settings. All personal memorised data (including access code) will be erased!

Reset

1. Switch the instrument ON while holding **MODE** pressed.
2. The display shows [rESEt] for a while and returns to the measurements.

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