CONSORT

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

P107

January 1997

Table of contents

Introduction	1
Manufacturer	1
Warranty	1
Servicing	1
Specifications	2
Error Codes	2
Inputs	2
Recorder output	2
pH-MEASUREMENT	3
Maintenance	3
pH-measurement	3
mV-MEASUREMENT	4
Maintenance	4
mV-measurement	4
Dead-stop titration	5
Temperature-measurement	5
TROUBLE SHOOTING	5
Simple test	5
Calibrate	5

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This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 4° and 40°C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

Connect the instrument only to an earthen power line. The required power source is indicated on the label at the back of the instrument. Do not cut and splice the power cord. When removing the power cord from the wall outlet, be sure to unplug by holding the plug attachment and not by pulling the cord. Do not hold the plug by wet hand.

CONSORT nv
Parklaan 36
B2300 Turnhout
Belgium

Tel (++32)(14)41 12 79
Fax (++32)(14)42 91 79

Introduction

Manufacturer

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

CONSORT is not liable for consequential damages arising out of the use or handling of its products.

Warranty

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply and **NOT** to send the following items unless they are suspect:

Manual
Cables
Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

Servicing

Specifications

Ranges	pH 0...+14 pH mV ±1000 mV °C 0...100°C
Resolution	0.01 pH, 0.5°C, 1 mV
Inputs	BNC-1 input for pH electrodes BNC-2 input for mV electrodes BANANA input for a Pt1000 temperature probe
Temp. Comp.	automatic with Pt1000 or manual
Analogue output	0...2 V, ca 8 kΩ, 8 bit D/A
Display	4 digit/character, 13 mm LED display, shows readings, °C & messages
Keys	6 tactile membrane keys
Ambient temp.	4...40°C
Rel. Humidity	0...90 % (non-condensing!)
Power supply	210-250 VAC (* 105-125 VAC), 50/60 Hz, max. 6 VA
Cabinet	rugged plastic/metal cabinet
Dimensions	300 x 180 x 120 mm
Weight	2 kg

Each time the instrument is switched on, it will check briefly its internal calibrations while showing all functions of the display and the software-version. An error message appears when something is wrong.

[Or.pH]	= Overrange pH (e.g. interrupted electrode cable).
[Or.mV]	= Overrange mV (e.g. interrupted electrode cable)
[Or.°C]	= Overrange °C (measured degrees out of range).
[MEM]	= Fatal memory error (all factory settings are cleared. Have your instrument serviced!).

Error codes

A recorder can be connected to the red (+) and black (-) REC. terminals. The 0 to 2V output range corresponds with 0 to 14 pH in the pH mode or -1000 to +1000 mV in the mV mode. Use only laboratory recorders with a high input impedance!

Recorder output

The pH electrode should be connected to the BNC-1 connector. If separate electrodes are used, connect the reference electrode to the REF. terminal. Automatic temperature compensation and temperature measurements are possible by plugging a Pt1000 thermocompensator into the °C terminals. You can also use a combination pH electrode with built in Pt1000. Its banana plug should be inserted in the upper °C terminal. Without compensator, the manual temperature compensation is automatically switched on. The mV electrode should be connected to the BNC-2 connector.

Inputs

pH-MEASUREMENT

A pH electrode is active and stable only after wetting! For this purpose it must be immersed for **at least ten hours** in a 3...4 M KCl solution. During short interruptions (e.g. storage) the electrode should be immersed in a 3...4 M KCl solution. In doing this it is always kept ready for use. When the interruption is longer than a month, refill the closing cap with 3...4 M KCl and plug it on the electrode tip in order to protect the glass bulb. Before use, ensure that the reference part of the electrode is topped up with a 3...4 M KCl solution.

Avoid a low pressure inside the electrode! Therefore always remove the closure from the refilling aperture during the measurements as well as during the standardisation. This allows the saltbridge solution to flow through the ceramic liquid junction and prevents contamination of the electrolyte. For the same reason, the inside level should always be higher than the outside level of the measuring solution. Close the refilling aperture again when storing the electrode.

A polluted electrode may be cleaned with a soft detergent or 0.1 M HCl. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

If the electrode is polluted by proteinaceous materials (such as blood), it should stand in a pepsin solution overnight and then be cleaned before use. The pH electrode wears away by being used. If the electrode tends to respond slower and calibration becomes difficult, even after cleaning, it should be replaced by a new one.

A pH measurement is always relative which makes it necessary to calibrate regularly.

Zero point:

Each electrode has its own specific pH value, corresponding to the zero potential point (ISO pH). This specific pH value can differ from its theoretical zero point (e.g. 7 pH) by a few tenths of a pH.

Slope:

A good electrode follows the Nernst equation within 90 to 100 %. If the slope lies outside these limits, the electrode becomes slow and may not be optimal for most applications.

Maintenance

pH measurement

1. Select the **pH** mode.
2. Adjust the **TEMPERATURE** compensation to the temperature of the solutions (not necessary when a Pt1000 is used).
3. After rinsing the electrodes with distilled water, immerse them in a buffer solution of pH as close as possible to the ISO pH of the electrode (e.g. 6.87 pH).
4. Standardise to the pH of the buffer with the **CALIBRATION** button.
5. After rinsing the electrodes with distilled water, immerse them in a second buffer solution of pH close to the pH to be measured (e.g. 4.01 pH).
6. Standardise to the pH of the buffer with the **SLOPE** button.
7. After rinsing the electrodes with distilled water, immerse them in the solution to be measured and read the display.
8. Rinse the electrodes always with distilled water after use and store them in a 3.5 M KCl solution.

- * Stirring the solution during the measurements promotes the homogeneity and is therefore always recommended.
- * See page 21 for the standard buffers.

mV-MEASUREMENT

Metal electrodes (Pt, Ag, Au): Metal electrodes are always ready for use. During short interruptions they are immersed in distilled water. **They should be cleaned regularly:**

- * Silver electrodes are immersed in a concentrated ammonia solution during one hour.
- * Platinum or gold electrodes are immersed in concentrated nitric acid during one hour.

1. Select the mV range by pressing **MODE**.
2. After rinsing the electrodes with distilled water, immerse them in the solution to be measured and read the display.
3. Rinse the electrodes always with distilled water after use and store them in a 3..4 M KCl solution.

Maintenance

***mV
measurement***

This model is provided with 2 output terminals from which a constant current of ca. 10 μ A can be drawn. In order to drive this current through the electrodes proceed as follows.

Dead-stop titration

1. Select the **mV** mode.
2. Make a connection with the optional adapter A4023 between the mV input and the blue POL. terminal.
3. Connect the platinum electrodes with the REF. input and the blue POL. terminal.
4. After rinsing the electrodes with distilled water, immerse them in the solution to be titrated and read the display.
5. Rinse the electrodes always with distilled water after use and store them in that solution.

Provided a Pt1000 thermocompensator is used, read temperature by selecting the °C mode.

Temperature measurement

TROUBLE SHOOTING

Try this procedure if you suspect the electronic performance of the instrument.

Simple test

1. Short-circuit all pH/mV inputs. Leave all other inputs open.
2. Switch the instrument ON.
3. Select the mV ranges. The display should show $[0] \pm 1$ mV. Run the **Calibrate** procedure if erroneous values are read.
4. Select the pH ranges. The display should show $[7.00] \pm 1$ pH. Run the **Calibrate** procedure if erroneous values are read.

This procedure permits to calibrate automatically all amplifiers of the instrument.

Calibrate

1. Short-circuit the pH and mV inputs.
2. Turn all three potentiometers fully to the right (clockwise) and select the °C mode.
3. Switch the instrument on.
4. The display shows [SELF] for a while and returns to the measurements.

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