CONSORT

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

P400

January 1997

Table of contents

Introduction	1
Manufacturer	1
Warranty	1
Servicing	1
Error Codes	1
Maintenance	2
pH-measurement	2

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This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 4° and 40°C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

The batteries can be replaced by opening the bottom compartment of the cabinet. Avoid to use rechargeable batteries!

CONSORT nv
Parklaan 36
B2300 Turnhout
Belgium

Tel (++32)(14)41 12 79
Fax (++32)(14)42 91 79

Introduction

Manufacturer

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

CONSORT is not liable for consequential damages arising out of the use or handling of its products.

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Manual
Cables
Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

[Or] = Overrange (e.g. interrupted electrode cable).
[CAL] = Calibration failed (e.g. defective electrode or used buffers!).
[MEM] = Fatal memory error (all factory settings are cleared. Have your instrument serviced!).

Warranty

Servicing

Error codes

A pH electrode is active and stable only after wetting! For this purpose it must be immersed for **at least ten hours** in a 3...4 M KCl solution. During short interruptions (e.g. storage) the electrode should be immersed in a 3...4 M KCl solution. In doing this it is always kept ready for use. When the interruption is longer than a month, refill the closing cap with 3...4 M KCl and plug it on the electrode tip in order to protect the glass bulb. Before use, ensure that the reference part of the electrode is topped up with a 3...4 M KCl solution.

Avoid a low pressure inside the electrode! Therefore always remove the closure from the refilling aperture during the measurements as well as during the standardisation. This allows the saltbridge solution to flow through the ceramic liquid junction and prevents contamination of the electrolyte. For the same reason, the inside level should always be higher than the outside level of the measuring solution. Close the refilling aperture again when storing the electrode.

A polluted electrode may be cleaned with a soft detergent or 0.1 M HCl. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

If the electrode is polluted by proteinaceous materials (such as blood), it should stand in a pepsin solution overnight and then be cleaned before use. The pH electrode wears away by being used. If the electrode tends to respond slower and calibration becomes difficult, even after cleaning, it should be replaced by a new one.

1. Select the pH range by pressing **MODE**. The display will immediately show the measured pH according to the previous standardization. Should you want to restandardise, press **CAL**.
2. After rinsing the electrodes with distilled water, immerse them in the first buffer solution.
3. The display shows one of the 9 buffers in memory, e.g. [6.87], while the upper indicator is blinking. Select the proper buffer with **MODE** and press **CAL**. The instrument shows the measured buffer and will standardize automatically when readings are stable ([c] stops blinking).
4. After rinsing the electrodes with distilled water, immerse them in the second buffer solution.
5. The display shows one of the other buffers in memory, e.g. [4.01], while the second indicator is blinking. Select the proper buffer with **MODE** and press **CAL**. The instrument shows the measured buffer and will standardize automatically when readings are stable ([c] stops blinking).
6. After rinsing the electrodes with distilled water, immerse them in the samples and read the display.
7. Rinse the electrodes always with distilled water after use and store them a 3...4 M KCl solution.

**pH
measurement**

- * See page 9 for the buffers in memory.
- * Stirring the solution during the measurements promotes the homogeneity and is therefore always recommended.
- * Adjust the manual temperature compensation by pressing **MODE**. Then press **CAL**, select the desired temperature with **MODE**, and press **CAL** again.
- * The instrument will refuse automatic standardization when the electrode is unstable. Insufficient stirring or a worn electrode may be the cause.

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