CONSORT

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

P601

January 1997

Table of contents

Introduction	1
Manufacturer	1
Warranty	1
Servicing	1
Specifications	2
Keyboard	2
Error Codes	3
SET-UP	3
Beeper	3
Alarm	3
Printer	4
Identification	4
Code	4
DATA ACQUISITION	4
Logging data	5
Adding data	5
Recalling data	5
Printing data	6
Inputs	6
pH-MEASUREMENT	6
Maintenance	6
pH-measurement	7
TEMPERATUREMEASUREMENT	8
Calibration	8
Measurement	8
mV-MEASUREMENT	8
Maintenance	8
mV-measurement	8
DATA COMMUNICATION	8
Recorder output	8
GLP-protocol	8
Digital output	9
Printer connection	10
Computer connection	10
TROUBLE SHOOTING	11
Simple test	11
Calibrate	11
Reset	11

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This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 4° and 40°C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

If it is equipped with rechargeable batteries, the instrument should be connected permanently to the mains in order to keep the batteries in good condition. Do not store the instrument for longer than 1 month without recharging the batteries. The batteries can be replaced by opening the bottom compartment of the cabinet.

CONSORT nv
Parklaan 36
B2300 Turnhout
Belgium

Tel (++32)(14)41 12 79
Fax (++32)(14)42 91 79

Introduction

Manufacturer

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

CONSORT is not liable for consequential damages arising out of the use or handling of its products.

Warranty

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Manual
Cables
Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

Servicing

Specifications

Ranges	pH 0...14 pH mV ±2000 mV °C 0..100°C
Resolution	0.01 pH, 1 mV, 0.1°C
Inputs	BNC input for pH/mV electrodes BANANA input for a Pt1000 temperature probe
Temp. Comp.	automatic with Pt1000 or manual
Digital output	programmable RS232, 2400 b/s, for bi-directional communication with a computer or printer
Analogue output	0...2 V, ca 8 kΩ, 8 bit D/A
Display	8 digit/character, 12 mm LCD display, shows readings, °C & messages
Keys	10 tactile membrane keys
Ambient temp.	4...40°C
Rel. Humidity	0...90 % (non-condensing!)
Power supply	210-250 V~ (* 105-125 V~), 50/60 Hz, max. 2 VA & 4 NiCd batteries, size AA
Cabinet	IP65 cabinet
Dimensions	252 x 121 x 50 mm
Weight	600 g

Each time the instrument is switched on, it will check briefly its internal calibrations while showing all functions of the display and the software-version. An error message appears when something is wrong.

MODE	= Selects all modes or escapes from error traps, calibration procedures, etc..by returning to the original mode.
°C	= Reads temperature, when a Pt1000 is present, or manual temperature compensation.
CAL	= Starts or proceeds a calibration or a function.
▲/▼	= Button for entering a value or for selecting a function.
HOLD	= Holds display when measuring.
RES	= Changes resolution/units of the display.
PRINT	= Displayed value is printed through the RS232 output.
ON/OFF	= Switches the instrument on or off.

Keyboard

Error codes

- [Or. pH] = Overrange pH (e.g. interrupted electrode cable).
[Or. mV] = Overrange mV (e.g. interrupted electrode cable).
[Or. °C] = Overrange °C (measured degrees out of range).
[no dAtA] = Data-logging memory empty.
[not.CAL] = Calibration procedure interrupted (the previous calibration remains valid!).
[Err.SLP] = Slope not within 80...120 % (e.g. worn pH electrode or used buffers).
[Err.pHo] = ISO pH not within 6...8 pH (e.g. defective electrode).
[Err.CAL] = Calibration error (check electrodes and standard solutions).
[Err.Stb] = Electrode not stable.
[Err.COM] = General communication error.
[Err.Prt] = Communication error with the attached printer.
[Err.bAt] = Battery charger defective (disconnect the meter immediately from the mains, remove the batteries and have it serviced!).
[Err.COd] = Invalid code (enter the proper sequence of keys!).
[Err.MEM] = Fatal memory error (all factory settings are cleared. Have your instrument serviced!).
[ILLEGAL] = An unauthorised person has failed to repair the unit.
Warranty is elapsed from now on!

SET-UP

The built-in beeper can be programmed to sound each time a key-path is pressed.

Beeper

1. Select the set-up procedure by pressing **MODE**.
2. Press **CAL** until you can choose with \blacktriangle or \blacktriangledown between [bBEEP on] or [bBEEP OFF]. Decide whether or not the beeper should sound each time the keyboard is touched.
3. Press **CAL** to continue or press **MODE** to cancel.

An audible alert can be given when the displayed values stray outside a programmable high- and low-level.

Alarm

1. Select the set-up procedure by pressing **MODE**.
2. Press **CAL** until the display lets you choose with \blacktriangle or \blacktriangledown between [ALrM OFF], [ALrM PH], or [ALrM MV]. Select [ALrM OFF] if no alarm is required. Press **CAL**.
3. The display shows e.g. [L=400MV] or [L=6.00PH]. Select the desired minimum level with \blacktriangle or \blacktriangledown and press **CAL**. Each time the measurements are below this level, the beeper will sound an alarm.
4. The display shows e.g. [H=600MV] or [H=8.00PH]. Select the desired maximum level with \blacktriangle or \blacktriangledown . Each time the measurements are above this level, the beeper will sound an alarm.
5. Press **CAL** to continue or press **MODE** to cancel.

Printer

The measured values can be printed automatically at timed intervals or manually on operator command.

1. Select the set-up procedure by pressing **MODE**.
2. Press **CAL** until the display shows e.g. [Pt=120]. Select the desired time interval (0...9999 s) between the transmitted data to a printer, with \blacktriangle or \blacktriangledown . Pre-set to zero if no automatic printing is required.
3. Press **CAL** to continue or press **MODE** to cancel.

The instrument can be identified by e.g. a computer when a specific number is allocated to it.

Identification number

1. Select the set-up procedure by pressing **MODE**.
2. Press **CAL** until the display shows e.g. [Id.no=5]. Select the desired identification number (0...999) for your meter with \blacktriangle or \blacktriangledown . Pre-set to zero if this number should not be mentioned while printing.
3. Press **CAL** to continue or press **MODE** to cancel.

A private code can be programmed to avoid undesired access to the instrument.

Secret code

1. Select the set-up procedure by pressing **MODE**.
2. Press **CAL** until the display shows [Code on], [Code OFF], or [Code CAL]. Decide, with \blacktriangle or \blacktriangledown , whether the access to all future calibrations should be by entering a code (on) or not (off). Press **CAL** to continue or press **MODE** to cancel. Select [Code CAL] to enter your personal code.
3. The display shows [E=-----] while [=] is blinking. Enter your secret sequence of 5 keys to proceed. Press **CAL** to continue or press **MODE** to cancel.

DATA-ACQUISITION

Up to 100 values can be stored in the non-volatile memory of the instrument.

- [no LOG] : No data-logging.
- [SET LOG] : Permits to pre-set the data-logging parameters.
- [LSt LOG] : Permits to view the stored values on the display.
- [Prt LOG] : Permits to send the stored values to a printer or computer.
- [Cnt LOG] : Permits to continue with a next series of measurements, without erasing the previous data-logging.

Set log: Procedure to pre-set the data-logging parameters.

Logging data

1. Select the set-up procedure by pressing **MODE** and press **CAL**.
2. Choose [SEt LOG] with \blacktriangle or \blacktriangledown . Select [no LOG] if no data-logging is required and press **MODE** to return to the original mode.
3. Press **CAL** until you can choose with \blacktriangle or \blacktriangledown between [LOG PH] or [LOG MV]. Select the desired mode and press **CAL**.
4. The display shows e.g. [dt=100]. Select the desired time interval (0...9999 s) between the data-logging, with \blacktriangle or \blacktriangledown and press **CAL**. Pre-set to zero for manual data-logging, each time **PRINT** is pressed (the next points 5 to 7 are skipped and press **CAL**).
5. The display shows e.g. [n=40]. Select the desired number of values (1...100) to be data-logged with \blacktriangle or \blacktriangledown and press **CAL**.
6. The display shows e.g. [d=31.12.91] while the year is blinking. Program the desired starting date with \blacktriangle or \blacktriangledown and press **CAL**. Program similarly the month and the day.
7. The display shows e.g. [t=23.59] while the minutes are blinking. Program the desired starting time with \blacktriangle or \blacktriangledown and press **CAL**. Program similarly the hours.
8. The display shows [StArt ?]. Press **CAL** to start the data-logging according the previous settings (a blinking [L] appears) or press **MODE** to cancel. After the last value has been catalogued the display will show [End LOG]. When manual data-logging has been selected, press **PRINT** to put a next measurement into memory. In the meantime the display shows the logging-number e.g. [LOG 7]. Press **MODE** to cancel.

[Cnt LOG]: Permits to continue with a next series of measurements, without erasing the previous data-logging

Adding data

1. Select the set-up procedure by pressing **MODE** and press **CAL**.
2. Select [Cnt LOG] with \blacktriangle or \blacktriangledown to continue with a next series of measurements, without erasing the previous data-logging, provided the memory has enough capacity left and press **CAL**.

[Lst LOG]: permits to view the stored values on the display.

Recalling data

1. Select the set-up procedure by pressing **MODE** and press **CAL**.
2. Select [Lst LOG] with \blacktriangle or \blacktriangledown and press **CAL**.
3. The display shows the first data-logged value in memory (a blinking [n] appears). Press \blacktriangle or \blacktriangledown to view all other values.
4. While pressing **HOLD** the corresponding date is shown (not relevant when a manual or continued data-logging has been performed).
5. While pressing **RES** the corresponding time shown (or: series + log-number, e.g. [2-017], when a continued data-logging has been performed).
6. While pressing $^{\circ}\text{C}$ the corresponding temperature is shown.
7. Press **PRINT** to send the recalled data to a printer.
8. Press **MODE** to return to cancel.

[Prt LOG]: Permits to send the stored values to a printer or computer.

Printing data

1. Select the set-up procedure by pressing **MODE** and press **CAL**.
2. Choose [Prt LOG] with \blacktriangle or \blacktriangledown and press **CAL**
3. The display shows [PrintinG] while all data-logged values, including corresponding dates and times, are sent to the RS232 output in the following format, compatible with EXCEL™, LOTUS 123™,... The columns are separated by a TAB character.

Log#	VALUE	UNIT	C	H	D
001	7.34	ppm	23.7	13:19:57	29/01/95
002	7.08	ppm	23.8	13:21:57	29/01/95
003	6.71	ppm	23.8	13:22:57	29/01/95
004	6.49	ppm	23.9	13:23:57	29/01/95

The measuring electrode should be connected to the coaxial INPUT connector. If separate electrodes are used, connect the reference electrode to the REF. terminal. Automatic temperature compensation and temperature measurements are possible by plugging a Pt1000 thermocompensator into the °C terminals. You can also use a combination pH electrode with built in Pt1000. Its banana plug should be inserted in the upper °C terminal. Without compensator, the manual temperature compensation is automatically switched on.

Inputs

pH-MEASUREMENT

A pH electrode is active and stable only after wetting! For this purpose it must be immersed for **at least ten hours** in a 3...4 M KCl solution. During short interruptions (e.g. storage) the electrode should be immersed in a 3...4 M KCl solution. In doing this it is always kept ready for use. When the interruption is longer than a month, refill the closing cap with 3...4 M KCl and plug it on the electrode tip in order to protect the glass bulb. Before use, ensure that the reference part of the electrode is topped up with a 3...4 M KCl solution.

Maintenance

Avoid a low pressure inside the electrode! Therefore always remove the closure from the refilling aperture during the measurements as well as during the standardisation. This allows the saltbridge solution to flow through the ceramic liquid junction and prevents contamination of the electrolyte. For the same reason, the inside level should always be higher than the outside level of the measuring solution. Close the refilling aperture again when storing the electrode.

A polluted electrode may be cleaned with a soft detergent or 0.1 M HCl. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

If the electrode is polluted by proteinaceous materials (such as blood), it should stand in a pepsin solution overnight and then be cleaned before use. The pH electrode wears away by being used. If the electrode tends to respond slower and calibration becomes difficult, even after cleaning, it should be replaced by a new one.

A pH measurement is always relative which makes it necessary to calibrate regularly.

Zero point:

Each electrode has its own specific pH value, corresponding to the zero potential point (ISO pH). This specific pH value can differ from its theoretical zero point (e.g. 7 pH) by a few tenths of a pH.

Slope:

A good electrode follows the Nernst equation within 90 to 100 %. If the slope lies outside these limits, the electrode becomes slow and may not be optimal for most applications.

1. Select the pH range by pressing **MODE**. The display will immediately show the measured pH according to the previous standardisation. Should you want to restandardise, press **CAL**.
2. Rinse the electrodes with distilled water and immerse them in the first buffer solution.
3. The display shows one of the 9 buffers in memory (b1) or the manual input (M1) of user specified buffers, e.g. [b1=4.01] while [b1] or [M1] is blinking. Select the proper buffer with ▲ or ▼ and press **CAL**. If you wish to change the manual buffer, first press **RES** to reset to the proper value with ▲ or ▼ and then press **CAL**.
4. The instrument shows the measured buffer and will standardise automatically when readings are stable ([=] stops blinking). If the buffer has a different temperature, compensate the indicated value manually with ▲ or ▼ (not necessary when using a Pt1000).
5. Go on in the same way with all the next buffers and always press **CAL**. Press **MODE** to start the measurements.
6. After rinsing the electrodes with distilled water, immerse them in the samples and read the pH-value from the display.
7. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.

- * See page 45 for the buffers in memory.
- * A blinking decimal point warns you for unstable measurements. Wait to read the display!
- * Stirring the solution during the measurements promotes the homogeneity and is therefore always recommended.
- * The instrument will refuse automatic standardisation when the electrode is unstable. Insufficient stirring or a worn electrode may be the cause. Select a lower resolution for a more rapid standardisation or replace the electrode!
- * Press **RES** to change the resolution from 0.1 to 0.01 pH.

TEMPERATURE MEASUREMENT

1. Read temperature by pressing °C.
2. Without Pt1000, adjust the manual temperature compensation with ▲ or ▼ and proceed by pressing **MODE**.

Temperature measurement

1. Press **CAL** to calibrate the °C scale with any reference. Immerse the electrode in a solution of known temperature.
2. The display shows e.g. [°C=22.3]. Calibrate to the proper value with ▲ or ▼ and press **CAL** again or press **RES** to reset the absolute readings.

Calibration

mV-MEASUREMENT

Metal electrodes (Pt, Ag, Au): Metal electrodes are always ready for use. During short interruptions they are immersed in distilled water. **They should be cleaned regularly:**

Maintenance

- * Silver electrodes are immersed in a concentrated ammonia solution during one hour.
- * Platinum or gold electrodes are immersed in concentrated nitric acid during one hour.

1. Select the mV range by pressing **MODE**.
2. After rinsing the electrodes with distilled water, immerse them in the solution to be measured and read the display.
3. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.

mV measurement

DATA COMMUNICATION

A recorder can be connected to the (+) and (-) terminals. The 0 to 2 V output range corresponds with 0 to 14 pH in the pH mode, -2000 to +2000 mV in the mV mode. Use only laboratory recorders with a high input impedance!

Recorder output

1. Select the pH range by pressing **MODE**.
2. Check the ISO pH and SLOPE of the electrode by pressing ▼ sequentially: e.g. [100.0 %] & [7.00 PHo].
3. Press **PRINT** to print a complete calibration report, e.g.:

GLP-protocol

```

CONSORT pH Meter P601
-----

SETTINGS
-----
Identification No      : 004
Software Version      : 3.0
User Code              : off
ATC                    (Pt1000) : yes
ATC Cal.Temp. (C)     : 30.8

Alarm pH              : on
Alarm High (pH)       : 12.86
Alarm Low (pH)        : 9.11

Batt.Capacity (%)     : 87

CALIBRATION
-----
Buffer-1 (pH)         : 6.87
Temperature (C)       : 19.4
Response Time (s)     : 9

Buffer-2 (pH)         : 4.01
Temperature (C)       : 19.5
Response Time (s)     : 13

Zero Point (pH)       : 6.91
Slope (%)             : 97.8
Average Resp. (s)     : 12

STATISTICS
-----
Zero Point (pH)       : -0.09
Slope (%)             : +2.8
Average Resp. (s)     : +1

```

A standard RS232 output terminal (DP9) is provided for interfacing the instrument with a printer or computer. Data is sent in the ASCII code at a BAUD rate of 2400 b/s (8 bit, no parity, 1 start and 2 stopbits, automatic line feed ON, data transfer protocol Xon/Xoff).

Digital output

Serial port pinout specifications:

pin 1 : RLSD, received line signal detector
pin 2 : TxD, transmit data
pin 3 : RxD, receive data
pin 4 : not connected
pin 5 : Gnd, signal ground
pin 6 : DSR, data set ready
pin 7 : not connected
pin 8 : CTS, clear to send
pin 9 : not connected

Interfacing with a printer:

Press **PRINT** each time you want to print the measured values. The columns are separated by a TAB character. Run the set-up procedure for automatic printing at timed intervals, (see 'Printer' page 4), while you are measuring. The **PRINT** keypath does not function in this case. It is also disabled during an automatic data-logging., e.g.:

Log#	VALUE	UNIT	C	H	D
001	7.34	ppm	23.7	13:19:57	29/01/95
002	7.08	ppm	23.8	13:21:57	29/01/95
003	6.71	ppm	23.8	13:22:57	29/01/95
004	6.49	ppm	23.9	13:23:57	29/01/95

Using WINDOWS™ to collect data:

1. Connect the computer with the RS232 terminal, using the optional AK2310 cable (eventually with adapter AB3190 for connection with a 25 pole RS232 terminal).
2. Start **WINDOWS™** and select *Terminal* in the *Accessories Group*.
3. Choose *Receive Text File* in the *Transfer Menu* and create a file name (TXT format) for the data to receive.
4. Select the *Settings Menu* and make sure the following parameters are programmed:

Baud Rate : 2400
Data bits : 8
Stop bits : 2
Parity : None
Flow control : Xon/Xoff

5. Start the instrument in the desired measuring mode for manual or automatic data collection (see 'Printer' page 4 for more details).
6. When finished, open the created file with a standard spreadsheet (e.g. **EXCEL™**) or wordprocessor (e.g. **WORD™**) for further processing.

Interfacing with a computer:

Send a single character to the instrument and it will execute the corresponding command as follows:

- 1 = press **MODE**.
- 2 = press **CAL**.
- 3 = press **▲**.
- 4 = press **▼**.
- 5 = press **SET**.
- 6 = press **HOLD**.
- 7 = press **RES**.
- 8 = press **PRINT**.
- ? = send display to computer.
- + = keyboard on (manual operation remains possible).
- = keyboard off (manual operation disabled).
- V = start the input of a value (only valid in certain routines).

Follow exactly the same measuring or calibration procedures, as described in this manual, to include any desired command in your computer data-acquisition program. This example continuously prints the measured values (pH and temperature) on the computer screen:

10	ON ERROR GOTO 70	'makes sure no error occurs
20	OPEN "COM1:2400,N,8,2" AS #1	'prepares computer input
30	PRINT #1,"8"	' PRINT button pressed
40	INPUT #1,V\$	'measured values transmitted
50	PRINT V\$	'received values printed
60	GOTO 30	'repeat procedure
70	RUN	

TROUBLE SHOOTING

Try this procedure if you suspect the electronic performance of the instrument.

Simple test

1. Short-circuit all pH/mV inputs. Leave all other inputs open.
2. Switch the instrument ON.
3. Select the mV ranges by pressing **MODE**. The display should show [0 MV] ± 1 mV. Run the **Calibrate** procedure if erroneous values are read.
4. Connect a fresh standard alkaline battery of 1.5 V to pH/mV input.
5. The display should now read approximately [1500 MV]. Run the **Reset** procedure if erroneous values are read.
6. Reverse the battery of 1.5 V and connect it to pH/mV input.
7. The display should now read approximately [-1500 MV]. Run the **Reset** procedure if erroneous values are read.

This procedure permits to calibrate automatically all amplifiers of the instrument.

Calibrate

1. Short-circuit all pH/mV inputs. Leave all other inputs open.
2. Switch the instrument ON while holding **CAL** pressed.
3. The display shows [SElFCAL] for a while and returns to the measurements.

This procedure resets the instrument to the original factory settings. All personal memorised data (including access code) will be erased!

Reset

1. Switch the instrument ON while holding **MODE** pressed.
2. The display shows [rESEt] for a while and returns to the measurements.

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