CONSORT

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

P603

January 1997

Table of contents

Introduction	1
Manufacturer	1
Warranty	1
Servicing	1
Specifications	2
Keyboard	2
Error Codes	3
SET-UP	3
Beeper	3
Recorder	3
Printer	4
Identification	4
Code	4
DATA ACQUISITION	4
Logging data	5
Adding data	5
Recalling data	5
Printing data	6
Inputs	6
pH-MEASUREMENT	6
Maintenance	6
pH-measurement	7
Relative pH-measurement	8
TEMPERATUREMEASUREMENT	8
Calibration	8
Measurement	8
mV-MEASUREMENT	9
Maintenance	9
Calibration	9
mV-measurement	9
Relative mV-measurement	9
ION-SELECTIVE MEASUREMENT	10
Calibration	11
Measurement	11
Direct concentration	11
Standard addition	11
Standard subtraction	12
Sample addition	13
Sample subtraction	13
DATA COMMUNICATION	14
Recorder output	14
GLP-protocol	14
Digital output	16
Printer connection	16
Computer connection	17
TROUBLE SHOOTING	18
Simple test	18
Calibrate	18
Reset	18

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This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 4° and 40°C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

If it is equipped with rechargeable batteries, the instrument should be connected permanently to the mains in order to keep the batteries in good condition. Do not store the instrument for longer than 1 month without recharging the batteries. The batteries can be replaced by opening the bottom compartment of the cabinet.

CONSORT nv
Parklaan 36
B2300 Turnhout
Belgium

Tel (++32)(14)41 12 79
Fax (++32)(14)42 91 79

Introduction

Manufacturer

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

CONSORT is not liable for consequential damages arising out of the use or handling of its products.

Warranty

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Manual
Cables
Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

Servicing

Specifications

Ranges	pH -2...+16 pH mV ±2000 mV Ion 10 ⁻¹¹ to 10 ⁺¹ , any unit °C -30...+130°C
Resolution	0.001 pH, 0.1 mV, 0.1°C, 1% Ion
Inputs	BNC input for pH/mV/Ion electrodes BANANA input for a Pt1000 temperature probe
Temp. Comp.	automatic with Pt1000 or manual
Digital output	programmable RS232, 2400 b/s, for bi-directional communication with a computer or printer
Analogue output	0...2 V, ca 8 kΩ, 8 bit D/A
Display	8 digit/character, 12 mm LCD display, shows readings, °C & messages
Keys	10 tactile membrane keys
Ambient temp.	4...40°C
Rel. Humidity	0...90 % (non-condensing!)
Power supply	210-250 V~ (* 105-125 V~), 50/60 Hz, max. 2 VA & 4 NiCd batteries, size AA
Cabinet	IP65 cabinet
Dimensions	252 x 121 x 50 mm
Weight	600 g

Each time the instrument is switched on, it will check briefly its internal calibrations while showing all functions of the display and the software-version. An error message appears when something is wrong.

MODE	= Selects all modes or escapes from error traps, calibration procedures, etc..by returning to the original mode.
°C	= Reads temperature, when a Pt1000 is present, or manual temperature compensation.
CAL	= Starts or proceeds a calibration or a function.
▲/▼	= Button for entering a value or for selecting a function.
HOLD	= Holds display when measuring.
RES	= Changes resolution/units of the display.
PRINT	= Displayed value is printed through the RS232 output.
ON/OFF	= Switches the instrument on or off.

Keyboard

Error codes

- [Or. pH] = Overrange pH (e.g. interrupted electrode cable).
[Or. mV] = Overrange mV (e.g. interrupted electrode cable).
[Or. °C] = Overrange °C (measured degrees out of range).
[no dAtA] = Data-logging memory empty.
[not.CAL] = Calibration procedure interrupted (the previous calibration remains valid!).
[Err.SLP] = Slope not within 80...120 % (e.g. worn pH electrode or used buffers).
[Err.pHo] = ISO pH not within 6...8 pH (e.g. defective electrode).
[Err.CAL] = Calibration error (check electrodes and standard solutions).
[Err.Stb] = Electrode not stable.
[Err.COM] = General communication error.
[Err.Prt] = Communication error with the attached printer.
[Err.bAt] = Battery charger defective (disconnect the meter immediately from the mains, remove the batteries and have it serviced!).
[Err.COd] = Invalid code (enter the proper sequence of keys!).
[Err.MEM] = Fatal memory error (all factory settings are cleared. Have your instrument serviced!).
[ILLEGAL] = An unauthorised person has failed to repair the unit.
Warranty is elapsed from now on!

SET-UP

The built-in beeper can be programmed to sound each time a key-path is pressed.

Beeper

1. Select the set-up procedure by pressing **MODE**.
2. Press **CAL** until you can choose with \blacktriangle or \blacktriangledown between [bBEEP on] or [bBEEP OFF]. Decide whether or not the beeper should sound each time the keyboard is touched.
3. Press **CAL** to continue or press **MODE** to cancel.

Before concentrations can be recorded properly, it is necessary to define the min./max. values corresponding to the left and right corners of the semi-logarithmic paper. A maximum of 4 decades is allowed.

Recorder

1. Select the set-up procedure by pressing **MODE**.
2. Press **CAL** until the display shows e.g. [rEc.L=1.-5]. Choose the lowest concentration corresponding to a 0 V recorder output with \blacktriangle or \blacktriangledown and press **CAL**.
3. The display shows e.g. [rEc.H=1.-2]. Choose the highest concentration corresponding to a 2 V recorder output with \blacktriangle or \blacktriangledown .
4. Press **CAL** to continue or press **MODE** to cancel.

Printer

The measured values can be printed automatically at timed intervals or manually on operator command.

1. Select the set-up procedure by pressing **MODE**.
2. Press **CAL** until the display shows e.g. [Pt=120]. Select the desired time interval (0...9999 s) between the transmitted data to a printer, with \blacktriangle or \blacktriangledown . Pre-set to zero if no automatic printing is required.
3. Press **CAL** to continue or press **MODE** to cancel.

The instrument can be identified by e.g. a computer when a specific number is allocated to it.

Identification number

1. Select the set-up procedure by pressing **MODE**.
2. Press **CAL** until the display shows e.g. [Id.no=5]. Select the desired identification number (0...999) for your meter with \blacktriangle or \blacktriangledown . Pre-set to zero if this number should not be mentioned while printing.
3. Press **CAL** to continue or press **MODE** to cancel.

A private code can be programmed to avoid undesired access to the instrument.

Secret code

1. Select the set-up procedure by pressing **MODE**.
2. Press **CAL** until the display shows [Code on], [Code OFF], or [Code CAL]. Decide, with \blacktriangle or \blacktriangledown , whether the access to all future calibrations should be by entering a code (on) or not (off). Press **CAL** to continue or press **MODE** to cancel. Select [Code CAL] to enter your personal code.
3. The display shows [E=-----] while [=] is blinking. Enter your secret sequence of 5 keys to proceed. Press **CAL** to continue or press **MODE** to cancel.

DATA-ACQUISITION

Up to 100 values can be stored in the non-volatile memory of the instrument.

- [no LOG] : No data-logging.
- [SET LOG] : Permits to pre-set the data-logging parameters.
- [LSt LOG] : Permits to view the stored values on the display.
- [Prt LOG] : Permits to send the stored values to a printer or computer.
- [Cnt LOG] : Permits to continue with a next series of measurements, without erasing the previous data-logging.

Set log: Procedure to pre-set the data-logging parameters.

Logging data

1. Select the set-up procedure by pressing **MODE** and press **CAL**.
2. Choose [SEt LOG] with \blacktriangle or \blacktriangledown . Select [no LOG] if no data-logging is required and press **MODE** to return to the original mode.
3. Press **CAL** until you can choose with \blacktriangle or \blacktriangledown between [LOG PH] or [LOG MV]. Select the desired mode and press **CAL**.
4. The display shows e.g. [dt=100]. Select the desired time interval (0...9999 s) between the data-logging, with \blacktriangle or \blacktriangledown and press **CAL**. Pre-set to zero for manual data-logging, each time **PRINT** is pressed (the next points 5 to 7 are skipped and press **CAL**).
5. The display shows e.g. [n=40]. Select the desired number of values (1...100) to be data-logged with \blacktriangle or \blacktriangledown and press **CAL**.
6. The display shows e.g. [d=31.12.91] while the year is blinking. Program the desired starting date with \blacktriangle or \blacktriangledown and press **CAL**. Program similarly the month and the day.
7. The display shows e.g. [t=23.59] while the minutes are blinking. Program the desired starting time with \blacktriangle or \blacktriangledown and press **CAL**. Program similarly the hours.
8. The display shows [StArt ?]. Press **CAL** to start the data-logging according the previous settings (a blinking [L] appears) or press **MODE** to cancel. After the last value has been catalogued the display will show [End LOG]. When manual data-logging has been selected, press **PRINT** to put a next measurement into memory. In the meantime the display shows the logging-number e.g. [LOG 7]. Press **MODE** to cancel.

[Cnt LOG]: Permits to continue with a next series of measurements, without erasing the previous data-logging

Adding data

1. Select the set-up procedure by pressing **MODE** and press **CAL**.
2. Select [Cnt LOG] with \blacktriangle or \blacktriangledown to continue with a next series of measurements, without erasing the previous data-logging, provided the memory has enough capacity left and press **CAL**.

[Lst LOG]: permits to view the stored values on the display.

Recalling data

1. Select the set-up procedure by pressing **MODE** and press **CAL**.
2. Select [Lst LOG] with \blacktriangle or \blacktriangledown and press **CAL**.
3. The display shows the first data-logged value in memory (a blinking [n] appears). Press \blacktriangle or \blacktriangledown to view all other values.
4. While pressing **HOLD** the corresponding date is shown (not relevant when a manual or continued data-logging has been performed).
5. While pressing **RES** the corresponding time shown (or: series + log-number, e.g. [2-017], when a continued data-logging has been performed).
6. While pressing **°C** the corresponding temperature is shown.
7. Press **PRINT** to send the recalled data to a printer.
8. Press **MODE** to return to cancel.

[Prt LOG]: Permits to send the stored values to a printer or computer.

Printing data

1. Select the set-up procedure by pressing **MODE** and press **CAL**.
2. Choose [Prt LOG] with \blacktriangle or \blacktriangledown and press **CAL**
3. The display shows [PrintinG] while all data-logged values, including corresponding dates and times, are sent to the RS232 output in the following format, compatible with EXCEL™, LOTUS 123™,... The columns are separated by a TAB character.

Log#	VALUE	UNIT	C	H	D
0001	7.348	pH	23.7	13:19:57	29/01/95
0002	7.215	pH	23.7	13:20:57	29/01/95
0003	7.088	pH	23.8	13:21:57	29/01/95
0004	6.712	pH	23.8	13:22:57	29/01/95
0005	6.491	pH	23.9	13:23:57	29/01/95

The measuring electrode should be connected to the coaxial INPUT connector. If separate electrodes are used, connect the reference electrode to the REF. terminal. Automatic temperature compensation and temperature measurements are possible by plugging a Pt1000 thermocompensator into the °C terminals. You can also use a combination pH electrode with built in Pt1000. Its banana plug should be inserted in the upper °C terminal. Without compensator, the manual temperature compensation is automatically switched on.

Inputs

pH-MEASUREMENT

A pH electrode is active and stable only after wetting! For this purpose it must be immersed for **at least ten hours** in a 3...4 M KCl solution. During short interruptions (e.g. storage) the electrode should be immersed in a 3...4 M KCl solution. In doing this it is always kept ready for use. When the interruption is longer than a month, refill the closing cap with 3...4 M KCl and plug it on the electrode tip in order to protect the glass bulb. Before use, ensure that the reference part of the electrode is topped up with a 3...4 M KCl solution.

Maintenance

Avoid a low pressure inside the electrode! Therefore always remove the closure from the refilling aperture during the measurements as well as during the standardisation. This allows the saltbridge solution to flow through the ceramic liquid junction and prevents contamination of the electrolyte. For the same reason, the inside level should always be higher than the outside level of the measuring solution. Close the refilling aperture again when storing the electrode.

A polluted electrode may be cleaned with a soft detergent or 0.1 M HCl. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

If the electrode is polluted by proteinaceous materials (such as blood), it should stand in a pepsin solution overnight and then be cleaned before use. The pH electrode wears away by being used. If the electrode tends to respond slower and calibration becomes difficult, even after cleaning, it should be replaced by a new one.

A pH measurement is always relative which makes it necessary to calibrate regularly.

Zero point:

Each electrode has its own specific pH value, corresponding to the zero potential point (ISO pH). This specific pH value can differ from its theoretical zero point (e.g. 7 pH) by a few tenths of a pH.

Slope:

A good electrode follows the Nernst equation within 90 to 100 %. If the slope lies outside these limits, the electrode becomes slow and may not be optimal for most applications.

Temperature correction:

The pH of each solution varies with the temperature. Therefore it is impossible to accurately compare pH readings unless they are corrected to a reference temperature. This model features a unique microprocessor controlled system, which helps you completely to overcome the problem of non linear temperature correction curves for the species to be measured. There is no need to keep sample temperatures close to each other, since a special temperature curve for any type of sample can be plotted in the non-volatile memory. Temperature corrected pH readings should not be confused with **the normal temperature compensation!**

1. Select the pH range by pressing **MODE**. The display will immediately show the measured pH according to the previous standardisation. Should you want to restandardise, press **CAL**.
2. Choose with \blacktriangle or \blacktriangledown between [Cor OFF], [Cor on] or [Cor CAL]. Decide whether or not (OFF) temperature corrected readings should be displayed (on), press **CAL** and proceed with point 4 (page 8). Select [Cor CAL] to plot the temperature curve of your sample in the memory of the instrument and press **CAL**.
3. The display shows e.g. [rEF=25°C] while [=] is blinking. Select the temperature to which all future measurements will be referred to with \blacktriangle or \blacktriangledown and press **CAL** to start the automatic plotting procedure or press **MODE** to cancel.

Prepare the sample solution to be measured in a thermostatic bath with magnetic stirrer and immerse the electrodes in it. Its temperature should be about 5°C lower than the lowest value you may possibly need during future measurements. The display shows e.g. [7.87cPH][18.9] which means that the sample has 7.87 pH at 18.9°C. A blinking [c], indicating that plotting has started, appears. Wait a while to enable the electrodes to adapt themselves to the new temperature.

Press **HOLD**, when readings are stable, to store the values in memory. Increase the temperature with at least 5°C and repeat the same procedure for as many temperatures as necessary. When the temperature has reached about 5°C more than the highest value you need for future measurements, press **CAL** to stop the temperature tracing. The display shows briefly [WAlt] while the curve is being computed.

4. Rinse the electrodes with distilled water and immerse them in the first buffer solution.
5. The display shows one of the 9 buffers in memory (b1) or the manual input (M1) of user specified buffers, e.g. [b1=4.01] while [b1] or [M1] is blinking. Select the proper buffer with ▲ or ▼ and press **CAL**. If you wish to change the manual buffer, first press **RES** to reset to the proper value with ▲ or ▼ and then press **CAL**.
6. The instrument shows the measured buffer and will standardise automatically when readings are stable ([=] stops blinking). If the buffer has a different temperature, compensate the indicated value manually with ▲ or ▼ (not necessary when using a Pt1000).
7. Go on in the same way with all the next buffers and always press **CAL**. Press **MODE** to start the measurements.
8. After rinsing the electrodes with distilled water, immerse them in the samples and read the pH-value from the display.
9. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.

1. While in the pH range, press ▲. At that moment the instrument stores the actual pH and automatically subtracts it from all next measurements. You can repeat this as many times as required. A blinking [r] in the display appears.
2. To cancel relative pH readings, press **MODE** and you will be in the normal measuring mode again.

Relative pH measurement

- * See page 73 for the buffers in memory.
- * A blinking decimal point warns you for unstable measurements. Wait to read the display!
- * Stirring the solution during the measurements promotes the homogeneity and is therefore always recommended.
- * The instrument will refuse automatic standardisation when the electrode is unstable. Insufficient stirring or a worn electrode may be the cause. Select a lower resolution for a more rapid standardisation or replace the electrode!
- * Press **RES** to change the resolution from 0.1 to 0.001 pH.

TEMPERATURE MEASUREMENT

1. Read temperature by pressing °C.
2. Without Pt1000, adjust the manual temperature compensation with ▲ or ▼ and proceed by pressing **MODE**.

Temperature measurement

1. Press **CAL** to calibrate the °C scale with any reference. Immerse the electrode in a solution of known temperature.
2. The display shows e.g. [°C=22.3]. Calibrate to the proper value with ▲ or ▼ and press **CAL** again or press **RES** to reset the absolute readings.

Calibration

mV-MEASUREMENT

Metal electrodes (Pt, Ag, Au): Metal electrodes are always ready for use. During short interruptions they are immersed in distilled water. **They should be cleaned regularly:**

- * Silver electrodes are immersed in a concentrated ammonia solution during one hour.
- * Platinum or gold electrodes are immersed in concentrated nitric acid during one hour.

1. Select the mV range by pressing **MODE**.
2. Press **CAL** to shift the mV scale with any reference. Immerse the electrodes in a standard solution of known potential .
3. The display shows e.g. [MV=476.3]. Calibrate to the proper value with \blacktriangle or \blacktriangledown and press **CAL** again or press **RES** to reset the absolute readings.

1. Select the mV range by pressing **MODE**.
2. After rinsing the electrodes with distilled water, immerse them in the solution to be measured and read the display.
3. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.

1. While in the mV range, press \blacktriangle . At that moment the instrument stores the actual mV and automatically subtracts it from all next measurements. You can repeat this as many times as required. A blinking [r] in the display appears.
2. To cancel relative mV readings, press **MODE** and you will be in the normal measuring mode again.

Maintenance

Calibration

mV measurement

Relative mV measurement

Press **RES** to change the resolution from 0.1 to 1 mV.

ION-SELECTIVE MEASUREMENT

Calibration

This instrument features a unique system requiring no conventional slope determination. It is performed in two phases, from which only one should be repeated regularly.

- * First of all, the relation between the electrode potential and a series of standards (up to 19 values) is plotted into the non volatile memory, after which a complex calibration curve is computed for recall later. At lower concentration levels ion selective electrodes tend to have a slow response making this procedure rather time consuming. However, because the curve remains in memory, this phase must not be repeated unless the electrode has been used for a longer period.
- * When using the direct concentration method, it is possible to calibrate on 2 points only, without changing the curve in memory. This quick procedure adapts the curve temporary to the actual measurements keeping the same curvature.

Always calibrate the electrode in the same conditions as future determinations will require and make sure temperature remains constant. The instrument only accepts the following standard values, in scientific notation: $1 \cdot 10^{-9}$ / $3 \cdot 10^{-9}$ / $1 \cdot 10^{-8}$ / $3 \cdot 10^{-8}$ / $1 \cdot 10^{-7}$ / $3 \cdot 10^{-7}$ / $1 \cdot 10^{-6}$ / $3 \cdot 10^{-6}$ / $1 \cdot 10^{-5}$ / $3 \cdot 10^{-5}$ / $1 \cdot 10^{-4}$ / $3 \cdot 10^{-4}$ / $1 \cdot 10^{-3}$ / $3 \cdot 10^{-3}$ / $1 \cdot 10^{-2}$ / $3 \cdot 10^{-2}$ / $1 \cdot 10^{-1}$ / $3 \cdot 10^{-1}$ / $1 \cdot 10^0$, any unit. Therefore, e.g. 1 mg/l is shown on the display as [1. 3] (= $1 \cdot 10^{-3}$).

1. Prepare 2 or more (up to 19) standards, covering your future measuring range. Each of them should be 3 to 10 times stronger than the previous one.
2. Press **MODE** until the display will show e.g. [EL03 PrG] or [EL03 CLr] depending on electrode 03 being programmed previously or being cleared. By pressing \blacktriangle or \blacktriangledown [EL01] through [EL10] can be recalled. Press **CAL** to start the automatic plotting procedure.
3. After rinsing the electrodes with distilled water, immerse them in the first standard solution. The display will show e.g. [3.-7 263.8] (= $3 \cdot 10^{-7}$ / 263.8 mV). Select with \blacktriangle or \blacktriangledown the first calibration value and press **CAL** when readings are stable.
4. After rinsing the electrodes with distilled water, immerse them in the second standard solution. The display will show e.g. [1.-6 206.5] (= $1 \cdot 10^{-6}$ / 206.5 mV). Select your next standard with \blacktriangle or \blacktriangledown and press **CAL** when readings are stable.
5. Go on in the same way with all the next standard solutions. After the last one has been entered, press **MODE** to stop the calibration procedure. The instrument will now compute all values to a curve. In the mean time the display shows briefly [WAIT].
6. After the display shows e.g. [EL03 PrG], which means that ion selective electrode no.03 is in the non volatile memory, select any mode to proceed.

Measurement

Pressing **RES** during the measuring phase, permits you to choose between readings in scientific notation or in ng/l, µg/l, mg/l or g/l.

- [Std.Add] : standard addition method
[Std.Sub] : standard subtraction method
[AnL Add] : sample addition method
[AnL Sub] : sample subtraction method

1. Press **MODE** until the display will show e.g. [EL04 PrG] or [EL04 CLr] depending on electrode 04 being programmed previously or being cleared. By pressing \blacktriangle or \blacktriangledown successively [EL01] through [EL10] can be recalled.
2. Select the desired ion selective electrode and verify if the display shows [PrG].
3. Press **CAL** and the display will immediately show the measured concentration according to the previous calibration.
4. Press **CAL** if a 2 point calibration is needed. This procedure adapts the memorised curve temporary to the actual measurements keeping the same curvature. The display will show e.g. [b1=10nG] or [b1=1.-8]. Press \blacktriangle or \blacktriangledown to recall the desired standard, then press **CAL** to proceed.
5. After rinsing the electrodes with distilled water, immerse them in the first standard solution. The display will show e.g. [E1=176.3]. When readings are stable press **CAL**.
6. The display will show e.g. [b2=100nG] or [b2=1.-7]. Press \blacktriangle or \blacktriangledown to recall the desired standard, then press **CAL** to proceed.
7. After rinsing the electrodes with distilled water, immerse them in the second standard solution. The display will show e.g. [E2=211.7]. When readings are stable press **CAL**.
8. Rinse the electrodes with distilled water, immerse them in the samples, and read the concentration on the display

Direct concentration

the standard addition method is used when the samples can be kept at the same temperature and a standard of 5 to 100 times the expected sample concentration is available. The most suitable increment can be found from: $\text{increment percentage} > 100 \times \text{expected sample value} / \text{standard value}$.

Standard addition

1. Press **MODE** until the display will show e.g. [EL04 PrG] or [EL04 CLr] depending on electrode 04 being programmed previously or being cleared. By pressing \blacktriangle or \blacktriangledown successively [EL01] through [EL10] can be recalled.
2. Select the desired ion selective electrode and verify if the display shows [PrG].
3. Press **MODE** and recall the desired method by pressing \blacktriangle or \blacktriangledown until the display shows [Std.Add]. Press **MODE** to proceed.
4. The display shows one of the possible increments (1, 2, 5, 10, or 20 %), e.g. [10% Incr]. Select the desired value with \blacktriangle or \blacktriangledown and press **CAL**.

5. Prepare a standard solution of the measured species. Its value should be 5 to 100 times higher than the expected sample concentration. The display shows one of the possible standards e.g. [Std 1.-3]. Select the desired value with \blacktriangle or \blacktriangledown and press **CAL**.
6. After rinsing the electrodes with distilled water, immerse them in e.g. 100 ml sample solution. The display will show the electrode potential e.g. [E=307.6]. When readings are stable press **CAL**.
7. Add a small increment corresponding to the chosen mode (e.g. 10 ml). The display will show the difference in electrode potential e.g. [dE=48.1]. When readings are stable press **CAL**. The instrument calculates the result and shows the concentration of the sample in the same unit as the standard used for the addition.
8. For a next measurement press **CAL** and the same procedure starts again from point 6). Press **MODE** to abort.

The standard subtraction method is used to measure species for which stable standards or ion selective electrodes do not exist. The procedure is similar to the standard addition method. The most suitable increment can be found from: increment percentage > 100 x expected sample value / standard value.

Standard subtraction

1. Calculate the stoichiometric ratio (n) of the reaction: $xS + yR \rightarrow S_xR_y$
 $n=x/y$ (S = sample, R = reagent).
2. Press **MODE** until the display will show e.g. [EL04 PrG] or [EL04 CLr] depending on electrode 04 being programmed previously or being cleared. By pressing \blacktriangle or \blacktriangledown successively [EL01] through [EL10] can be recalled.
3. Select the desired ion selective electrode and verify if the display shows [PrG].
4. Press **MODE** and recall the desired method by pressing \blacktriangle or \blacktriangledown until the display shows [Std.Sub]. Press **MODE** to proceed.
5. The display shows one of the possible increments (1, 2, 5, 10, or 20 %), e.g. [10% Incr]. Select the desired value with \blacktriangle or \blacktriangledown and press **CAL**.
6. Prepare a standard solution of the measured species. Its value should be 5 to 100 times higher than the expected sample concentration. The display shows one of the possible standards e.g. [Std 1.-5]. Select the desired value with \blacktriangle or \blacktriangledown and press **CAL**.
7. After rinsing the electrodes with distilled water, immerse them in e.g. 100 ml sample solution. The display will show the electrode potential e.g. [E=178.2]. When readings are stable press **CAL**.
8. Add a small increment corresponding to the chosen mode (e.g. 5 ml). The display will show the difference in electrode potential e.g. [dE=34.4]. When readings are stable press **CAL**. The instrument calculates the result and shows the concentration of the sample in the same unit as the standard used for the addition.
9. Calculate the real concentration by multiplying the result with the stoichiometric ratio (n). For a next measurement press **CAL** and the same procedure starts again from point 7). Press **MODE** to abort.

Sample addition

The analate addition method is used when the samples vary widely in temperature or have a 5 to 100 times higher concentration than the available standard. The most suitable increment can be found from: increment percentage < 100 x standard value / expected sample value.

1. Press **MODE** until the display will show e.g. [EL04 PrG] or [EL04 CLr] depending on electrode 04 being programmed previously or being cleared. By pressing \blacktriangle or \blacktriangledown successively [EL01] through [EL10] can be recalled.
2. Select the desired ion selective electrode and verify if the display shows [PrG].
3. Press **MODE** and recall the desired method by pressing \blacktriangle or \blacktriangledown until the display shows [Anl.Add]. Press **MODE** to proceed.
4. The display shows one of the possible increments (1, 2, 5, 10, or 20 %), e.g. [10% Incr]. Select the desired value with \blacktriangle or \blacktriangledown and press **CAL**.
5. Prepare a standard solution of the measured species. Its value should be 5 to 100 times lower than the expected sample concentration. The display shows one of the possible standards e.g. [Std 3.-6]. Select the desired value with \blacktriangle or \blacktriangledown and press **CAL**.
6. After rinsing the electrodes with distilled water, immerse them in e.g. 100 ml standard solution. The display will show the electrode potential e.g. [E= 281.5]. When readings are stable press **CAL**.
7. Add a small increment corresponding to the chosen mode (e.g. 20 ml). The display will show the difference in electrode potential e.g. [dE=23.9]. When readings are stable press **CAL**. The instrument calculates the result and shows the concentration of the sample in the same unit as the standard used for the addition.
8. For a next measurement press **CAL** and the same procedure starts again from point 6). Press **MODE** to abort.

The analate subtraction method is used to measure species for which stable standards or ion selective electrodes do not exist. The procedure is similar to the analate addition method. The most suitable increment can be found from: increment percentage < 100 x standard value / expected sample value.

Sample subtraction

1. Calculate the stoichiometric ratio (n) of the reaction: $xS + yR \rightarrow S_xR_y$
 $n=x/y$ (S = sample, R = reagent).
2. Press **MODE** until the display will show e.g. [EL04 PrG] or [EL04 CLr] depending on electrode 04 being programmed previously or being cleared. By pressing \blacktriangle or \blacktriangledown successively [EL01] through [EL10] can be recalled.
3. Select the desired ion selective electrode and verify if the display shows [PrG].
4. Press **MODE** and recall the desired method by pressing \blacktriangle or \blacktriangledown until the display shows [Anl.Sub]. Press **MODE** to proceed.
5. The display shows one of the possible increments (1, 2, 5, 10, or 20 %), e.g. [10% Incr]. Select the desired value with \blacktriangle or \blacktriangledown and press **CAL**.

6. Prepare a standard solution of the measured species. Its value should be 5 to 100 times lower than the expected sample concentration. The display shows one of the possible standards e.g. [Std 1.-6]. Select the desired value with ▲ or ▼ and press **CAL**.
7. After rinsing the electrodes with distilled water, immerse them in e.g. 100 ml standard solution. The display will show the electrode potential e.g. [E= 129.3]. When readings are stable press **CAL**.
8. Add a small increment corresponding to the chosen mode (e.g. 1 ml). The display will show the difference in electrode potential e.g. [dE= 41.7]. When readings are stable press **CAL**. The instrument calculates the result and shows the concentration of the sample in the same unit as the standard used for the addition.
9. Calculate the real concentration by multiplying the result with the stoichiometric ratio (n). For a next measurement press **CAL** and the same procedure starts again from point 7). Press **MODE** to abort.

DATA COMMUNICATION

A recorder can be connected to the (+) and (-) terminals. The 0 to 2 V output range corresponds with -2 to +16 pH in the pH mode, -2000 to +2000 mV in the mV mode or can be programmed for 1 to 4 decades in the concentration mode (run the set-up procedure 'Recorder' page 3). Use only laboratory recorders with a high input impedance!

Recorder output

1. Select the pH range by pressing **MODE**.
2. Check the ISO pH and SLOPE of the electrode by pressing ▼ sequentially: e.g. [100.0 %] & [7.00 PHo].
3. Press **PRINT** to print a complete calibration report, e.g.:

GLP-protocol

CONSORT Ion Analyser P603

SETTINGS

Identification No : 004
Software Version : 3.0
User Code : off
Relative pH : off

pH Correction : off
ATC (Pt1000) : no
Manual Temp. (C) : 23.5

Batt.Capacity (%) : 87

CALIBRATION

Buffer-1 (pH) : 6.87
Temperature (C) : 19.4
Response Time (s) : 9

Buffer-2 (pH) : 4.01
Temperature (C) : 19.5
Response Time (s) : 13

Zero Point (pH) : 6.91
Slope (%) : 97.8
Average Resp. (s) : 12

STATISTICS

Zero Point (pH) : -0.09
Slope (%) : +2.8
Average Resp. (s) : +1

Digital output

A standard RS232 output terminal (DP9) is provided for interfacing the instrument with a printer or computer. Data is sent in the ASCII code at a BAUD rate of 2400 b/s (8 bit, no parity, 1 start and 2 stopbits, automatic line feed ON, data transfer protocol Xon/Xoff).

Serial port pinout specifications:

pin 1 : RLSD, received line signal detector
pin 2 : TxD, transmit data
pin 3 : RxD, receive data
pin 4 : not connected
pin 5 : Gnd, signal ground
pin 6 : DSR, data set ready
pin 7 : not connected
pin 8 : CTS, clear to send
pin 9 : not connected

Interfacing with a printer:

Press **PRINT** each time you want to print the measured values. The columns are separated by a TAB character. Run the set-up procedure for automatic printing at timed intervals, (see 'Printer' page 4), while you are measuring. The **PRINT** keypath does not function in this case. It is also disabled during an automatic data-logging., e.g.:

Log#	VALUE	UNIT	C	H	D
0001	7.348	pH	23.7	13:19:57	29/01/95
0002	7.215	pH	23.7	13:20:57	29/01/95
0003	7.088	pH	23.8	13:21:57	29/01/95
0004	6.712	pH	23.8	13:22:57	29/01/95
0005	6.491	pH	23.9	13:23:57	29/01/95

Printer connection

Using WINDOWS™ to collect data:

1. Connect the computer with the RS232 terminal, using the optional AK2310 cable (eventually with adapter AB3190 for connection with a 25 pole RS232 terminal).
2. Start **WINDOWS™** and select *Terminal* in the *Accessories Group*.
3. Choose *Receive Text File* in the *Transfer Menu* and create a file name (TXT format) for the data to receive.
4. Select the *Settings Menu* and make sure the following parameters are programmed:

Baud Rate : 2400
Data bits : 8
Stop bits : 2
Parity : None
Flow control : Xon/Xoff

5. Start the instrument in the desired measuring mode for manual or automatic data collection (see 'Printer' page 4 for more details).
6. When finished, open the created file with a standard spreadsheet (e.g. **EXCEL™**) or wordprocessor (e.g. **WORD™**) for further processing.

Interfacing with a computer:

Send a single character to the instrument and it will execute the corresponding command as follows:

- 1 = press **MODE**.
- 2 = press **CAL**.
- 3 = press **▲**.
- 4 = press **▼**.
- 5 = press **SET**.
- 6 = press **HOLD**.
- 7 = press **RES**.
- 8 = press **PRINT**.
- ? = send display to computer.
- + = keyboard on (manual operation remains possible).
- = keyboard off (manual operation disabled).
- V = start the input of a value (only valid in certain routines).

Follow exactly the same measuring or calibration procedures, as described in this manual, to include any desired command in your computer data-acquisition program.

This example continuously prints the measured values (pH and temperature) on the computer screen:

10	ON ERROR GOTO 70	'makes sure no error occurs
20	OPEN "COM1:2400,N,8,2" AS #1	'prepares computer input
30	PRINT #1,"8"	' PRINT button pressed
40	INPUT #1,V\$	'measured values transmitted
50	PRINT V\$	'received values printed
60	GOTO 30	'repeat procedure
70	RUN	

TROUBLE SHOOTING

Try this procedure if you suspect the electronic performance of the instrument.

Simple test

1. Short-circuit all pH/mV inputs. Leave all other inputs open.
2. Switch the instrument ON.
3. Select the mV ranges by pressing **MODE**. The display should show [0 MV] ± 1 mV. Run the **Calibrate** procedure if erroneous values are read.
4. Connect a fresh standard alkaline battery of 1.5 V to pH/mV input.
5. The display should now read approximately [1500 MV]. Run the **Reset** procedure if erroneous values are read.
6. Reverse the battery of 1.5 V and connect it to pH/mV input.
7. The display should now read approximately [-1500 MV]. Run the **Reset** procedure if erroneous values are read.

This procedure permits to calibrate automatically all amplifiers of the instrument.

Calibrate

1. Short-circuit all pH/mV inputs. Leave all other inputs open.
2. Switch the instrument ON while holding **CAL** pressed.
3. The display shows [SELFCAL] for a while and returns to the measurements.

This procedure resets the instrument to the original factory settings. All personal memorised data (including access code) will be erased!

Reset

1. Switch the instrument ON while holding **MODE** pressed.
2. The display shows [rESEt] for a while and returns to the measurements.

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