# CONSORT

---

MANUAL

HANDLEIDING


MODE D'EMPLOI

ANLEITUNG

**P800**  
**K810**  
**C830**

---

April 2001

QMS Quality Management Systems Limited

## Registration Certificate

**T**his document certifies that the quality administration systems of  
**Consort NV**

Parklaan 36  
B-2300 Turnhout

have been assessed and approved by  
QMS Quality Management Systems Limited to the following  
quality administration systems, standards and guidelines:-  
NEN-EN-ISO 9001 : 1994

The approved quality administration systems apply to the following:-

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
Certificate Expiry: 24 Oktober 2005

Certificate Number: BE 4007

On behalf of  
**QMS Quality Management Systems Limited**

Issued by:  
Stichting Quality Management  
Systems Nederland

This Certificate remains valid while the holder maintains their quality administration systems in accordance with the standards and guidelines above, which will be audited by QMS Quality Management Systems Limited

This Certificate is the property of QMS Quality Management Systems Limited and must be returned in the event of cancellation.

**Instruments manufactured by CONSORT**

Electrophoresis power supplies

•

Ion analysers

•

pH meters

•

Conductometers

•

Dissolved oxygen meters

•

Thermometers

•

Industrial controllers

**CONSORT nv**  
Parklaan 36  
B2300 Turnhout  
Belgium

*Find out more at [www.consort.be](http://www.consort.be)*

Tel (+32)(14)41 12 79

Fax (+32)(14)42 91 79

E-mail: [info@consort.be](mailto:info@consort.be)

**DECLARATION OF CONFORMITY**

We declare under our sole responsibility that the product

**Multi-channel analysers**  
**content of the type numbers**  
**P800, K810, C830**

to which this declaration relates is in conformity

with the following standards

**EN61010**

**LOW VOLTAGE DIRECTIVE 73/23/EEG**

**EN50081-1**

**EN50082-1**

**EN60555-2**

**EMC DIRECTIVE 89/336/EEG**

Turnhout, January 7, 1999

A handwritten signature in black ink, appearing to be 'L. G. G.', is written below the date.

on behalf of **CONSORT nv**

**Table of contents**

|                               |   |
|-------------------------------|---|
| Introduction.....             | 5 |
| Manufacturer .....            | 5 |
| Warranty .....                | 5 |
| Servicing.....                | 5 |
| Specifications .....          | 6 |
| Keyboard .....                | 6 |
| Error codes .....             | 6 |
| AC adaptor .....              | 6 |
| Inputs.....                   | 6 |
| pH measurement.....           | 7 |
| Conductivity measurement..... | 7 |
| mV measurement.....           | 7 |
| Temperature measurement ..... | 7 |
| Maintenance .....             | 8 |

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## ***Introduction***

---

This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 4° and 40°C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

## ***Manufacturer***

---

### **CONSORT nv**

Parklaan 36  
B2300 Turnhout  
Belgium

Tel (++32)(14)41 12 79

Fax (++32)(14)42 91 79

E-mail: info@consort.be

## ***Warranty***

---

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

**CONSORT** is not liable for consequential damages arising out of the use or handling of its products.

## ***Servicing***

---

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

*Manual*

*Cables*

*Accessories*

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

---

|                      |  |               |                   |
|----------------------|--|---------------|-------------------|
| <b>Ranges</b>        | <b>pH</b>  | 0...14 pH     | <i>(not K810)</i> |
|                      | <b>mV</b>  | ±1000 mV      | <i>(not K810)</i> |
|                      | <b>Conductivity</b>                                | 0...100 mS/cm | <i>(not P800)</i> |
|                      | <b>°C</b>  | 0..100°C      |                   |
| <b>Resolution</b>    | 0.01 pH, 0.1 µS/cm, 1 mV, 0.5°C <i>(C830: 1°C)</i> |               |                   |
| <b>Inputs</b>        | BNC input for pH/mV electrodes                     |               |                   |
|                      | BNC input for conductivity cell                    |               |                   |
|                      | BANANA input for Pt1000 temperature probe          |               |                   |
| <b>Temp. Comp.</b>   | automatic with Pt1000 or manual                    |               |                   |
| <b>Display</b>       | 3½ digit, 12 mm LED display                        |               |                   |
| <b>Keys</b>          | 5 tactile membrane keys                            |               |                   |
| <b>Ambient temp.</b> | 4...40°C   |               |                   |
| <b>Rel. Humidity</b> | 0...90 % (non-condensing!)                         |               |                   |
| <b>Power supply</b>  | AC adapter for 230 V: model A4070                  |               |                   |
|                      | AC adapter for 115 V: model A4071                  |               |                   |
| <b>Cabinet</b>       | IP65 cabinet                                       |               |                   |
| <b>Dimensions</b>    | 13 x 18 x 10 cm                                    |               |                   |
| <b>Weight</b>        | 600 g  |               |                   |

---

**Keyboard**

- MODE** = Selects all modes or escapes from error traps, calibration procedures, etc. by returning to the original mode.
- CAL** = Starts or proceeds a calibration or a function.
- ▲/▼** = Button for entering a value or for selecting a function.
- ON/OFF** = Switches the instrument on or off.

---

**Error codes**

- [or]** = Overrange.
- [cc]** = Cell constant out of range.
- [CAL]** = Calibration failed.
- [MEM]** = Fatal memory error. Have your instrument serviced!

---

**AC adaptor**

Connect the jack of a suitable AC adapter (model **A4070** for 230 V~ or model **A4071** for 115 V~) to the DC socket. Do not hold the adapter by wet hand.

---

**Inputs**

The pH/mV electrode should be connected to the first BNC connector. The conductivity cell should be connected to the second BNC connector. Automatic temperature compensation and temperature measurements are possible by plugging a Pt1000 thermocompensator into the two black terminals. Without compensator, the manual temperature compensation is automatically switched on.

1. Select the pH range by pressing **MODE**. The display will immediately show the measured pH according to the previous calibration. Should you want to recalibrate, press **CAL**.
2. Rinse the electrodes with distilled water and immerse them in the first buffer solution.
3. The display shows one of the 9 buffers in memory, e.g. [4.01] while the pH indicator is blinking. Select the proper buffer with  $\blacktriangle$  or  $\blacktriangledown$  and press **CAL**. The instrument shows the measured buffer and will standardise automatically when readings are stable ([CAL] stops blinking).
4. After rinsing the electrodes with distilled water, immerse them in the second buffer solution.
5. The display shows one of the 9 buffers in memory, e.g. [9.18], while the pH indicator is blinking. Select the proper buffer with  $\blacktriangle$  or  $\blacktriangledown$  and press **CAL**. The instrument shows the measured buffer and will standardise automatically when readings are stable ([CAL] stops blinking).
6. After rinsing the electrodes with distilled water, immerse them in the samples and read the display.
7. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.

1. Select the conductivity-range by pressing **MODE**.
2. After rinsing the cell several times with a 0.01 M KCl standard solution (1413  $\mu\text{S}/\text{cm}$ ), immerse it in that standard. The solution temperature is not so critical but should lie between 0°C and 30°C. When no Pt1000 is used, do not forget to compensate manually first! Press **CAL**.
3. The display shows the reference temperature [r.20] or [r.25]. Select the desired value with  $\blacktriangle$  or  $\blacktriangledown$ . Press **CAL** to continue.
4. The display will show the computed cell constant e.g. [1.045] and will calibrate automatically when readings are stable ([CAL] stops blinking).
5. Rinse the cell several times with the sample, immerse it in that solution and read the display.
6. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

1. Select the mV range by pressing **MODE**.
2. After rinsing the electrodes with distilled water, immerse them in the solution to be measured and read the display.
3. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.

1. Read temperature by pressing **MODE**. Without Pt1000, adjust the manual temperature compensation with  $\blacktriangle$  or  $\blacktriangledown$  and proceed by pressing **MODE**.

**A pH electrode is active and stable only after wetting!** For this purpose it must be immersed for **at least ten hours** in a 3...4 M KCl solution. During short interruptions (e.g. storage) the electrode should be immersed in a 3...4 M KCl solution. In doing this it is always kept ready for use. When the interruption is longer than a month, refill the closing cap with 3...4 M KCl and plug it on the electrode tip in order to protect the glass bulb. Before use, ensure that the reference part of the electrode is topped up with a 3...4 M KCl solution.

**Avoid a low pressure inside the electrode!** Therefore always remove the closure from the refilling aperture during the measurements as well as during the calibration. This allows the saltbridge solution to flow through the ceramic liquid junction and prevents contamination of the electrolyte. For the same reason, the inside level should always be higher than the outside level of the measuring solution. Close the refilling aperture again when storing the electrode.

A polluted electrode may be cleaned with a soft detergent or 0.1 M HCl. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

If the electrode is polluted by proteinaceous materials (such as blood), it should stand in a cleaning solution overnight and then be cleaned with distilled water before use. The pH electrode wears away by being used. If the electrode tends to respond slower and calibration becomes difficult, even after cleaning, it should be replaced by a new one.

**Metal electrodes (Pt, Ag, Au):** Metal electrodes are always ready for use. During short interruptions they are immersed in distilled water. **They should be cleaned regularly:**

- Silver electrodes are immersed in a concentrated ammonia solution during one hour.
- Platinum or gold electrodes are immersed in concentrated nitric acid during one hour.

**A conductivity cell is active and stable only after wetting!** For this purpose it must be immersed for **at least one hour** in distilled water. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

A polluted cell may be cleaned with a soft detergent or diluted nitric acid. Greasy substances may be removed with acetone or alcohol (never do this with plastic electrodes!).

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