CONSORT

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

R321

January 1997

Index

Introduction	1
Manufacturer	1
Warranty	1
Servicing	1
Specifications	2
Keyboard	3
Error Codes	3
Inputs	3
Sample point	3
Connections	4
SET-UP	4
Channels	4
Relay	5
Recorder	5
Electrode test	5
Parameter	6
RS232	7
Identification	7
Code	7
MEASUREMENT	8
Theory	8
Maintenance	9
Measurements with 1-point calibration	9
Measurements with 2-point calibration	10
Manual O ₂ correction	10
Rapid °C shift	11
DATA COMMUNICATION	11
Diagnosis	11
Analogue output	11
Digital output	12
Computer connection	12
Printer connection	14
TROUBLE SHOOTING	14
Reset	14

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Introduction

This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 4° and 40°C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

Connect the instrument only to an earthed power line. The required power source is indicated on the label at the back of the instrument. Do not cut and splice the power cord. When removing the power cord from the wall outlet, be sure to unplug by holding the plug attachment and not by pulling the cord. Do not hold the plug by wet hand.

CONSORT nv
Parklaan 36
B2300 Turnhout
Belgium

Tel (++32)(14)41 12 79
Fax (++32)(14)42 91 79

Manufacturer

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

CONSORT is not liable for consequential damages arising out of the use or handling of its products.

Warranty

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Manual
Cables
Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

Servicing

Specifications

Ranges	O₂ 0...60 ppm Saturation 0...600 % Temperature 0...50°C
Resolution	0.1 ppm, 1%, 0.1°C or 0.1% of reading
Inputs	2 inputs for O ₂ -probe and Pt1000 (or Pt100) <ul style="list-style-type: none">- independent mode (operates as 2 independent controllers)- AND mode (relays will only close when both inputs exceed a pre-set alarm level)- OR mode (relays will close when one of the inputs exceeds a pre-set alarm level)- differential mode (reads and controls on the difference between both inputs)
Temp. Comp.	automatic with NTC (100 kΩ or 22 kΩ), built-in O ₂ -probe
Digital output	programmable RS232, 150...4800 b/s, for bi-directional communication between other R300's, computer or printer
Analogue output	channel-1: 4...20 mA, max. load 350 Ω, programmable 8 bit D/A with adjustable span channel-2: 4...20 mA, max. load 350 Ω, programmable 8 bit D/A with adjustable span
Relays	3 programmable solid-state relays for high /low/fail alarm, voltage-free contacts max. 250 V~/0.5 A
Display	8 digit/character, 12.5 mm LCD display, shows readings & messages
Keys	4 tactile membrane keys
Ambient temp.	4...40°C
Rel. Humidity	0...90 % (non-condensing!)
Power supply	210-250 V~, 50/60 Hz, max. 4 VA
Cabinet	splashproof (IP65) cabinet, for wall mounting with 4 screws
Dimensions	200 x 120 x 90 mm
Weight	1.2 kg
Safety	<ul style="list-style-type: none">- non volatile memory, protected against any kind of interference- resets automatically to measurement after power down, interference or unauthorised manipulation- returns automatically to measurement after power down, interference, unauthorised manipulation or when no key is pressed for 5 minutes

Each time the instrument is switched on, it will check briefly its internal calibrations while showing all functions of the display and the software-version. An error message appears when something is wrong.

- SET** = Installs the meter according your application.
▲/▼ = Button for entering a value or for selecting a function.
CAL = Starts or proceeds a calibration or a function.

Keyboard

The display alternates automatically (± 4 s interval) between the selected measuring ranges by pressing ▲ and ▼ together until [ALT] is shown.

- [Or. ppm]** = Overrange ppm O₂ (e.g. defective membrane).
[Or. °C] = Overrange °C (measured degrees out of range).
[not.CAL] = Calibration procedure interrupted (the previous calibration remains valid!).
[Err.CAL] = Calibration error (check electrodes and standard solutions).
[Err.COd] = Invalid code (enter the proper sequence of keys!).
[Err.COM] = General communication error.
[Err.Prb] = Temperature probe is not connected.
[EnGiNEEr] = Code access not possible (call a service engineer to reset the controller!).
[Err.MEM] = Fatal memory error (all factory settings are cleared. Have your instrument serviced!).
[ILLEGAL] = An unauthorised person has failed to repair the unit.
Warranty is elapsed from now on!

Error codes

Before starting a measurement, the O₂ electrodes should be connected to the proper inputs. An automatic temperature compensation is provided for each channel. Only electrodes, equipped with a built-in NTC thermistor of 100 k Ω or 22 k Ω (at 25°C) and supplying a saturation current between 30 & 100 nA (at 800 mV), are accepted.

Inputs

Before the exact location of the electrodes is finalised, it is worthwhile to carry out a measuring profile of the process stream to ensure that the required homogeneity is being achieved. As a general rule the electrodes should protrude into the process stream by one-third of the pipe width and be positioned so as to avoid gas pockets or sediments. Generally the measuring electrodes are inserted into the flow with the sample pick-up facing downstream in order to reduce the risk of blockage by suspended material. All sensors are usually positioned vertically at a suitable point downstream of the dosing point, where the measurements are stable and reliable. The response of the dosing system will be affected by the time required to disperse the reagents and their arrival at the measuring point.

Sample point

- O₂-1** Pt-cathode (inner conductor) to terminal (+) and Ag-anode (outer braid) to terminal (-) of (CHAN.1). NTC-thermistor (external wires) to terminals (+) and (-) of (°C1).
- O₂-2** Pt-cathode (inner conductor) to terminal (+) and Ag-anode (outer braid) to terminal (-) of (CHAN.2). NTC-thermistor (external wires) to terminals (+) and (-) of (°C2).

AC line:

phase to terminal (P).

neutral to terminal (N).

earth to earthing terminal.

Since the instrument is directly wired to the AC power line, an external mains switch should be provided when recommended by the local legislation. If a mains failure has occurred, it will always automatically proceed with the process control when the AC power comes back on.

Relays:

There are three relays (REL1) through (REL3) with voltage-free contacts, which may be programmed to be activated by the high/low levels or a pre-set alarm. The ratings of the connected load should lie between 100...250 V~, maximum 0.5 A.

SET-UP

Procedure to pre-set the different channels and their possibilities according your applications.

Channels

1. Press **SET**.
2. Press **▲** or **▼** until the display shows [SEt CHAn] and press **CAL** to continue or **SET** to cancel.
3. The display shows [CH1=on] or [CH1=OFF]. Select the desired mode for channel-1 with **▲** or **▼** and press **CAL** to continue or **SET** to cancel.
4. The display shows [CH2=on] or [CH2=OFF]. Select the desired mode for channel-2 with **▲** or **▼** and press **CAL** to continue or **SET** to cancel.
5. The display shows [Atc=100] or [Atc=22]. Select the proper thermistor (100 kΩ or 22 kΩ) of the electrodes with **▲** or **▼** and press **CAL** to continue or **SET** to cancel.
6. This step is skipped if one of the channels is not used. The display shows [CH.ind.], [CH.diF.], [CH.And], or [CH.or]. Select the desired mode with **▲** or **▼** and press **CAL** to continue or **SET** to cancel. Mode [CH.ind.] means that both channels operate independently. Mode [CH.diF.] means that the relays are only activated when a pre-set difference between both channels is exceeded. Mode [CH.And] means that both channels must exceed a pre-set level before the relays are activated. Mode [CH.or] means that the relays are activated when one of both channels exceeds a pre-set level.

Relay

All relays can be programmed for controlling any of the high or low level alarms or for connection in a general alarm situation (ALrM), e.g. when an electrode breaks. The displayed possibilities depend from the channel set-up.

1. Press **SET**.
2. Press \blacktriangle or \blacktriangledown until the display shows [SEt rEL] and press **CAL** to continue or **SET** to cancel.
3. The display shows e.g. [r1=PPM1(H)], [r1=PPM1(L)], [r1=PPM2(H)], [r1=PPM2(L)], or [r1=ALrM]. Select the desired mode for relay-1 with \blacktriangle or \blacktriangledown and press **CAL** to continue with all the next relays (r2...r3) or press **SET** to cancel.

Before measurements can be recorded properly, it is necessary to define the recorder-span (values corresponding to the left and right corners of the graphic paper). The displayed possibilities depend from the channel set-up.

Recorder

1. Press **SET**.
2. Press \blacktriangle or \blacktriangledown until the display shows [SEt rEc] and press **CAL** to continue or **SET** to cancel.
3. The display shows e.g. [rc1=PPM1], [rc1=PPM2],[rc1=°C1], or [rc1=°C2]. Select the desired channel & mode for recorder-1 with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to cancel.
4. The display shows e.g. [rc=3.0(L)]. Choose the value corresponding to a 4 mA recorder-1 output with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to cancel.
5. The display shows e.g. [rc=10.0(H)]. Choose the value corresponding to a 20 mA recorder-1 output with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to cancel.
6. The display shows e.g. [rc2=PPM1], [rc2=PPM2], [rc2=°C1], or [rc2=°C2]. Select the desired channel & mode for recorder-2 with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to cancel.
7. The display shows e.g. [rc=7.0(L)]. Choose the value corresponding to a 4 mA recorder-2 output with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to cancel.
8. The display shows e.g. [rc=13.0(H)]. Choose the value corresponding to a 20 mA recorder-2 output with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to cancel.

Allows to check the polarisation-current of the connected electrode.

Electrode test

1. Press **SET**.
2. Press \blacktriangle or \blacktriangledown until the display shows [SEtUP] and press **CAL** to continue or **SET** to cancel.
3. Press **CAL** until the display shows [teSt Cur]. Press **SET** if you want to test the polarisation-current of your oxygen electrode.
4. The current is shown in nA (800 mV applied to the cell) while a blinking [i] appears.
5. Press **SET** to cancel.

The displayed possibilities depend from the channel set-up and their logic function (independent, differential, AND, OR).

1. Press **SET**.
2. Press \blacktriangle or \blacktriangledown until the display shows e.g. [SEt PAR1] and press **CAL** to continue or **SET** to cancel. Select [SEt PAR2] to program the parameters of channel-2. Indicator (1) will now be replaced by (2). If both channels have been linked by one of the logical functions, both indicators will blink.
3. The display shows e.g. [1.7PPM] while (1) and (L) are blinking. Pre-set the desired low level for channel-1 with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to cancel.
4. The display shows e.g. [8.3PPM] while (1) and (H) are blinking. Pre-set the desired high level for channel-1 with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to cancel.
5. The display shows the hysteresis between the ON and OFF state of the relays e.g. [Hy=0.10] while (1) and (=) are blinking. Pre-set the desired hysteresis for channel-1 with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to cancel.
6. The display shows the delay (s) after which the relays should be activated, when one of the pre-set levels is exceeded, e.g. [dt=120] while (1) and (=) are blinking. Pre-set the desired delay for channel-1 with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to cancel.
7. The display shows the max. level difference from which proportional regulation should start, e.g. [dO2=0.5] while (1) and (=) are blinking. Pre-set the desired difference for channel-1 with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to cancel. Pre-set this value to zero for simple ON/OFF regulations.
8. The display shows the pulse time (s) of the relays, e.g. [Pt=15] while (1) and (=) are blinking. The relays will pulse at a rate proportional to the regulation difference. When the difference is superior to the pre-set max. value, the relay is continuously activated. However, the wait-time between the pulses increases near the high or low pre-set level. Pre-set the desired pulse-time for channel-1 with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to cancel. This step is skipped for simple ON/OFF regulations.
9. The display shows the alarm time (s) after which the process control should be interrupted in case of a failure, e.g. [At=360] while (1) and (=) are blinking. Pre-set the desired alarm-time for channel-1 with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to cancel. Pre-set this value to zero if no alarm is required.
10. When both channels are linked to each other by an AND or an OR function, the display shows the maximum allowed difference between the 2 measurements, e.g. [dA=0.3] while (1), (2), and (=) are blinking. Pre-set the desired value with \blacktriangle or \blacktriangledown and press **CAL** to continue or **SET** to cancel. The process control is immediately stopped and a general alarm is given should this value be exceeded.

The differential function between both channels permits to control a process on the difference between two measured values. A high level pre-set of e.g. 3 ppm causes the relevant relay to be activated only when the difference in ppm between both channels exceeds 3 ppm and NOT when one of the channels exceeds 3 ppm!

It is recommended to decrease the BAUD rate in case of difficulties with the digital connections, especially for long distances.

1. Press **SET**.
2. Press \blacktriangle or \blacktriangledown until the display shows [SEtuP] and press **CAL** to continue or **SET** to cancel.
3. Press **CAL** until the display shows e.g. [br=2400]. Pre-set the desired baud-rate (150...4800 b/s) with \blacktriangle or \blacktriangledown and press **CAL**.
4. The display shows e.g. [rS=60] which means that data is transmitted every 60 s during a normal measuring phase (no pre-set levels are exceeded). Pre-set the desired time interval (0...9999 s) with \blacktriangle or \blacktriangledown . Pre-set to zero if not required.
5. The display shows e.g. [rS=15], while the (H) & (L) indicators are blinking, which means that data is transmitted every 15 s when a pre-set level is exceeded. Pre-set the desired time interval (0...9999 s) with \blacktriangle or \blacktriangledown . Pre-set to zero if not required.
6. Press **CAL** to continue or **SET** to cancel.

The instrument can be identified by e.g. a computer when a specific number is allocated to it.

Identification

1. Press **SET**.
2. Press \blacktriangle or \blacktriangledown until the display shows [SEtuP] and press **CAL** to continue or **SET** to cancel.
3. Press **CAL** until the display shows e.g. [Id.no=5]. Pre-set the desired identification number (0...999) for your meter with \blacktriangle or \blacktriangledown . Pre-set to zero when not relevant.
4. Press **CAL** to continue or **SET** to cancel.

A private code can be programmed to avoid undesired access to the instrument.

Code

1. Press **SET**.
2. Press \blacktriangle or \blacktriangledown until the display shows [SEtuP] and press **CAL** to continue or **SET** to cancel.
3. Press **CAL** until the display shows [Code on], [Code OFF], or [Code CAL]. Decide whether the access to all future calibrations should be by entering a code (on) or not (off) and press **CAL**. Select [Code CAL] to enter your personal code.
4. The display shows [E=-----] while [=] is blinking. Enter your secret sequence of 5 keys to proceed.

If [Code on] has been programmed after a general reset without entering first a personal code, the factory pre-set (5 x SET) is automatically valid.

MEASUREMENT

Theory

Principle:

The oxygenmeter and its electrode function according to the Clark principle with platinum as cathode and silver as anode in an electrolyte cell to which a polarisation voltage of 800 mV is applied. Oxygen gas present in the electrolyte is reduced to OH ions at the platinum cathode. The resulting current is diffusion limited and therefore proportional to the oxygen concentration in the sample solution. This current is amplified, corrected, and displayed in ppm dissolved oxygen.

Electrode:

The anode and cathode of the electrode are immersed in a special buffer solution inside the electrode, to give high stability readings and short polarisation times. The electrolyte compartment is separated by a PTFE film from the sample solution. Only gases, but no ions, can diffuse through this film from the sample solution into the internal electrolyte.

Temperature compensation:

Speed of electrode reactions depends on temperature. Therefore, for a given oxygen concentration, the electrode current becomes greater at higher temperature, and conversely. This temperature effect is compensated with the NTC resistor incorporated in the electrode. For aqueous solutions the oxygen meter therefore always gives correct readings in ppm, independent of the actual temperature.

Interferences:

All substances which can diffuse through the membrane and for which 800 mV potential suffices for polarographic reduction, will be reduced in the electrode. This will give a corresponding current contribution, if they are present. Interference can be caused by ions entering the electrode through porous or mechanically damaged membranes and by diffusion of other reactive gases apart from oxygen, e.g. CO₂, Cl₂, SO₂, and H₂S. These substances react in undesired manner with the electrode.

Acidic or basic gases change the pH value of the electrolyte solution and thus disturb the reading, particularly when measuring small oxygen concentrations. High salt concentrations in the sample solution can falsify readings too. The oxygen meter therefore is provided with a salinity compensation for salt concentrations up to 4 %.

Always store the electrode in distilled water after use! The membrane may cease to function correctly after prolonged use or mechanical damage. This will falsify the readings and the electrolyte may drain out. To replace the membrane (see page 57), unscrew the measuring head and dismantle it further to remove the old membrane. Take a new membrane and mount it carefully, fill the measuring head with electrolyte solution, and then slowly screw it onto the electrode. The electrode is now ready for use again.

The electrode body consisting of the electrode shaft and the plug head has almost unlimited life and requires no maintenance. After prolonged use of the electrode in solutions containing high concentrations of oxygen, the silver anode may become deactivated by undesired deposits. An indication is that the electrode no longer responds correctly to calibration. In this case, screw off the measuring head and cautiously remove the precipitated deposits from the silver electrode surface with the aid of filter paper. The electrode is then ready again for operation after refilling the measurement head with new electrolyte solution.

Start this method only with a well polarised electrode!

1. Before making measurements, check that the membrane is undamaged and that there are no gas bubbles present in the electrolyte space.
2. Connect the electrode and switch the instrument on. After a few minutes, the electrode is polarised and ready for making measurements.
3. Select the ppm (mg/l) or % range by pressing **MODE**. The display will immediately show the measured oxygen according to the previous calibration. Should you want to restandardise, press **CAL**.
4. The display shows the salinity correction, e.g. [SAL=2.1%]. Leave this value to zero unless you are going to measure in heavily salted solutions such as e.g. sea-water (3.5%). Pre-set the proper value with ▲ or ▼ and press **CAL** or **MODE** to cancel.
5. The display shows the pressure correction, e.g. [P=1023hP]. Pre-set the actual air pressure (hP) with ▲ or ▼ and press **CAL** or **MODE** to cancel.
6. Press ▲ or ▼ until the display shows [1.Point]. Press **CAL** to start the calibration procedure.
7. The electrode exposed in this manner to the air reaches an equilibrium corresponding to the partial pressure of oxygen and thus to saturation in water at the given temperature. Alternatively, rinse the electrode with distilled water and then dip it into freshly air saturated distilled water for calibration. The advection rate must be at least 10 cm/s (**stir!**). Air saturated water can be prepared by shaking or vigorously stirring distilled water for 10 minutes with free access of air.
8. The instrument shows the measured saturation, e.g. [SAT.=8.3] and will calibrate automatically when readings are stable ([=] stops blinking).
9. The display shows [InSErt]. After rinsing the electrode with distilled water, insert it back in the system, press **CAL** and start the O₂ control.

Measurements with 1-point calibration

Measurements with 2-point calibration

1. Before making measurements, check that the membrane is undamaged and that there are no gas bubbles present in the electrolyte space.
2. Connect the electrode and switch the instrument on. After a few minutes, the electrode is polarised and ready for making measurements.
3. Select the ppm (mg/l) or % range by pressing **MODE**. The display will immediately show the measured oxygen according to the previous calibration. Should you want to restandardise, press **CAL**.
4. The display shows the salinity correction, e.g. [SAL=2.1%]. Leave this value to zero unless you are going to measure in heavily salted solutions such as e.g. sea-water (3.5%). Pre-set the proper value with \blacktriangle or \blacktriangledown and press **CAL** or **MODE** to cancel.
5. The display shows the pressure correction, e.g. [P=1023hP]. Pre-set the actual air pressure (hP) with \blacktriangle or \blacktriangledown and press **CAL** or **MODE** to cancel.
6. Press \blacktriangle or \blacktriangledown until the display shows [2.Point]. Press **CAL** to start the calibration procedure.
7. Dip the electrode into a saturated sodium sulphite solution (without stirring). The instrument shows the measured zero-current (nA) of the electrode, e.g. [O2.o=1.3] and will calibrate automatically when readings are stable ([=] stops blinking).
8. Rinse the electrode with distilled water and cautiously dry it without rubbing (use cotton wool or soft filter paper). The electrode exposed in this manner to the air reaches an equilibrium corresponding to the partial pressure of oxygen and thus to saturation in water at the given temperature. Alternatively, rinse the electrode with distilled water and then dip it into freshly air saturated distilled water for calibration. The advection rate must be at least 10 cm/s (**stir!**). Air saturated water can be prepared by shaking or vigorously stirring distilled water for 10 minutes with free access of air.
9. The instrument shows the measured saturation, e.g. [SAT.=8.3] and will calibrate automatically when readings are stable ([=] stops blinking).
10. The display shows [InSErt]. After rinsing the electrode with distilled water, insert it back in the system, press **CAL** and start the O₂ control.

Start this method only with a well polarised electrode!

1. Before making measurements, check that the membrane is undamaged and that there are no gas bubbles present in the electrolyte space.
2. Connect the electrode and switch the instrument on. After a few minutes, the electrode is polarised and ready for making measurements.
3. Select the ppm (mg/l) or % range by pressing **MODE**. The display will immediately show the measured oxygen according to the previous calibration. Should you want to restandardise, press **CAL**.
4. The display shows the salinity correction, e.g. [SAL=2.1%]. Leave this value to zero unless you are going to measure in heavily salted solutions such as e.g. sea-water (3.5%). Pre-set the proper value with \blacktriangle or \blacktriangledown and press **CAL** or **MODE** to cancel.
5. The display shows the pressure correction, e.g. [P=1023hP]. Pre-set the actual air pressure (hP) with \blacktriangle or \blacktriangledown and press **CAL** or **MODE** to cancel.

Manual O₂ correction

6. Press \blacktriangle or \blacktriangledown until the display shows [MAn.Cor]. Press **CAL** to start the calibration procedure.
7. The display shows e.g. [PPM=8.1] while [=] is blinking as long as readings are not stable. Shift to the real value with \blacktriangle or \blacktriangledown when [=] stops blinking, and press **CAL** to return to the measurements.

1. Select the desired °C-channel with \blacktriangle or \blacktriangledown and press **CAL**.
2. The display shows e.g. [°C=27.3] while [=] is blinking as long as readings are not stable. Shift to the real value with \blacktriangle or \blacktriangledown when [=] stops blinking, and press **CAL** to return to the measurements. Press **SET** here if you wish to return to the original readings. The display shows [rESEt?]. Press **CAL** to reset the shifted values or press **SET** to cancel.

Rapid °C shift

- * Read temperature by pressing \blacktriangle or \blacktriangledown .
- * Read % saturation by pressing \blacktriangle or \blacktriangledown .
- * The instrument will refuse automatic calibration when readings are unstable.

DATA COMMUNICATION

This procedure permits to recall some important values recorded during the measuring & calibration phases of the controller. These values are reset after each new calibration!

Diagnosis

1. Press **SET**.
2. Press \blacktriangle or \blacktriangledown until the display shows [diAGn] and press **CAL** to continue or **SET** to cancel.
3. The display shows the saturation-current (nA) of the electrode connected to channel-1, e.g. [Cur=86]. Press **CAL** to continue or **SET** to cancel.
4. The display shows the saturation-point (ppm) measured by the electrode connected to channel-1, e.g. [SAt=9.7]. Press **CAL** to continue or **SET** to cancel.
5. The display shows the saturation-current (nA) of the electrode connected to channel-2, e.g. [Cur=73]. Press **CAL** to continue or **SET** to cancel.
6. The display shows the saturation-point (ppm) measured by the electrode connected to channel-2, e.g. [SAt=9.8]. Press **CAL** to continue or **SET** to cancel.
7. The display shows the lowest value ever measured by channel-1, e.g. [Mi=6.3]. Press **CAL** to continue or **SET** to cancel.
8. The display shows the highest value ever measured by channel-1, e.g. [MA=8.0]. Press **CAL** to continue or **SET** to cancel.
9. The display shows the lowest value ever measured by channel-2, e.g. [Mi=3.1]. Press **CAL** to continue or **SET** to cancel.
10. The display shows the highest value ever measured by channel-2, e.g. [MA=3.8]. Press **CAL** to continue or **SET** to cancel.

Two recorders, one for each channel, can be connected to the (+) and (-) terminals (OUT1) and (OUT2). The 4 to 20 mA output range corresponds with the pre-set values (see 'Recorder' page 5).

Analogue output

A standard RS232 connection is provided for interfacing the instrument with a computer, printer or other R300's. Data is sent or received in the ASCII code at a programmable BAUD rate of 150...4800 b/s (8 bit, no parity, 1 start and 2 stopbits, automatic line feed ON, data transfer protocol X-ON/X-OFF).

TxD, transmit data to terminal (+) of (IN)(RS232).

TxD, signal ground to terminal (-) of (IN)(RS232).

RxD, receive data to terminal (+) of (OUT)(RS232).

RxD, signal ground to terminal (-) of (OUT)(RS232).

Interfacing with a computer:

Computer connection

Any desired command can be included in a computer data acquisition program. Follow exactly the same measuring or calibration procedures as described in this manual. The instrument is able to recognise commands in two ways:

Simple: the simple method can be used when only one controller is connected to the computer (identification number should be zero!). Send a single character to the instrument and it will immediately execute the corresponding command as follows:

- 1 = press **SET**.
- 2 = press **▲**.
- 3 = press **▼**.
- 4 = press **CAL**.
- ? = send display to computer.
- + = keyboard on (manual operation remains possible).
- = keyboard off (manual operation disabled).
- V = start the input of a value (only valid in certain routines).

This example makes the computer press **SET** in the simple mode:

```
10 OPEN "COM1:2400,N,8,2" AS #1 'prepares communication port
20 PRINT #1,"1"                'SET key pressed
```

Sophisticated: this method is used for all types of communications such as a network of several R300's connected to the computer. The commands should always be preceded by the identification number of the controller concerned (see Set-up):

#xxx = "#" + identification number (xxx) in ASCII code. All blank (ASCII-32) or CR (ASCII-13) characters following the identification number are ignored!

- C1 = press **SET**.
- C2 = press **▲**.
- C3 = press **▼**.
- C4 = press **CAL**.
- C? = send display to computer.
- C+ = keyboard on (manual operation remains possible).
- C- = keyboard off (manual operation disabled).
- V = start the input of a value (only valid in certain routines).
- LF = linefeed (ASCII-10), closes a command or input and permits the instrument to receive a next command.

This example makes the computer press **SET** on controller #7, in the sophisticated mode:

```
10 OPEN "COM1:2400,N,8,2" AS #1 'prepares communication port
20 PRINT #1,"#007";           'identification number 7
30 PRINT #1,"C1";             'SET key pressed
40 PRINT #1,CHR$(10);         'LF, command closed
```

Input of a value: the direct input of values ("V" command) is only possible when the instrument is running a routine where also a manual input of these values could be made. If not, it will ignore any "V" command from the computer!

Examples of allowed routines are:

- * temperature input (e.g. manual temperature compensation).
- * control parameter input such as Low, High, ...
- * time settings.
- * input of manual buffer values.

Examples of NOT allowed routines are:

- * during measurements or controls.
- * when values have to be selected rather than changed (e.g. buffer values in memory).
- * while an electrode is being calibrated.

A 16-bit value (2's complement) should be transmitted in the following sequence:

```
1st character = "V" (start a value input)
2nd character = highest byte of the value in ASCII
3rd character = lowest byte of the value in ASCII
4th character = checksum of 2nd and 3rd character in ASCII
5th character = LF, linefeed (ASCII-10)
```

When a correct checksum has been received, the instrument will send a confirmation character "!" to the computer. If not, a question mark "?" is sent. Both are eventually preceded by the identification number.

This example makes the computer to enter the value "1000" on controller #7, in the sophisticated mode:

$$1000 = (\text{H-byte } 3) \times 256 + (\text{L-byte } 232)$$

```
10 OPEN "COM1:2400,N,8,2" AS #1 'prepares communication port
20 PRINT #1,"#007";           'identification number 7
30 PRINT #1,"V";             'start a value input
40 PRINT #1,CHR$(3);         'highest byte of value
50 PRINT #1,CHR$(232);       'lowest byte of value
60 PRINT #1,CHR$(235);       'checksum of both bytes
70 PRINT #1,CHR$(10);        'LF, closes command
```

Interfacing with a printer:

Printer connection

Run the RS232 set-up procedure for automatic printing at timed intervals during the process control. The result looks e.g. as follows:

#005

7.5 ppm1, 21.4 °C (/H) (/P) (H=high-alarm, P=proportional)

8.5 ppm2, 54.3 °C (/L) (/O) (L=low-alarm, O=on/off)

#006

3.6 ppm1, 30.9 °C (/H) (/W) (H=high-alarm, W=wait-phase)

5.7 ppm2, 15.0 °C (/L) (/A) (L=low-alarm, A=general alarm)

TROUBLE SHOOTING

This procedure resets the instrument to the original factory settings.

Reset

1. Switch the instrument ON while pressing **SET**.
2. The display shows [Ec=----] while [=] is blinking. Enter the engineer-code (**SET**, **CAL**, **▼**, **SET**) to reset all parameters in the memory to the original manufacturer settings.
3. The display shows [rESEt] for a while and returns to the measurements.

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