Consort

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

R361
R362
R363

Table of contents

Introduction	1
Manufacturer	1
Warranty	1
Servicing.....	1
Keyboard	1
Mains connections	2
Outputs	2
Inputs	2
Language	3
Contrast	3
Channels.....	3
Date and time	4
Password	4
Data-logging	4
Analog output	4
RS485	4
Sample point	5
pH measurement.....	5
mV measurement	5
Temperature measurement	5
Free chlorine measurement.....	5
Conductivity measurement.....	6
Dissolved oxygen measurement	6
Air pressure measurement	6
GLP.....	7
Calibration reminder	7
Maintenance of electrodes	7

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Introduction

This instrument is manufactured with the latest technology and needs no particular maintenance. **Consort** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 0° and 40° C. Never use the instrument in a room with high humidity (>95 %) or at very low temperatures (condensation water!).

Manufacturer

Consort bvba

Parklaan 36

B2300 Turnhout

Belgium

Tel (+32)(14)41 12 79

Fax (+32)(14)42 91 79

E-mail: info@consort.be

Warranty

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **Consort** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **Consort** before returning any product for warranty repair on a freight prepaid basis!

Consort is not liable for consequential damages arising out of the use or handling of its products.

Servicing

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Manual

Cables

Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **Consort** dealer.

Keyboard


SET = Selects all settings or escapes from error traps, calibration procedures, etc..by returning to the original setting.

←→↑↓ = Button for entering a value or for selecting a function.

CAL = Starts or proceeds a calibration or a function.

STOP = Stops the control without interrupting the measurements.

HELP = Built-in manual.

Mains

Connect phase to terminal (P), neutral to terminal (N), and earth to earthing terminal.

Relays

There are four relays (REL1...REL4) with voltage-free contacts. Connect their two contacts in series with the load. The ratings of the connected load should lie between 12...250 VAC, maximum 1 A.


Inputs

The measuring electrodes should be connected to the coaxial connectors CH1 or CH2. Automatic temperature compensation and measurement are possible by connecting a Pt1000 thermocompensator to the coaxial connectors °C1 or °C2. Without compensator, the manual temperature compensation is automatically switched on.

- Avoid ground loops by immersing an earthing rod in the process liquid close to the electrodes. Connect it to terminal (MASS)!
- Always cover the unused input connectors with the supplied protection hoods!

Outputs

(R362, R363)


Recorder

Two recorders (4-20 mA), one for each channel, can be connected to the (+) and (-) terminals of (OUT1) and (OUT2). Max. impedance = 300 Ω.

RS485

Only use STP cable to connect up to 31 controllers to a computer!

Language

1. Press **SET**.
2. Select **[SYSTEM]** and press **CAL**.
3. Select **[LANGUAGE]** and press **CAL**.
4. Select the desired language and press **CAL**.

Contrast

1. Press **SET**.
2. Select **[SYSTEM]** and press **CAL**.
3. Select **[CONTRAST]** and press **CAL**.
4. Adjust the contrast and press **CAL**.

Channels

1. Press **SET**.
2. Select **[CHANNELS]** and press **CAL**.
3. Select the desired ranges **[INPUT]** and press **CAL**.
4. Select the desired inputs **[ATC]** (automatic temperature compensation) and press **CAL**.
5. Select **[CONTROL]** (yes or no) and press **CAL**.
6. Select the control type (on/off or proportional) and press **CAL**.
7. Pre-set the desired low level (Lo) and press **CAL**.
8. Select the corresponding relay (REL) and press **CAL**.
9. Pre-set the desired high level (Hi) and press **CAL**.
10. Select the corresponding relay (REL) and press **CAL**.
11. Pre-set the desired hysteresis (Hy) between the on and off state of the relays and press **CAL**.
12. Pre-set the desired delay (dt) after which the relays should be activated, when one of the pre-set levels is exceeded and press **CAL**.

Proportional control only:

13. Pre-set the desired the proportional control area (dP), in which proportional regulation should start and press **CAL**.
14. Pre-set the desired pulse-time (Pt) of the relays and press **CAL**. The relays will pulse at a rate proportional to the regulation difference.

Optional alarm procedure:

15. Select **[ALARM]** and press **CAL**.
16. The display shows the alarm time after which the process control should be interrupted in case of a failure. Pre-set the desired alarm-time and press **CAL**.
17. Select the corresponding relay (REL) to be closed in case of an alarm and press **CAL**.

Optional maintenance procedure:

18. Select **[MAINTENANCE]** and press **CAL**.
19. The display shows the interval for a regular maintenance. Pre-set the desired interval and press **CAL**.
20. Pre-set the desired waiting period before the beginning of the maintenance and press **CAL**.
21. Select the corresponding relay (REL) and press **CAL**.
22. Pre-set the maintenance time and press **CAL**.
23. Pre-set the desired waiting period after the finish of the maintenance and press **CAL**.
24. Pre-set the desired date and time when the first maintenance should start and press **CAL**.

- *The process control will automatically be interrupted during each maintenance period!*

1. Press **SET**.
2. Select [SYSTEM] and press **CAL**.
3. Select [DATE/TIME] and press **CAL**.
4. Select [MODIFY] and press **CAL**.
5. Adjust date and time and press **CAL**.

Password

1. Press **SET**.
2. Select [SYSTEM] and press **CAL**.
3. Select [PASSWORD] and press **CAL**.
4. A private code can be programmed to avoid undesired access to the instrument. Choose [YES] and press **CAL** to enter your secret sequence of 5 keys.

Data-logging

(R362, R363)

1. Press **SET**.
2. Select [DATA] and press **CAL**.
3. Select [INTERVAL] and press **CAL**.
4. Select the desired time interval between the data-logging and press **CAL**.

View the stored values:

1. Select [TABLE] and press **CAL**.
2. Select the desired channel by pressing **CAL**.
3. Select the first line to view and press **CAL**.
4. Select [SHOW TABLE] and press **CAL**.

Send the stored values:

1. Select [SEND] and press **CAL** to send the stored data to the digital port.

Erase the stored values:

1. Select [ERASE] and press **CAL**. Follow the instructions on the screen.

Analog output

(R362, R363)

1. Press **SET**.
2. Select [4-20 mA] and press **CAL**.
3. Select the desired range and press **CAL**.
4. Select the desired minimum level corresponding to a 4 mA recorder output and press **CAL**.
5. Select the desired maximum level corresponding to a 20 mA recorder output and press **CAL**.

RS485

(R362, R363)

1. Press **SET**.
2. Select [RS485] and press **CAL**.
3. Select [INTERVAL] and press **CAL**.
4. Select the desired time interval between the data transmission and press **CAL**.
5. Select [ID#] and press **CAL**.
6. Select the desired identification number (0...999) and press **CAL**. The instrument can be identified by e.g. a computer when a specific number is allocated to it.
7. Select [BAUD] and press **CAL**.
8. Select the desired baudrate and press **CAL**.

Sample point

Before the exact location of the electrodes is finalised, it is worthwhile to carry out a measuring profile of the process stream to ensure that the required homogeneity is being achieved. As a general rule the electrodes should protrude into the process stream by one-third of the pipe width and be positioned so as to avoid gas pockets or sediments.

Generally the measuring electrodes are inserted into the flow with the sample pick-up facing downstream in order to reduce the risk of blockage by suspended material.

All sensors are usually positioned vertically at a suitable point downstream of the dosing point, where the measurements are stable and reliable. The response of the dosing system will be affected by the time required to disperse the reagents and their arrival at the measuring point.

pH measurement

1. Select the desired mode with \leftrightarrow . The display will immediately show the measured value according to the previous calibration. Should you want to recalibrate, press **CAL**.
 2. The display shows five of the 9 buffers in memory (1.68, 4.00, 4.01, 6.87, 7.00, 9.18, 9.21, 10.01, 12.45). Select the proper values and press **CAL**. The unused buffers should be switched off.
 3. Rinse the electrodes with distilled water and immerse them in the first buffer solution. Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
 4. After rinsing the electrodes with distilled water, immerse them in the samples and read the display.
 5. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.
- Press \downarrow to change the resolution from 0.1 to 0.001 pH.

mV measurement

1. Select the desired mode with \leftrightarrow . The display will immediately show the measured value according to the previous calibration. Should you want to recalibrate, press **CAL**.
 2. Rinse the electrodes with distilled water and immerse them in a standard solution of known potential. Select [CALIBRATE], press **CAL**, adjust to the correct value and press **CAL**.
- Press \downarrow to change the resolution from 1 to 0.1 mV.
 - Select [RESET] and press **CAL** to reset the calibration.

Temperature measurement

1. Select the desired mode with \leftrightarrow . Without Pt1000, adjust the manual temperature compensation and press **CAL**. Should you want to recalibrate, press **CAL**.
 2. Immerse the Pt1000 in a standard solution of known temperature. Select [CALIBRATE], press **CAL**, adjust to the proper value and press **CAL**.
- Select [RESET] and press **CAL** to reset the calibration.

Free chlorine measurement

(R363)

1. Select the desired mode with \leftrightarrow . The display will immediately show the measured value according to the previous calibration. Should you want to recalibrate, press **CAL**.
2. Let the electrodes rest in the water stream for at least 24 hours (all dosing systems switched off). Measure the free chlorine accurately with a suitable photometer immediately before calibration. Select [CALIBRATE], press **CAL**, adjust to the proper value and press **CAL**.

Measurement:

1. Select the desired mode with \leftrightarrow . The display will immediately show the measured value according to the previous calibration. Should you want to recalibrate, press **CAL**.
2. The display shows the three standards in memory (0.01, 0.1, 1 M KCl). Select the proper values and press **CAL**. The unused standards should be switched off.
3. Select the temperature to which all future measurements will be referred to (REF: 25°C) and press **CAL**.
4. After rinsing the electrode several times with the first standard solution, immerse it in that standard. The solution temperature is not so critical but should lie between 0°C and 30°C. When no Pt1000 is used, do not forget to compensate manually first! Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
5. Rinse the electrode several times with the sample, immerse it in that solution and read the display.
6. Rinse the electrode always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

Capacitive compensation:

1. Capacitive compensation increases accuracy in the very low conductivity ranges (<10 $\mu\text{S}/\text{cm}$). Select whether or not this compensation should be applied. Verify if the attached electrode is completely dry and press **CAL**. Select [CAP.COMP] (yes or no) and press **CAL**. Follow the instructions on the screen.

Dissolved oxygen measurement

1. Select the desired mode with \leftrightarrow . The display will immediately show the measured value (in ppm or %O₂) according to the previous calibration. Should you want to recalibrate, press **CAL**.
 2. The display shows the salinity correction. Leave salinity correction to zero unless you are going to measure in heavily salted solutions such as e.g. seawater (35 g/l). Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
 3. The electrode exposed to the air reaches an equilibrium corresponding to the partial pressure of oxygen and thus to saturation in water at the given temperature. The instrument shows the measured saturation, current, temperature, and will calibrate automatically when readings are stable.
 4. After rinsing the electrodes with distilled water, immerse them in the samples and read the display. Stirring the solution during the measurements promotes the homogeneity and is **obligatory!** The advection rate must be at least 10 cm/s.
 5. Rinse the electrode always after use and store it in distilled water.
- Press \downarrow to change the resolution from 0.1 to 0.01 ppm or from 1 to 0.1 %.

Air pressure measurement:

1. Select the desired mode with \leftrightarrow . The display will immediately show the measured value according to the previous calibration. Should you want to recalibrate, press **CAL**.
 2. Select [CALIBRATE], press **CAL**, adjust to the proper value and press **CAL**.
- Select [RESET] and press **CAL** to reset the calibration.

1. Select the desired mode with \leftrightarrow and then press **CAL**.
2. Select [GLP] and press **CAL**.
3. Select [SHOW REPORT] and press **CAL**. Browse with \leftrightarrow \updownarrow to show a complete calibration report.
4. Select [SEND REPORT] and press **CAL** to send the report to the digital port.

Calibration reminder

1. Select the desired mode with \leftrightarrow and then press **CAL**.
2. Select [GLP] and press **CAL**.
3. Select [INTERVAL] and press **CAL**.
4. Select the desired time interval between each automatic warning for a new calibration of the electrodes and press **CAL**.

Maintenance of electrodes

A pH electrode is active and stable only after wetting! For this purpose it must be immersed for **at least ten hours** in a 3...4 M KCl solution. During short interruptions (e.g. storage) the electrode should be immersed in a 3...4 M KCl solution. In doing this it is always kept ready for use. When the interruption is longer than a month, refill the closing cap with 3...4 M KCl and plug it on the electrode tip in order to protect the glass bulb. Before use, ensure that the reference part of the electrode is topped up with a 3...4 M KCl solution. A polluted electrode may be cleaned with a soft detergent or 0.1 M HCl. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**). The pH electrode wears away by being used. If the electrode tends to respond slower and calibration becomes difficult, even after cleaning, it should be replaced by a new one.

Metal electrodes (Pt, Ag, Au): Metal electrodes are always ready for use. During short interruptions they are immersed in distilled water. **They should be cleaned regularly:**

- Silver electrodes are immersed in a concentrated ammonia solution during one hour.
- Platinum or gold electrodes are immersed in concentrated nitric acid during one hour.

A conductivity cell is active and stable only after wetting! For this purpose it must be immersed for **at least one hour** in distilled water. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition). A polluted cell may be cleaned with a soft detergent or diluted nitric acid. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

A dissolved oxygen electrode is active and stable only when polarised! The electrode body has almost unlimited life and requires no maintenance. After prolonged use of the electrode, it may become deactivated. An indication is that the electrode no longer responds correctly to calibration. In this case:

1. Carefully screw off the measuring head.
 2. Cautiously remove the precipitated deposits with the aid of filter paper. **Do not use grinding paper or a glass fibre!**
 3. Rinse several times with distilled water and shake carefully off the water drops.
 4. Take the measuring head, fill it with new electrolyte solution and very slowly screw it onto the electrode while holding it vertically.
 5. **Let the electrode rest for a few hours!** The electrode is now ready for use again.
- *Replace the membrane only when damaged, not when calibration is no longer possible!*

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Specifications	R361	R362	R363
pH	Range	-2...+16 pH	-2...+16 pH
	Resolution	0.001 pH	0.001 pH
	Accuracy	0.1% ± 1 digit	0.1% ± 1 digit
	Calibration	1...5 points	1...5 points
	Buffers	9 pre-programmed	9 pre-programmed
	Temperature Compensation	-30...+130 °C	-30...+130 °C
	ISO-pH	6...8 pH	6...8 pH
Slope	80...120%	80...120%	
mV	Range	±2000 mV	±2000 mV
	Resolution	0.1 mV	0.1 mV
	Accuracy	0.1% ± 1 digit	0.1% ± 1 digit
	Calibration	1 point	1 point
CONDUCTIVITY	Range	0...2000 mS/cm	0...2000 mS/cm
	Resolution	0.001 µS/cm	0.001 µS/cm
	Accuracy	0.5% f.s. of range	0.5% f.s. of range
	Calibration	1...3 points	1...3 points
	Standards	0.01/0.1/1M KCl	0.01/0.1/1M KCl
	Cell Constant	0.05...15 cm ⁻¹	0.05...15 cm ⁻¹
	Temperature Compensation	-30...+130 °C	-30...+130 °C
	Reference Temperature	20° or 25 °C	20° or 25 °C
	Temperature Coefficient	natural waters (EN27888)	natural waters (EN27888)
	Capacitive Compensation	✓	✓
DISSOLVED	Range	0...60 mg/l (0...600%)	0...60 mg/l (0...600%)
	Resolution	0.01 mg/l (0.1%)	0.01 mg/l (0.1%)
	Accuracy	1% ± 1 digit	1% ± 1 digit
	Calibration	1 point	1 point
	Temperature Compensation	0...50 °C	0...50 °C
	Salinity Compensation	0...40 ppt	0...40 ppt
	Air Pressure Compensation	600...1300 hPa	600...1300 hPa
OXYGEN	Range		0...10 mg/l
	Resolution		0.01 mg/l
	Accuracy		5% ± 1 digit
	Calibration		1 point
	pH Compensation		5...9 pH
FREE CHLORINE	Range		10...40 °C
	Resolution		
	Accuracy		
	Calibration		
	Temperature Compensation		
TEMPERATURE	Range	-30...+130 °C	-30...+130 °C
	Resolution	0.1 °C	0.1 °C
	Accuracy	0.3 °C	0.3 °C
	Calibration	1 point	1 point
AIR PRESSURE	Range	600...1300 hPa	600...1300 hPa
	Calibration	1 point	1 point
CHANNELS	Measurement	2	2
	Free Chlorine		1
	Temperature	2	1
INPUTS	Measurement	2 BNC	3 BNC
	Temperature	2 BNC, for Pt1000	1 BNC, for Pt1000
CALIBRATION	Reminder	0...999 h	0...999 h
	GLP	✓	✓
CONTROL	On/Off	✓	✓
	Proportional	✓	✓
	Wash program	✓	✓
	Alarm timer	✓	✓
DISPLAY	LCD	128x64 pixels	128x64 pixels
	White Backlight	✓	✓
	Real Time Clock	✓	✓
RECORDER	Two Outputs		4...20 mA, max. 300 Ω load
RS485	Baud Rate		300...19200 b/s
	Computer		via RS485/RS232 interface
DATA-LOGGING	Values		12000 + °C/date/time
	Modes		all
	Interval		1 s...4 h
RELAY OUTPUTS	Relays	4 solid state (*)	4 solid state (*)
	Voltage	12...250 VAC min. 1 mA/ max. 1 A	12...250 VAC min. 1 mA/ max. 1 A
SECURITY	Identification Number		✓
	Password Protection	✓	✓
POWER SUPPLY	Mains	210...250 VAC, 50/60 Hz	210...250 VAC, 50/60 Hz
DIMENSIONS	WxDxH	28x17x6 cm	28x17x6 cm
WEIGHT	Meter	1.3 kg	1.3 kg

DECLARATION OF CONFORMITY

We declare under our sole responsibility that the product

Multi-channel analysers

content of the type numbers

R361, R362,R363

to which this declaration relates is in conformity

with the following standards

EN61010

LOW VOLTAGE DIRECTIVE 73/23/EEG

EN50081-1

EN50082-1

EN60555-2

EMC DIRECTIVE 89/336/EEG

Turnhout, September 18, 2008

on behalf of CONSORT nv

WARRANTY CERTIFICATE

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory.

Accessories and breakable items such as electrodes are not warranted unless proven to be defective before shipment.

The original purchase order numbers, Consort invoice numbers and serial numbers of the products must be provided.

CONSORT will repair all defective equipment returned to it during the warranty period without charge (CIF Turnhout prepaid by sender), provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures described in the instruction manual and when not having been subject to accident, alteration, misuse or abuse.

If the products have been used with or have come into contact with fluids, an MSDS (material safety data sheet) must be supplied prior to issuing a return authorisation.

A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

Consort

www.consort.be

Parklaan 36 • B-2300 Turnhout • Belgium

Tel: (+32) (0)14 41 12 79

Fax: (+32) (0)14 42 91 79

Sales: sales@consort.be

Support: support@consort.be

Information: info@consort.be