Consort

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

R305
R315

October 2006

Table of contents

Introduction	5
Manufacturer	5
Warranty	5
Servicing.....	5
Keyboard	5
Contrast	5
Connections.....	6
Sample point	7
Language	7
Date and Time.....	7
Recorder.....	7
Password	7
Channels.....	8
Data-logging	8
RS485	8
Important	9
Measurements	9
Measurements	10
Important	10
Maintenance of pH electrodes.....	11
Maintenance of metal electrodes	11
Maintenance of conductivity electrodes	11
Maintenance of oxygen electrodes.....	11
GLP.....	12
Calibration reminder	12

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This instrument is manufactured with the latest technology and needs no particular maintenance. **Consort** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 0° and 40°C. Never use the instrument in a room with high humidity (>95 %) or at very low temperatures (condensation water!).

Manufacturer

Consort nv

Parklaan 36
B2300 Turnhout
Belgium

Tel (+32)(14)41 12 79
Fax (+32)(14)42 91 79
E-mail: info@Consort.be

Warranty

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **Consort** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **Consort** before returning any product for warranty repair on a freight prepaid basis!

Consort is not liable for consequential damages arising out of the use or handling of its products.

Servicing

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:


- Manual
- Cables
- Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **Consort** dealer.


Keyboard

- MODE** = Selects all modes or escapes from error traps, calibration procedures, etc..by returning to the original mode.
- ↔↑↓** = Button for entering a value or for selecting a function.
- CAL** = Starts or proceeds a calibration or a function.
- STOP** = Stops the control without interrupting the measurements.
- HELP** = Built-in manual.

Contrast

You can adjust the display contrast here!

pH/mV/O₂ electrode

Both channels are provided with floating inputs. The maximum common mode voltage (potential between the (-) terminals of CHAN.1 and CHAN.2) should not exceed ±1 V. Connect the internal wire of the electrode cable to terminal (+) and the screening of the cable to terminal (-) of (CHAN1 or CHAN2). Ground loops can be avoided by immersing an earthing rod in the process liquid close to the electrodes. Connect it only to terminal (M)!

Temperature probe

A separate automatic temperature compensation is provided for each channel but a common Pt1000 (or Pt100) thermocompensator for both channels can also be used. Connect the internal wire of the electrode cable to terminal (+) and the screening of the cable to terminal (-) of (C1 or C2).

Recorder

Two recorders, one for each channel, can be connected to the (+) and (-) terminals of (OUT1) and (OUT2). Max. impedance = 300 Ω.

Relays

There are three (R315) or four (R305) relays (REL1...REL4) with voltage-free contacts. Connect their two contacts in series with the load. The ratings of the connected load should lie between 12...250 V~, maximum 0.5 A.

Mains


Connect phase to terminal (P), neutral to terminal (N), and earth to earthing terminal.

- Always short circuit the unused pH/mV/O₂ input terminals!
- Special care should be taken when connecting a typical pH/mV/O₂ coaxial cable. Make sure that all parts are dry and clean and do not forget to remove the black antistatic foil between the inner lead and the outer screening!

RS485

Only use twisted pair UTP cable to connect up to 31 controllers to a computer!

Before the exact location of the electrodes is finalised, it is worthwhile to carry out a measuring profile of the process stream to ensure that the required homogeneity is being achieved. As a general rule the electrodes should protrude into the process stream by one-third of the pipe width and be positioned so as to avoid gas pockets or sediments.

Generally the measuring electrodes are inserted into the flow with the sample pick-up facing downstream in order to reduce the risk of blockage by suspended material.

All sensors are usually positioned vertically at a suitable point downstream of the dosing point, where the measurements are stable and reliable. The response of the dosing system will be affected by the time required to disperse the reagents and their arrival at the measuring point.

Language

1. Press **SET**.
2. Select [**LANGUAGE**] and press **CAL**.
3. Select the desired language and press **CAL**.

Date and Time

1. Press **SET**.
2. Select [**DATE/TIME**] and press **CAL**.
3. Adjust the date (respectively year, month and day) by pressing **CAL** each time.
4. Adjust the time (respectively hour, minutes and seconds) by pressing **CAL** each time.

Recorder

1. Press **SET**.
2. Select [**4-20 mA**] and press **CAL**.
3. Select the desired range and press **CAL**.
4. Select the desired minimum level corresponding to a 4 mA recorder output and press **CAL**.
5. Select the desired maximum level corresponding to a 20 mA recorder output and press **CAL**.

Adjustment:

1. Select [**ADJUST**] and press **CAL**.
2. Connect a mA-meter to the 4-20 mA outputs.
3. Select **OUT1 / 4 mA**, press **CAL** and read the mA-meter.
4. Adjust to the correct current and press **CAL**.
5. Repeat point 3-4 for all other currents.

Password

1. Press **SET**.
2. Select [**PASSWORD**] and press **CAL**.
3. A private code can be programmed to avoid undesired access to the instrument. Enter your secret sequence of 5 keys and press **CAL**.

1. Press **SET**.
2. Select [CHANNELS] and press **CAL**.
3. Select [INPUT] and press **CAL**.
4. Select the desired range and press **CAL**.
5. Select [ATC] and press **CAL**.
6. Select the desired temperature compensation input and press **CAL**.
7. Select [CONTROL] and press **CAL**.
8. Select the control type (on/off or proportional).
9. Pre-set the desired low level (Lo).
10. Select the corresponding relay (REL).
11. Pre-set the desired high level (Hi).
12. Select the corresponding relay (REL).
13. Pre-set the desired hysteresis (Hy) between the on and off state of the relays.
14. Pre-set the desired delay (dt) after which the relays should be activated, when one of the pre-set levels is exceeded.

Proportional control only:

15. Pre-set the desired the proportional control area (dP), in which proportional regulation should start.
16. Pre-set the desired pulse-time (Pt) of the relays. The relays will pulse at a rate proportional to the regulation difference.

Optional alarm procedure:

17. Select [ALARM AFTER] and press **CAL**.
18. The display shows the alarm time after which the process control should be interrupted in case of a failure. Pre-set the desired alarm-time and press **CAL**.
19. Select [REL] and press **CAL**.
20. Select the desired relay to be closed in case of an alarm and press **CAL**.

Data-logging

1. Press **SET**.
2. Select [DATA LOGGER] and press **CAL**.
3. Select [INTERVAL] and press **CAL**.
4. Select the desired time interval between the data-logging and press **CAL**.

View the stored values:

1. Select [TABLE], [GRAPH] or [SEND] and press **CAL**.
2. Select the desired form to display the stored data or send them to a computer and press **CAL** to continue. The scales of a graph can be modified by pressing **CAL**. Follow the instructions on the screen.

Erase the stored values:

1. Select [ERASE] and press **CAL**. Follow the instructions on the screen.

RS485

1. Press **SET**.
2. Select [RS485] and press **CAL**.
3. Select [BAUDRATE] and press **CAL**.
4. Select the desired baudrate and press **CAL**.
5. Select [IDENTIF.NUMBER] and press **CAL**.
6. Select the desired identification number (0...999) and press **CAL**. The instrument can be identified by e.g. a computer when a specific number is allocated to it.

- A blinking decimal point warns you for unstable measurements. Wait to read the display!
- Stirring the solution during the measurements promotes the homogeneity and is obligatory!
- The instrument will refuse automatic calibration when the electrode is unstable. Insufficient stirring or a worn electrode may be the cause.

Measurements

1. Should you want to recalibrate, select the desired measurement and press **CAL**. Follow the instructions on the screen. The displayed possibilities depend on the channel set-up and selected measurement.

(R305)

pH calibration:

2. Select the proper buffer from one of the values in memory or enter manually a special value and press **CAL**. The unused buffers should be switched off.
3. Remove the electrodes from their fittings, rinse them with distilled water and immerse them in the first buffer solution. Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
4. After rinsing the electrodes with distilled water, insert them back in their fittings and start the control.

mV calibration:

2. Remove the electrodes from their fittings, rinse them with distilled water, immerse them in a standard solution of known potential, adjust to the proper value and press **CAL**.
3. After rinsing the electrodes with distilled water, insert them back in their fittings and start the control.

Oxygen calibration:

2. The display shows the salinity correction. Select the proper values and press **CAL**. Leave salinity correction to zero unless you are going to measure in heavily salted solutions such as e.g. sea-water (35 g/l). Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
3. Remove the electrode from its fitting and rinse it with distilled water. The electrode exposed to the air reaches an equilibrium corresponding to the partial pressure of oxygen and thus to saturation in water at the given temperature. The instrument shows the measured saturation, current, temperature, and will calibrate automatically when readings are stable.
4. After rinsing the electrode with distilled water, insert it back in its fitting and start the control. Stirring the solution during the measurements promotes the homogeneity and is **obligatory!** The advection rate must be at least 10 cm/s.

Temperature calibration:

1. Without Pt1000, adjust the manual temperature compensation and proceed by pressing **SET**. Should you want to recalibrate, press **CAL**.
2. Immerse the Pt1000 in a standard solution of known temperature. Calibrate to the proper value and press **CAL**.

Adjustment:

2. Measure the real value of the processed liquid with an other electrode and meter. Adjust to the real value and press **CAL**. Press **STOP** here if you wish to return to the original readings.
- Press **↑** to maximum measurement.
 - Press **↓** to minimum measurement.
 - Press **STOP** to stop the control without interrupting the measurements.

- A blinking decimal point warns you for unstable measurements. Wait to read the display!
- Stirring the solution during the measurements promotes the homogeneity and is obligatory!
- The instrument will refuse automatic calibration when the electrode is unstable. Insufficient stirring or a worn electrode may be the cause.

Measurements

1. Should you want to recalibrate, select the desired measurement and press **CAL**. Follow the instructions on the screen. The displayed possibilities depend on the channel set-up and selected measurement. (R315)

Capacitive compensation:

2. Capacitive compensation increases accuracy in the very low conductivity ranges (<10 $\mu\text{S}/\text{cm}$). Select whether this compensation should be applied or not. Verify if the attached electrode is completely dry. Select [COMP.ON] or [COMP.OFF] and press **CAL**. Follow the instructions on the screen.

Conductivity calibration:

2. Select the proper standard from one of the values in memory or enter manually a special value and press **CAL**. Also select the temperature to which all future measurements will be referred to.
3. Remove the electrode from its fitting.
4. After rinsing the electrode several times with the standard solution, immerse it in that standard. The solution temperature is not so critical but should lie between 0°C and 30°C. Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
5. After rinsing the electrode with distilled water, insert it back in its fitting and start the control.

Temperature calibration:

2. Without Pt1000, adjust the manual temperature compensation and proceed by pressing **SET**. Should you want to recalibrate, press **CAL**.
3. Immerse the Pt1000 in a standard solution of known temperature. Calibrate to the proper value and press **CAL**.

Adjustment:

2. Measure the real value of the processed liquid with an other electrode and meter. Adjust to the real value and press **CAL**. Press **STOP** here if you wish to return to the original readings.
- Press **↑** to show the maximum measurement.
 - Press **↓** to show the minimum measurement.
 - Press **STOP** to stop the control without interrupting the measurements.

A pH electrode is active and stable only after wetting! For this purpose it must be immersed for **at least ten hours** in a 3...4 M KCl solution. During short interruptions (e.g. storage) the electrode should be immersed in a 3...4 M KCl solution. In doing this it is always kept ready for use. When the interruption is longer than a month, refill the closing cap with 3...4 M KCl and plug it on the electrode tip in order to protect the glass bulb. Before use, ensure that the reference part of the electrode is topped up with a 3...4 M KCl solution.

Avoid a low pressure inside the electrode! Therefore always remove the closure from the refilling aperture during the measurements as well as during the calibration. This allows the saltbridge solution to flow through the ceramic liquid junction and prevents contamination of the electrolyte. For the same reason, the inside level should always be higher than the outside level of the measuring solution. Close the refilling aperture again when storing the electrode.

A polluted electrode may be cleaned with a soft detergent or 0.1 M HCl. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

If the electrode is polluted by proteinaceous materials (such as blood), it should stand in a cleaning solution overnight and then be cleaned with distilled water before use. The pH electrode wears away by being used. If the electrode tends to respond slower and calibration becomes difficult, even after cleaning, it should be replaced by a new one.

Maintenance of metal electrodes

Metal electrodes (Pt, Ag, Au): Metal electrodes are always ready for use. During short interruptions they are immersed in distilled water. **They should be cleaned regularly:**

- Silver electrodes are immersed in a concentrated ammonia solution during one hour.
- Platinum or gold electrodes are immersed in concentrated nitric acid during one hour.

Maintenance of conductivity electrodes

A conductivity cell is active and stable only after wetting! For this purpose it must be immersed for **at least one hour** in distilled water. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

A polluted cell may be cleaned with a soft detergent or diluted nitric acid. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

Maintenance of oxygen electrodes

A dissolved oxygen electrode is active and stable only when polarised! The electrode body has almost unlimited life and requires no maintenance. After prolonged use of the electrode, it may become deactivated. An indication is that the electrode no longer responds correctly to calibration. In this case:

1. Carefully screw off the measuring head.
 2. Cautiously remove the precipitated deposits with the aid of filter paper. **Do not use grinding paper or a glass fibre!**
 3. Rinse several times with distilled water and shake carefully off the water drops.
 4. Take the measuring head, fill it with new electrolyte solution and very slowly screw it onto the electrode while holding it vertically.
 5. **Let the electrode rest for a few hours!** The electrode is now ready for use again.
- *Replace the membrane only when damaged, not when calibration is no longer possible!*

1. Press **SET**.
2. Select [GLP] and press **CAL**.
3. Select [SHOW REPORT] and press **CAL**. Browse with **←→↑↓** to show a complete calibration report.
4. Select [SEND] and press **CAL** to send the report to a computer.

pH/mV/O₂/°C-controller R305

```
-----
Date           : 15/09/1998
Time          : 11:32:04
Version       : 1.0
```

SETTINGS

```
-----
Identification No : 003
Password         : OFF
GLP reminder     (h) : ---

Channel1        : pH
ATC             : °C1
Manual temp.    (°C) : 23.4
```

```
Channel2        : mV
```

pH1

CALIBRATION

```
-----
Date           : 15/09/1998
Time          : 11:20:06
```

Buffer 4.01/6.87

```
Slope          (%) : 98.4
Zero point     (pH) : 6.87
```

AVERAGE VALUES

```
-----
Response time  (s) : 11
Temperature    (°C) : 22.6
```

STATISTICS

```
-----
Slope          (%) : +0.4
Zero point     (pH) : +0.02
Response time  (s) : -1
```

mV2

```
Date           : 15/09/1998
Time          : 11:20:06
Calibr. at     (mV) : 235
```

Calibration reminder

1. Select [INTERVAL] and press **CAL**.
2. Select the desired time interval between each automatic warning for a new calibration of the electrodes and press **CAL**.

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Specifications		R305
pH	Range	0...14 pH
	Resolution	0.01 pH
	Accuracy	0.1% ± 1 digit
	Calibration	1...2 points
	Buffers	9 + 2 user specified
	Temperature Compensation	-30...+130 °C
	ISO-pH	6...8 pH
	Slope	80...120%
mV	Range	±2000 mV
	Resolution	1 mV
	Accuracy	0.1% ± 1 digit
	Calibration	1 point
DISS. OXYGEN	Range	0...60 mg/L (0...600%)
	Resolution	0.05 mg/L (0.5%)
	Accuracy	1% ± 1 digit
	Calibration	1 point
	Temperature Compensation	0...50 °C
TEMPERATURE	Range	-30...+130 °C
	Resolution	0.1 °C
	Accuracy	0.3 °C
	Calibration	±10 °C
CHANNELS	pH/mV/Dissolved Oxygen	2
	Temperature	2 for Pt1000 (or Pt100)
INPUTS		screw terminals
CALIBRATION	Reminder	0...999 h
	GLP	✓
CONTROL	On/Off	✓
	Proportional	✓
DISPLAY	Set-up/Calibration	128x64 pixels
	Real Time Clock	✓
	Languages	EN, NL, FR, DE
RECORDER	Built-in	1200 points
	Interval	1...9999 s
	Two Outputs	4...20 mA, max. 300 Ω load
RS485	Baud Rate	300...9600 b/s
	Computer	via RS485/RS232 interface
DATA-LOGGING	Values (computer)	1200 + °C/date/time
	Interval	1...9999 s
RELAY OUTPUTS	Relays	4 solid state
	Mains	12...250 VAC/ max. 0.5 A
SECURITY	Identification Number	✓
	Password Protection	✓
POWER SUPPLY	Mains	210...250 VAC
DIMENSIONS	WxDxH	20x9x12 cm
WEIGHT	Meter	1.2 kg

Specifications		R315
CONDUCTIVITY	Range	0...2000 mS/cm
	Resolution	0.01 µS/cm
	Accuracy	0.5% f.s. of range
	Calibration	1 point
	Standards	0.01/0.1/1 M KCl
	Cell Constant	0.1/1/10 cm ⁻¹ ± 50%
	Temperature Compensation	-30...+130 °C
	Reference Temperature	20° or 25 °C
	Temperature Coefficient	natural waters (EN27888)
	Capacitive Compensation	✓
RESISTIVITY	Range	0...20 MΩ.cm
	Resolution	1 Ω.cm
TDS	Range	0...100 g/l
	Resolution	0.01 mg/l
TEMPERATURE	Range	-30...+130 °C
	Resolution	0.1 °C
	Accuracy	0.3 °C
	Calibration	±10 °C
CHANNELS	pH/mV/Dissolved Oxygen	2
	Temperature	2 for Pt1000 (or Pt100)
INPUTS		screw terminals
CALIBRATION	Reminder	0...999 h
	GLP	✓
CONTROL	On/Off	✓
	Proportional	✓
DISPLAY	Set-up/Calibration	128x64 pixels
	Real Time Clock	✓
	Languages	EN, NL, FR, DE
RECORDER	Built-in	1200 points
	Interval	1...9999 s
	Two Outputs	4...20 mA, max. 300 Ω load
RS485	Baud Rate	300...9600 b/s
	Computer	via RS485/RS232 interface
DATA-LOGGING	Values (computer)	1200 + °C/date/time
	Interval	1...9999 s
RELAY OUTPUTS	Relays	3 solid state
	Mains	12...250 VAC/ max. 0.5 A
SECURITY	Identification Number	✓
	Password Protection	✓
POWER SUPPLY	Mains	210...250 VAC
DIMENSIONS	WxDxH	20x9x12 cm
WEIGHT	Meter	1.2 kg

WARRANTY CERTIFICATE

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory.

CONSORT will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse.

A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

DECLARATION OF CONFORMITY

We declare under our sole responsibility that the product

**Multi-channel analysers
content of the type numbers
R305, R315**

to which this declaration relates is in conformity
with the following standards

**EN61010
LOW VOLTAGE DIRECTIVE 73/23/EEG
EN50081-1
EN50082-1
EN60555-2
EMC DIRECTIVE 89/336/EEG**

Turnhout, September 8, 1998

A handwritten signature in black ink, appearing to be 'J. G. ...', is written over a horizontal line.

on behalf of **CONSORT** nv

CERTIFICATE OF REGISTRATION

This is to certify that:

Consort NV
Turnhout

Has been assessed by QualityMasters in respect of its Quality Management System and found to comply with:

NBN EN ISO 9001:2000

Approval is hereby granted for registration providing the rules and conditions relating to certification are observed at all times.

Scope: Production of high-end equipment for electrochemistry and electrophoresis.

Original approval: 2005-08-29
Current certificate: 2005-09-17

Certificate expiry: 2010-08-29
Certificate number: BE-0211-317

Authorized signatory

For QualityMasters Belgium

This certificate remains the property of QualityMasters.