Consort

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

R735

October 2006

Table of contents

Manufacturer	3
Warranty	3
Introduction	3
Servicing.....	3
Keyboard	3
Error codes.....	4
Inputs	4
Recorder output.....	4
Digital output.....	4
Date	5
Time.....	5
Temperature	5
Recorder.....	5
Channels.....	5
Relais	6
Alarm	6
Specific features.....	6
mV measurement	6
pH measurement.....	7
Conductivity measurement.....	7
Maintenance of pH electrodes.....	8
Maintenance of metal electrodes	8
Maintenance of conductivity electrodes	8
Data-logging	9

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This instrument is manufactured with the latest technology and needs no particular maintenance. **Consort** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 0° and 40°C. Never use the instrument in a room with high humidity (>95 %) or at very low temperatures (condensation water!).

Connect the instrument only to an earthed power line. The required power source is indicated on the label at the back of the instrument. Do not cut and splice the power cord. When removing the power cord from the wall outlet, be sure to unplug by holding the plug attachment and not by pulling the cord. Do not hold the plug by wet hand.

Manufacturer

Consort nv

Parklaan 36
B2300 Turnhout
Belgium

Tel (+32)(14)41 12 79
Fax (+32)(14)42 91 79
E-mail: info@Consort.be

Warranty

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **Consort** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **Consort** before returning any product for warranty repair on a freight prepaid basis!

Consort is not liable for consequential damages arising out of the use or handling of its products.

Servicing

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

- Manual
- Cables
- Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **Consort** dealer.

Keyboard

- MODE** = Selects all modes or escapes from error traps, calibration procedures, etc..by returning to the original mode.
- CAL** = Starts or proceeds a calibration or a function.
- ▲/▼** = Button for entering a value or for selecting a function.
- SET** = Installs the meter according your application.
- HOLD** = Holds display when measuring.
- RES** = Changes resolution/units of the display.
- PRINT** = Displayed value is printed through the RS232 output.
- ON/OFF** = Switches the instrument on or off.

[Or. pH]	=	overrange pH
[Or. mV]	=	overrange mV
[Or. μ S]	=	overrange μ S
[Or. MS]	=	overrange mS
[Or. OhM]	=	overrange Ω
[Or. °C]	=	Overrange °C
[Or.PPt]	=	overrange ppt (TDS)
[Or.SAL]	=	overrange salinity
[no dAtA]	=	data-logging memory empty
[not.CAL]	=	calibration procedure interrupted
[Err.SLP]	=	slope not within 80...120 %
[Err.pHo]	=	ISO pH not within 6...8 pH
[Err.cc]	=	cell constant not within 70...130% of the nominal value
[Err.CAL]	=	calibration error
[Err.COd]	=	invalid code
[Err.COM]	=	general communication error
[Err.MEM]	=	fatal memory error; have your instrument serviced!
[ILLEGAL]	=	an unauthorised person has failed to repair the unit; warranty is elapsed from now on!

Inputs

The pH electrode should be connected to the coaxial pH connector. The mV electrode should be connected to the coaxial mV connector. If separate electrodes are used, connect the reference electrodes to the REF. terminals. Connect a conductivity cell to the coaxial EC input. Automatic temperature compensation and temperature measurements are possible by plugging a Pt1000 thermocompensator into the °C terminals. Without compensator, the manual temperature compensation is automatically switched on.

Recorder output

A recorder can be connected to the red (+) and black (-) REC terminals. Use only laboratory recorders with a high input impedance!

Digital output

A standard RS232 output terminal (DP9) is provided for interfacing the instrument with a printer or computer (using the optional AK2310 cable). Data is sent in the ASCII code at a BAUD rate of 150...4800 bps (8 bit, no parity, 1 stopbit).

Serial port pinout specifications:

pin 1 : connected to pin 4 and pin 6
pin 2 : TxD, transmit data
pin 3 : RxD, receive data
pin 4 : connected to pin 1 and pin 6
pin 5 : Gnd, signal ground
pin 6 : connected to pin 1 and pin 4
pin 7 : connected to pin 8
pin 8 : connected to pin 7
pin 9 : not connected

1. Press **SET** until the display shows [SEt dAtE] and press **CAL**.
2. The display shows [dAtE=on] or [dAtE=OFF]. Select the desired mode and press **CAL**.
3. Program the date (respectively year, month and day) by pressing **CAL** each time.

1. Press **SET** until the display shows [SEt tiME] and press **CAL**.
2. The display shows [tiME=on] or [tiME=OFF]. Select the desired mode and press **CAL**.
3. Program the time (respectively hour, minutes and seconds) by pressing **CAL** each time.

1. Press **SET** until the display shows [SEt °C] and press **CAL**.
2. The display will show the measured temperature e.g. [°C=23.7][Aut]. If necessary, calibrate your probe and press **CAL**. Press **RES** to reset the temperature calibration to the original factory calibration.
3. When no Pt1000 temperature probe is connected to the instrument the display will show the manual temperature e.g. [°C=23.7][MAn]. Shift to the proper temperature and press **CAL**.

1. Press **SET** until the display shows [SEt rEc] and press **CAL**.
2. The display shows for all channels e.g. [rEc=PH], [rEc=MV], [rEc=EC], [rEc=°C]...etc. Select the desired mode and press **CAL**.
3. This step is **skipped** unless [rEc=EC] has been chosen. Press **CAL** and the display shows the conductivity measuring range and units, e.g. [rEc 0.01MS][Unit]. Select the desired range and press **CAL**.
4. The display shows e.g. [rEc=6.00][LoW]. Select the desired minimum level corresponding to a 0 V recorder output and press **CAL**.
5. The display shows e.g. [rEc=8.00][HiGH]. Select the desired maximum level corresponding to a 2 V recorder output and press **CAL**.

1. Press **SET** until the display shows [SEt CHAn] and press **CAL**.
2. The display shows e.g. [PH=on] or [PH=OFF]...etc. Select the desired mode (respectively pH, mV and EC) by pressing **CAL** each time.
3. The display shows e.g. [rEF=25°C]. Select the temperature to which all future conductivity measurements will be referred and press **CAL**.

1. Press **SET** until the display shows [SEt rEL] and press **CAL**.
2. The display shows e.g. [r1=PH.Hi], [r1=PH.Lo], [r1=ALrM]...etc. Select the desired mode for relay-1 and press **CAL** each time.
3. The display shows the delay (s) after which the relays should be activated, when one of the pre-set levels is exceeded, e.g. [r1.dt=120]. Pre-set the desired delay for the relays and press **CAL**.

1. Press **SET** until the display shows [SEt ALArM] and press **CAL**.
2. The display shows e.g. [ALrM PH], [ALrM MV], [ALrM EC]...etc. Select the desired mode and press **CAL**.
3. The display shows [AL.= on] or [AL.= OFF]. Select the desired mode and press **CAL**.
4. The display shows e.g. [AL.=6.00][LoW]. Select the desired minimum level and press **CAL**.
5. The display shows e.g. [AL.=8.00][HiGH]. Select the desired maximum level and press **CAL**.
6. The display shows the hysteresis between the ON and OFF state of the relays e.g. [AL.=0.10][HySt]. Pre-set the desired hysteresis and press **CAL**.

1. Press **SET** until the display shows [SEtUP] and press **CAL**.
2. The display lets you choose between [bEEP on] or [bEEP OFF]. Decide whether or not the beeper should sound each time the keyboard is touched. Press **CAL**.
3. The display shows [ScAn OFF] or [ScAn on]. Decide whether the display should alternate between the selected modes with an interval of 4 s or not and press **CAL**. Pressing **PRINT** in this mode will print all values of the selected channels.
4. The display shows e.g. [br=2400]. Pre-set the desired baud-rate (150...4800 b/s) and press **CAL**.
5. The display shows e.g. [rS=120]. Select the desired time interval (0...9999 s) between the transmitted data (RS232). Pre-set to zero if no automatic printing is required. Press **CAL**.
6. The display shows e.g. [Id.no=5]. Select the desired identification number (0...999) for your meter and press **CAL**. The instrument can be identified by e.g. a computer when a specific number is allocated to it. Pre-set to zero to ignore this number.
7. The display shows [Code on], [Code OFF] or [Code CAL]. Decide whether all future access to the instrument should be by entering a code (on) or not (off). Press **CAL**. Select [Code CAL] to enter your personal code to avoid undesired access to the instrument. The display shows [E=-----] while [=] is blinking. Enter your secret sequence of 5 keys to proceed. Press **CAL**.

1. Select the mV range by pressing **MODE**. The display will immediately show the measured mV according to the previous calibration. Should you want to recalibrate, press **CAL**.
2. Immerse the electrodes in a standard solution of known potential.
3. The display shows e.g. [MV=476.3]. Adjust to the proper value and press **CAL** again or press **RES** to reset the absolute readings.
4. After rinsing the electrodes with distilled water, immerse them in the solution to be measured and read the display.
5. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.

1. Select the pH range by pressing **MODE**. The display will immediately show the measured pH according to the previous calibration. Should you want to recalibrate, press **CAL**.
2. Rinse the electrodes with distilled water and immerse them in the first buffer solution.
3. The display shows one of the 9 buffers in memory (b1) or the manual input (M1) of user specified buffers, e.g. [b1=4.01] while [b1] or [M1] is blinking. Select the proper buffer and press **CAL**. If you wish to change the manual buffer, first press **RES** to reset to the proper value and then press **CAL**.
4. The instrument shows the measured buffer and will calibrate automatically when readings are stable ([=] stops blinking). If the buffer has a different temperature, compensate the indicated value manually (not necessary when using a Pt1000).
5. Go on in the same way with all the next buffer (b2) and always press **CAL**. Press **MODE** to start the measurements.
6. After rinsing the electrodes with distilled water, immerse them in the samples and read the pH-value from the display.
7. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.

Conductivity measurement

1. Select the conductivity range by pressing **MODE** (the display shows [μ S] or [MS]) and press **CAL**.
 2. The display will show the cell constant (e.g. [cc=1.073]). Press **CAL** to start the automatic determination procedure.
 3. The display shows one of the KCl standard solutions in memory [St=0.01M][cc] (1413 μ S/cm), [St=0.1M][cc] (12.88 mS/cm), [St=1M][cc] (111.8 mS/cm) or the manual input [St=MAn]. Select the proper standard.
 4. This step is **skipped** unless the manual input has been chosen. The display shows e.g. [MS=10.00][cc] while the units are blinking. Choose between 1.000 (units), 10.00 (tens) or 100.0 (hundreds) the desired calibration range and press **CAL** again. The display shows e.g. [MS=1.234][cc] while [=] is blinking. Select the desired standard value.
 5. After rinsing the cell several times with the standard solution, immerse it the standard solution. The solution temperature is not so critical but it should lie between 0°C and 30°C. If the standard has a different temperature, compensate the indicated value manually (not necessary when using a Pt1000). Press **CAL**.
 6. The display will show the computed cell constant e.g. [cc=1.045] and will calibrate automatically when readings are stable ([=] stops blinking). Press **MODE** to start the measurements.
 7. Rinse the electrodes with the sample and immerse them in it.
 8. Read the display (remember: compensated to the reference temperature!). The autoranging feature gives you the largest possible number of digits + their unit.
 9. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).
- Press **RES** to switch between ppm (TDS: ppm = mg/l, ppt = g/l), SAL (salinity), S/cm (conductivity) or Ω .cm (resistivity: E6 = M Ω , E3 = k Ω , E0 = Ω).

A pH electrode is active and stable only after wetting! For this purpose it must be immersed for **at least ten hours** in a 3...4 M KCl solution. During short interruptions (e.g. storage) the electrode should be immersed in a 3...4 M KCl solution. In doing this it is always kept ready for use. When the interruption is longer than a month, refill the closing cap with 3...4 M KCl and plug it on the electrode tip in order to protect the glass bulb. Before use, ensure that the reference part of the electrode is topped up with a 3...4 M KCl solution.

Avoid a low pressure inside the electrode! Therefore always remove the closure from the refilling aperture during the measurements as well as during the calibration. This allows the saltbridge solution to flow through the ceramic liquid junction and prevents contamination of the electrolyte. For the same reason, the inside level should always be higher than the outside level of the measuring solution. Close the refilling aperture again when storing the electrode.

A polluted electrode may be cleaned with a soft detergent or 0.1 M HCl. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

If the electrode is polluted by proteinaceous materials (such as blood), it should stand in a cleaning solution overnight and then be cleaned with distilled water before use. The pH electrode wears away by being used. If the electrode tends to respond slower and calibration becomes difficult, even after cleaning, it should be replaced by a new one.

Metal electrodes (Pt, Ag, Au): Metal electrodes are always ready for use. During short interruptions they are immersed in distilled water. **They should be cleaned regularly:**

- Silver electrodes are immersed in a concentrated ammonia solution during one hour.
- Platinum or gold electrodes are immersed in concentrated nitric acid during one hour.

A conductivity cell is active and stable only after wetting! For this purpose it must be immersed for **at least one hour** in distilled water. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

A polluted cell may be cleaned with a soft detergent or diluted nitric acid. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

Pre-set and start the data-logging:

1. Press **SET** until the display shows [SEt LOG] and press **CAL**.
2. Select [ConFiG][LOG] and press **CAL**.
3. Choose between [tiMEr], [MANuAL] or [OFF]. Select the desired mode and press **CAL**.
4. This step is **skipped** unless automatic data-logging is chosen (timer). The display shows e.g. [dt=100]. Select the desired time interval (0...9999 s) between the automatic data-logging and press **CAL**.
5. When manual data-logging has been selected, press **PRINT** to put a next measurement into memory. In the meantime the display shows the storage number e.g. [LOG 7].

View the stored values on the display:

1. Press **SET** until the display shows [SEt LOG] and press **CAL**.
2. Select [LSt dAtA] and press **CAL**.
3. The display shows the first data-logged value in memory (a blinking [n] appears). Press \blacktriangle or \blacktriangledown to view all other values.
4. While pressing **HOLD** the corresponding date is shown.
5. While pressing **RES** the corresponding time is shown.
6. While pressing **CAL** the corresponding number is shown.
7. Press **PRINT** to send the recalled data to a printer.
8. Press **SET** to return.

Erase the the stored values:

1. Press **SET** until the display shows [SEt LOG] and press **CAL**.
2. Select [CLr LOG] and press **CAL**.

Send the stored values to a printer or computer:

1. Press **SET** until the display shows [SEt LOG] and press **CAL**.
2. Select [Prt dAtA] and press **CAL** to print the data.
3. The display shows [PrintinG] while all data-logged values, including corresponding dates and times, are sent to the RS232 output in the following format, compatible with EXCEL™, LOTUS 123™,....The columns are separated by a TAB character:

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Specifications		R735
pH	Range	0...14 pH
	Resolution	0.01 pH
	Accuracy	0.1% ± 1 digit
	Calibration	1...2 points
	Buffers	9 pre-programmed 5 user specified
	Temperature Compensation	0...100°C
	ISO-pH	6...8 pH
Slope		80...120%
mV	Range	±2000 mV
	Resolution	1 mV
	Accuracy	0.1% ± 1 digit
	Calibration	1 point
CONDUCTIVITY	Range	0...200 mS/cm
	Resolution	0.01 µS/cm
	Accuracy	0.5% f.s. of range
	Calibration	1 point
	Standards	0.01/0.1/1 M KCl 3 user specified
	Cell Constant	0.5...1.5 cm ⁻¹
	Temperature Compensation	0...100°C
	Reference Temperature	20° or 25°C
	Temperature Coefficient	natural waters (EN27888)
RESISTIVITY	Range	0...20 MΩ.cm
	Resolution	10 Ω.cm
SALINITY	Range	0...70 ppt
	Resolution	0.1 ppt
TDS	Range	0...100 g/l
	Resolution	0.01 mg/l
TEMPERATURE	Range	0...100°C
	Resolution	0.1°C
	Accuracy	0.3°C
	Calibration	±10°C
CHANNELS	pH	1
	mV	1
	Conductivity	1
	Temperature	1
INPUTS	pH	1 BNC, 10 ¹² Ω
	mV	1 BNC, 10 ¹² Ω
	Conductivity	1 BNC
	Temperature	2 banana, for Pt1000
CALIBRATION	GLP	✓
DISPLAY	Readings	8 digits LED
	Temperature	4 digits LED
	Alternating	✓
	Hold Function	✓
	Selectable Resolution	✓
	Real Time Clock	✓
RECORDER	Output	0...2 V programmable
RS232	Baud Rate	300...19200 b/s
	Interval	1...9999 s
DATA-LOGGING	Values	500 + °C/date/time
	Modes	pH, mV, conductivity
	Manual or Timed	✓
	Interval	1...9999 s
ALARM	Audible	✓
	Relay	✓
RELAY OUTLETS	Relays	2 mechanical
	Mains	2 IEC, max. 5 A, fused
CONTROL	On/Off	✓
	Proportional	✓
SECURITY	Identification Number	✓
	Password Protection	✓
POWER SUPPLY	Mains	210...250 VAC
DIMENSIONS	WxDxH	30x18x12 cm
WEIGHT	Meter	2.4 kg

WARRANTY CERTIFICATE

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory.

CONSORT will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse.

A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

DECLARATION OF CONFORMITY

We declare under our sole responsibility that the product

**Multi-parameter analyser
content of the type numbers
R735**

to which this declaration relates is in conformity
with the following standards

**EN61010
LOW VOLTAGE DIRECTIVE 73/23/EEG
EN50081-1
EN50082-1
EN60555-2
EMC DIRECTIVE 89/336/EEG**

Turnhout, May 25, 1999

A handwritten signature in black ink, appearing to be 'J. G. G.', is written over a horizontal line.

on behalf of **CONSORT** nv

CERTIFICATE OF REGISTRATION

This is to certify that:

Consort NV
Turnhout

Has been assessed by QualityMasters in respect of its Quality Management System and found to comply with:

NBN EN ISO 9001:2000

Approval is hereby granted for registration providing the rules and conditions relating to certification are observed at all times.

Scope: Production of high-end equipment for electrochemistry and electrophoresis.

Original approval: 2005-08-29

Certificate expiry: 2010-08-29

Current certificate: 2005-09-17

Certificate number: BE-0211-317

Authorized signatory

For QualityMasters Belgium

This certificate remains the property of QualityMasters.