CONSORT

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

T851
T852

August 2004

QUALITYMASTERS

CERTIFICATE OF REGISTRATION

This is to certify that:

Consort NV

Parklaan 36
B-2300 Turnhout

Has been assessed by QualityMasters in respect of its Quality Management System and found to comply with:

NEN-EN-ISO 9001 : 2000

Approval is hereby granted for registration providing the rules and conditions relating to certification are observed at all times.

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Lack of fulfillment of conditions as set forth in the certification regulations may render this certificate invalid. The use of the accreditation mark indicates accreditation in respect to the activities covered by the accreditation number EU9706010F.

Q

M

REGISTRATION

CERTIFICATE OF

Authorized signatory

A. LEUWATY
QualityMasters B.V.

This certificate remains the property of QualityMasters B.V.

Instruments manufactured by CONSORT

Electrophoresis power supplies

•

Ion analysers

•

pH meters

•

Conductometers

•

Dissolved oxygen meters

•

Thermometers

•

Industrial controllers

CONSORT nv
Parklaan 36
B2300 Turnhout
Belgium

Find out more at www.consort.be

Tel (+32)(14)41 12 79

Fax (+32)(14)42 91 79

E-mail: info@consort.be

DECLARATION OF CONFORMITY

We declare under our sole responsibility that the product

Thermometer

content of the type numbers

T851, T852

to which this declaration relates is in conformity

with the following standards

EN61010

LOW VOLTAGE DIRECTIVE 73/23/EEG

EN50081-1

EN50082-1

EN60555-2

EMC DIRECTIVE 89/336/EEG

Turnhout, December 11, 1996

A handwritten signature in black ink, appearing to be 'J. G. G.', is written over a horizontal line.

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This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 4° and 40° C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

Connect the jack of a suitable AC adapter (model **A4072**) to the DC socket. Do not hold the adapter by wet hand.

Manufacturer

CONSORT nv

Parklaan 36
B2300 Turnhout
Belgium

Tel (+32)(14)41 12 79
Fax (+32)(14)42 91 79
E-mail: info@consort.be

Warranty

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

CONSORT is not liable for consequential damages arising out of the use or handling of its products.

Servicing

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Manual

Cables

Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

Ranges	J: -200 ... +900°C (-328 ...+1652°F) K: -200 ... +1370°C (-328 ...+2498°F) T: -200 ... +600°C (-328 ...+1112°F) E: 0 ... +900°C (+32 ...+1652°F) Pt100: -200... +850°C (-328 ...+1562°F)
Resolution	<1000°C:0.1 °C (0.1 °F) >999.9°C: 1°C (1°F)
Inputs	4 x terminal block
Digital output	programmable RS232, 150...4800 b/s, for bi-directional communication with a computer or printer
Analogue output	0...2 V, ca 8 kΩ, programmable 8 bit D/A with adjustable span
Display	8+4 digit/character, 13 mm LED display, shows readings, °C and messages
Keys	9 tactile membrane keys
Ambient temp.	4...40°C
Rel. humidity	0...90 % (non-condensing!)
Power supply	210-250 V~, 50/60 Hz, max. 4 VA
Cabinet	IP65 cabinet
Dimensions	260 x 180 x 90 mm
Weight	1 kg

Keyboard

MODE	= Selects all modes or escapes from error traps, calibration procedures, etc..by returning to the original mode.
CAL	= Starts or proceeds a calibration or a function.
▲/▼	= Button for entering a value or for selecting a function.
SET	= Installs the meter according your application.
HOLD	= Holds display when measuring.
RES	= Changes resolution/units of the display.
PRINT	= Displayed value is printed through the RS232 output.
ON/OFF	= Switches the instrument on or off.

Error codes


[no dAtA]	= Data-logging memory empty.
[Err.COd]	= Invalid code.
[Err.CAL]	= Calibration error.
[not.CAL]	= Calibration procedure interrupted.
[Err.COM]	= General communication error.
[Err.MEM]	= Fatal memory error. Have your instrument serviced!.
[ILLEGAL]	= An unauthorised person has failed to repair the unit. Warranty is elapsed from now on!

T851:

Connect the temperature probes with the terminal blocks and insert them in the instrument. The channel numbers start from the left to the right (backview). The even numbers are the negative poles. The uneven numbers are the positive poles. **Important:** connect the earth terminal to the earthing to avoid ground loop errors!

T852:

Connect the temperature probes with the terminal blocks and insert them in the instrument. The channel numbers start from the left to the right (backview).

**Recorder output**

A recorder can be connected to the red (+) and black (-) terminals. Use only laboratory recorders with a high input impedance!

Digital output

A standard RS232 output terminal (DP9) is provided for interfacing the instrument with a printer or computer. Data is sent in the ASCII code at a BAUD rate of 150...4800 bps (8 bit, no parity, 1 stopbit).

Serial port pinout specifications:

- pin 1 : connected to pin 4 and pin 6
- pin 2 : TxD, transmit data
- pin 3 : RxD, receive data
- pin 4 : connected to pin 1 and pin 6
- pin 5 : Gnd, signal ground
- pin 6 : connected to pin 1 and pin 4
- pin 7 : connected to pin 8
- pin 8 : connected to pin 7
- pin 9 : not connected

Channels

1. Press **SET** until the display shows [SEt CHAn] and press **CAL**.
2. The display shows [CH.01=on], [CH.01=OFF]. Select the desired mode for channel-1 and press **CAL**.
3. Go on in the same way with all next channels by repeating the procedure from point 2). Always press **CAL** to continue. Press **RES** to choose simultaneously the same mode for all channels.

1. Press **SET** until the display shows [SEt dAtE] and press **CAL**.
2. The display shows [dAtE=on] or [dAtE=OFF]. Select the desired mode and press **CAL**.
3. Program the date (respectively year, month and day) by pressing **CAL** each time.

1. Press **SET** until the display shows [SEt tiME] and press **CAL**.
2. The display shows [tiME=on] or [tiME=OFF]. Select the desired mode and press **CAL**.
3. Program the time (respectively hour, minutes and seconds) by pressing **CAL** each time.

1. Press **SET** until the display shows [SEt tyPE] and press **CAL**.
2. The display shows e.g. for channel-1 [Cr-Al][CH.01] (type-K), [Fe-Co][CH.01] (type-J), [Cu-Co][CH.01] (type-T) or [Ni-Co][CH.01] (type-E). Select the desired thermocouple for this channel and press **CAL**.
3. Go on in the same way with all next channels by repeating the procedure from point 2). Always press **CAL** to continue. Press **RES** to choose simultaneously the same type of thermocouple for all channels.

1. Press **SET** until the display shows [SEt rEC] and press **CAL**.
2. The display shows for all channels e.g. [rEc=CH.01], [rEc=CH.02], ...etc. Select the desired channel and press **CAL**.
3. The display shows e.g. [Lo=-200°][rEc] while [=] is blinking. Select the desired minimum level corresponding to a 0 V recorder output and press **CAL**.
4. The display shows e.g. [Hi=1370°][rEc] while [=] is blinking. Select the desired maximum level corresponding to a 2 V recorder output and press **CAL**.
5. Go on in the same way with all next channels by repeating the procedure from point 2). Always press **CAL** to continue.

1. Press **SET** until the display shows [SEtUP] and press **CAL**.
2. The display lets you choose between [bBEEP on] or [bBEEP OFF]. Decide whether or not the beeper should sound each time the keyboard is touched. Press **CAL**.
3. The display lets you choose between [POWER][on] or [POWER][OFF]. Choose [POWER][on] to inactivate the ON/OFF-key (meter is always on and cannot be switched off). After an interruption of the mains power, the meter will resume measurements automatically. Choose [POWER][OFF] to activate the ON/OFF-key (meter can be switched on or off). Press **CAL**.
4. The display shows e.g. [br=2400]. Pre-set the desired baud-rate (150...4800 b/s) and press **CAL**.
5. The display shows e.g. [rS=120]. Select the desired time interval (0...9999 s) between the transmitted data (RS232). Pre-set to zero if no automatic printing is required and press **CAL**.
6. The display shows e.g. [Id.no=5]. Select the desired identification number (0...999) for your meter and press **CAL**. The instrument can be identified by e.g. a computer when a specific number is allocated to it. Pre-set to zero to ignore this number.
7. The display shows [Code on], [Code OFF] or [Code CAL]. Decide whether all future access to the instrument should be by entering a code (on) or not (off). Press **CAL**. Select [Code CAL] to enter your personal code to avoid undesired access to the instrument. The display shows [E=-----] while [=] is blinking. Enter your secret sequence of 5 keys to proceed. Press **CAL**.

Alarm

1. Press **SET** until the display shows [SEt ALrM] and press **CAL**.
2. The display shows [bBEEp on] or [bBEEp OFF]. Select the desired mode and press **CAL**.
3. The display shows the delay (s) if one of the thresholds is passed, e.g.[dt=0] [ALrM]. Select the delay and press **CAL**.
4. The display shows e.g. [CH.01=on][ALrM] while [=] is blinking. Select the desired channel and press **CAL**.
5. The display shows e.g. [CH.01=on] or [CH.01=OFF] while [=] is blinking. Select the desired mode and press **CAL**.
6. The display shows e.g. [Lo=2°][ALrM] while [=] is blinking. Select the desired minimum level and press **CAL**. Each time the measurements are below this level, the beeper will sound an alarm.
7. The display shows e.g. [Hi=900°][ALrM] while [=] is blinking. Select the desired maximum level and press **CAL**. Each time the measurements are above this level, the beeper will sound an alarm.
8. Go on in the same way with all next channels by repeating the procedure from point 4). Always press **CAL**.

1. Press **SET** until the display shows [SEt Unit] and press **CAL**.
2. Choose between °C, °F or At (absolute temperature, K) and press **SET** to return.

Calibration

1. Verify that the instrument is switched on in a stable temperature environment for at least 1 hour.
2. Select the temperature range by pressing **MODE**. The display immediately shows the measured value according to the last calibration. Press **MODE** to see the values of the other channels. Press **CAL** to calibrate.
3. The display shows [ALL CHAn] (simultaneous calibration of all channels) or e.g. [CH 01] (calibration of the selected channel). Select the desired mode and press **CAL**.
4. The display shows [rEF=AUT] (channel-1 is the reference for all the other channels) or [rEF=MAAn] (manual calibration). Select the desired mode and press **CAL**.
5. The display shows e.g. [°C1=18.9][CAL] while [=] is blinking. Select the desired value (manual calibration only). Wait at least 10 seconds, after stable readings, and press **CAL**.
6. The display shows e.g. [°C2=99.7][CAL] while [=] is blinking. Select the desired value ($>^{\circ}\text{C}1+50$) (manual calibration only). Wait at 10 seconds, after stable readings, and press **CAL** to continue or **MODE** to return (1-point calibration).

Measurement

1. Verify that the instrument is switched on in a stable temperature environment for at least 1 hour.
 2. Select the temperature range by pressing **MODE**. The display immediately shows the measured value according to the last calibration. Press **MODE** to see the values of the other channels.
- Press **RES** to change the resolution from 0.1 to 1.
 - Press **▲** to display the maximum temperature ever measured by the selected channel (after switching on), e.g. [56.3⁻°C]. Press **▲** again to return to the normal measurements.
 - Press **▼** to display the minimum temperature ever measured by the selected channel (after switching on), e.g. [-4.5₋°C]. Press **▼** again to return to the normal measurements.
 - Any alarm will close the relay, sound an alert and display [ALrM]. Press any key to stop the beeper and to open the relay again. The display will show e.g. [AL.07] to indicate the channel in a alarm.

Alternating measurements

1. Press **SET** until the display shows [SCAn] and press **MODE** to start the alternating mode. The display alternates between the selected modes with an interval of 4 s. Pressing **PRINT** in this mode will print all values of the selected channels.
2. Press **MODE** again to stop this mode and continue with the actual channel.

1. Press **CAL** while measuring.
2. The display shows [ALL CHAN] or e.g. [CH 01]. Select the desired mode and press **PRINT** to print a complete calibration report of the selected channel, e.g.:

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Slope (d°C/°C) : +0.0007
Offset (°C) : -0.1

Pre-set and start the data-logging:

1. Press **SET** until the display shows [SEt LoG] and press **CAL**.
2. Select [LOG dAtA] and press **CAL**.
3. Choose between [tiMEr], [MAnuAL], [ALArM] or [OFF] and press **CAL**. Select [OFF] if no data-logging is needed.
4. Select [MAnuAL] for manual data-logging and press **MODE**. Press **PRINT** to put a next measurement into memory. In the mean time the display shows the logging-number e.g. [LOG 357].
5. Select [ALArM] (collects up to 1200 values with 2 different time intervals: alarm on or off) or [tiMEr] (collects up to 3600 values). Press **CAL**.
6. The display shows e.g. [dt=100]. Select the desired time interval (1...9999 s) between a normal data-logging and press **MODE**. A blinking [L] appears while measuring. Press **SET** to end. The display shows [End LOG].
7. The display shows e.g. [dAt=20]. Select the desired time interval (1...9999 s) between a data-logging, when the alarm is activated, and press **MODE**. A blinking [L] appears while measuring. Press **SET** to end. The display shows [End LOG].

Clear the memory:

1. Press **SET** until the display shows [SEt LoG] and press **CAL**.
2. Select [CLr dAtA] and press **CAL** to clear the memory.

View the stored values on the display:

1. Press **SET** until the display shows [SEt LoG] and press **CAL**.
2. Select [LSt dAtA] and press **CAL**.
3. The display shows the first data-logged value in memory (a blinking [n] appears). Press \blacktriangle or \blacktriangledown to view all other values.
4. While pressing **HOLD** the corresponding date is shown.
5. While pressing **RES** the corresponding time shown.
6. Press **PRINT** to send the recalled data to a printer.
7. Press **SET** to return.

Send the stored values to a printer or computer:

1. Press **SET** until the display shows [SEt LoG] and press **CAL**.
2. Select [Prt dAtA] and press **CAL** to print the data.
3. The display shows [PrintinG] while all data-logged values, including corresponding dates and times, are sent to the RS232 output in the following format, compatible with EXCEL™, LOTUS 123™,....The columns are separated by a TAB character:

Identification No: 037

D	H	CH02	CH03	CH04	CH05
29/01/95	13:19:57	25.6	200.4	45.8	587.6
29/01/95	13:20:57	25.4	200.7	45.3	587.2
29/01/95	13:21:57	25.7	200.9	45.4	587.4
29/01/95	13:22:57	25.8	201.2	45.5	588.3
29/01/95	13:23:57	25.3	200.7	45.8	588.1

Interfacing with a computer:

Send a single character to the instrument and it will execute the corresponding command as follows:

- 1 = press **MODE**.
- 2 = press **CAL**.
- 3 = press **▲**.
- 4 = press **▼**.
- 5 = press **SET**.
- 6 = press **HOLD**.
- 7 = press **RES**.
- 8 = press **PRINT**.
- ? = send display to computer.
- + = keyboard on (manual operation remains possible).
- = keyboard off (manual operation disabled).
- V = start the input of a value (only valid in certain routines).

Follow exactly the same measuring or calibration procedures, as described in this manual, to include any desired command in your computer data-acquisition program. This example continuously prints the measured values on the computer screen:

```

10  ON ERROR GOTO 70           'makes sure no error occurs
20  OPEN "COM1:2400,N,8,2" AS #1 'prepares computer input
30  PRINT #1,"8"              'PRINT button pressed
40  INPUT #1,V$               'measured values transmitted
50  PRINT V$                  'received values printed
60  GOTO 30                   'repeat procedure
70  RUN

```

Input of a value:

The direct input of values ("V" command) is only possible when the instrument is running a routine where also a manual input of these values could be made. If not, it will ignore any "V" command from the computer!

Examples of allowed routines are:

- temperature input (e.g. manual temperature compensation).
- control parameter input such as Low, High, ...
- time settings.
- input of manual buffer values.

Examples of NOT allowed routines are:

- during measurements.
- when values have to be selected rather than changed (e.g. buffer values in memory).
- while an electrode is being calibrated.

A 16-bit value (2's complement) should be transmitted in the following sequence:

```

1st character = "V" (start a value input)
2nd character = highest byte of the value in ASCII
3rd character = lowest byte of the value in ASCII
4th character = checksum of 2nd and 3rd character in ASCII
5th character = LF, linefeed (ASCII-10)

```

When a correct checksum has been received, the instrument will send a confirmation character "!" to the computer. If not, a question mark "?" is sent. Both are eventually preceded by the identification number. This example makes the computer to enter the value "1000" on instrument #7, in the sophisticated mode:

```

1000 = (H-byte 3)(x 256) + (L-byte 232)
10  OPEN "COM1:2400,N,8,2" AS #1   'prepares communication port
20  PRINT #1,"V";                  'start a value input
30  PRINT #1,CHR$(3);              'highest byte of value
40  PRINT #1,CHR$(232);           'lowest byte of value
50  PRINT #1,CHR$(235);           'checksum of both bytes
60  PRINT #1,CHR$(10);            'LF, closes command

```

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We would like to know your opinion, provide us with your remarks and fax a copy of this page to Consort: +32/14/429179.

1. purchase date:

2. purchased from:

2. indicate model:

4. serial no:

5. organisation:

6. name:

7. division:

8. address:

9. postal code:

10. city:

11. telephone:

12. telefax:

13. e-mail:

- Are you generally satisfied with the performance of the instrument ?
- Are you satisfied with the service of the company you purchased from ?
- What is the degree of satisfaction with:
 - the design of the instrument ?
 - the operating instructions ?
 - the supplied accessories ?
 - the delivery ?
 - the price/quality ratio ?
- What are your suggestions:

☺	☹	☹
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....


.....

Data acquisition software to transfer data to a personal computer:

- For Windows 95 or higher.
- Allows to log, view and print the incoming values in real time or to retrieve data-logged values from the memory of the instrument.
- All channels of multi-channel instruments can be processed simultaneously.
- Values are shown in separate windows, tables or graphs.
- A graph window also allows to set alarm limits and to view several statistic values.

To receive your free copy of the DIS-1 software, you must registrate first!

Make sure to fill in all lines, sign and send, fax or e-mail this page to Consort nv, Parklaan 36, B-2300 Turnhout, Belgium. Fax.: +32/14/429179. E-mail: info@consort.be

Please PRINT in capitals!

1. purchase date:

2. purchased from:

2. indicate model:

T851 - T852

4. serial no:

5. organisation:

6. name:

7. division:

8. address:

9. postal code:

10. city:

11. telephone:

12. telefax:

13. e-mail:

14. I understand that I may use this software only for evaluation and that any other use requires payment of the registration fee.

Signature: