Consort

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

T871
T872

Table of contents

Introduction	1
Manufacturer	1
Warranty	1
Servicing.....	1
AC adaptor	1
Keyboard	1
Inputs	2
Recorder output.....	2
Digital output.....	2
System	2
Recorder.....	2
Password	2
Channels.....	3
Alarm	3
Calibration	3
Measurement	3
Alternating measurements	4
GLP.....	4
Data-logging	5
RS232	5

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Introduction

This instrument is manufactured with the latest technology and needs no particular maintenance. **Consort** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 0° and 40°C. Never use the instrument in a room with high humidity (>95 %) or at very low temperatures (condensation water!).

Manufacturer

Consort bvba

Parklaan 36
B2300 Turnhout
Belgium

Tel (+32)(14)41 12 79
Fax (+32)(14)42 91 79
E-mail: info@Consort.be

Warranty

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **Consort** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **Consort** before returning any product for warranty repair on a freight prepaid basis!

Consort is not liable for consequential damages arising out of the use or handling of its products.

Servicing

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Manual

Cables

Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **Consort** dealer.

AC adaptor

Connect the jack of a suitable AC adaptor (model **A4072**) to the DC socket. Do not hold the adaptor by wet hand.

Keyboard

MODE = Selects all modes or escapes from error traps, calibration procedures, etc..by returning to the original mode.

↔↑↓ = Button for entering a value or for selecting a function.

CAL = Starts or proceeds a calibration or a function.

STORE = Stores the displayed value or send it to a computer.

HOLD = Holds display when measuring.


ON/OFF = Switches the instrument on or off.

T871:

Connect the temperature probes with the terminal blocks and insert them in the instrument. The channel numbers start from the left to the right (backview). The even numbers are the negative poles. The uneven numbers are the positive poles. **Important:** connect the earth terminal to the earthing to avoid ground loop errors!

T872:

Connect the temperature probes with the terminal blocks and insert them in the instrument. The channel numbers start from the left to the right (backview).

**Recorder output**

A recorder can be connected to the red (+) and black (-) terminals. Use only laboratory recorders with a high input impedance!

Digital output

A standard RS232 output terminal (DP9) is provided for interfacing the instrument with a printer or computer. Data is sent in the ASCII code at a BAUD rate of 300...19200 bps (8 bit, no parity, 1 stopbit).

Serial port pinout specifications:

- pin 1: connected to pin 4 and pin 6
- pin 2: TxD, transmit data
- pin 3: RxD, receive data
- pin 4: connected to pin 1 and pin 6
- pin 5: Gnd, signal ground
- pin 6: connected to pin 1 and pin 4
- pin 7: connected to pin 8
- pin 8: connected to pin 7
- pin 9: not connected

System

1. Select the menu by pressing **MODE**.
2. Select [SYSTEM] and press **CAL**.
3. Follow the instructions on the screen to adjust language, contrast, date and time.

Recorder

1. Select the menu by pressing **MODE**.
2. Select [RECORDER] and press **CAL**.
3. Select the desired channel and press **CAL**.
4. Select the desired minimum level corresponding to a 0 V recorder output and press **CAL**.
5. Select the desired maximum level corresponding to a 2 V recorder output and press **CAL**.

Password

1. Select the menu by pressing **MODE**.
2. Select [PASSWORD] and press **CAL**.
3. A private code can be programmed to avoid undesired access to the instrument. Choose [YES] and press **CAL** to enter your secret sequence of 5 keys.

Channels

1. Select the menu by pressing **MODE**.
2. Select [CHANNELS] and press **CAL**.
3. Select [ACTIVATE] and press **CAL**.
4. Select the desired channel and press **CAL** to activate or deactivate.
5. Select [THERMOCOUPLE] and press **CAL**.
6. Select the desired channel and press **CAL**.
7. Select the type of thermocouple and press **CAL**.
8. Select [UNIT] and press **CAL**.
9. Select the desired unit and press **CAL**.

Alarm

1. Select the menu by pressing **MODE**.
 2. Select [ALARM] and press **CAL**.
 3. Select [BEEPER] and press **CAL** to activate or deactivate.
 4. Select the desired channel and press **CAL**.
 5. Select [ACTIVATE] and press **CAL** to activate or deactivate.
 6. Select [TEMP. >] and press **CAL**.
 7. Select the desired maximum level and press **CAL**. Each time the measurements are above this level, the beeper will sound an alarm.
 8. Select [TEMP. <] and press **CAL**.
 9. Select the desired minimum level and press **CAL**. Each time the measurements are below this level, the beeper will sound an alarm.
- *Alarm hysteresis = 1 °C.*

Calibration

1. Verify that the instrument is switched on in a stable temperature environment.
2. Press **CAL** while measuring.
3. Select [CHANNEL] and press **CAL**.
4. Select the desired channel and press **CAL**.
5. Select [METHOD] and press **CAL**.
6. Select the desired calibration method and press **CAL**.
7. Select [TYPE] and press **CAL**.
8. Select the 1- or 2-point calibration and press **CAL**.
9. Select [CALIBRATE] and press **CAL**. Follow the instructions on the screen.

Measurement

1. Verify that the instrument is switched on in a stable temperature environment. The display immediately shows the measured value according to the last calibration.
- *Press \leftrightarrow to change between the channels.*
 - *Press \uparrow to display the maximum temperature ever measured by the selected channel (after switching on). Press \uparrow again to return to the normal measurements. Press **CAL** to reset the minimum/maximum memory.*
 - *Press \downarrow to display the minimum temperature ever measured by the selected channel (after switching on). Press \downarrow again to return to the normal measurements. Press **CAL** to reset the minimum/maximum memory.*
 - *Any alarm will close the relay and sound an alert. Press **MODE** to stop the beeper and to open the relay again.*

1. Select the menu by pressing **MODE**.
2. Select [ALTERNATE] and press **CAL** to activate or deactivate. The display alternates between the selected modes with an interval of ca. 5 s.
3. Press **MODE** to return to the measurements.

GLP

1. Press **CAL** while measuring.
2. Select [GLP] and press **CAL**.
3. Browse with **←→↑↓** to show a complete calibration report. Press **STORE** to print the report.

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1. Select the menu by pressing **MODE**.
2. Select [DATALOGGER] and press **CAL**.

Start the data-logging:

1. Select [STORE] and press **CAL**.
2. Select [ACTIVATE] and press **CAL**.
3. Select YES or NO and press **CAL**.
4. Select [INTERVAL] and press **CAL**.
5. Select the desired time interval between the data-logging (or MANUAL) and press **CAL**.
6. Select [ALARM] and press **CAL**.
7. Select the desired time interval between the data-logging (or MANUAL), when the alarm is activated, and press **CAL**.
8. Select [NUMBER] and press **CAL**.
9. Select the desired number of values to be data-logged and press **CAL**.
10. Press **MODE** to return to the measurements and to start the data-logging according to the previous settings. When manual data-logging has been selected, press **STORE** to put a next measurement into memory. In the meantime the display shows the logging-number e.g. [#00027].

View the stored values on the display:

1. Select [RECALL] and press **CAL**.
2. Select [TABLE] or [PRINT] to display or print the stored data and press **CAL** to continue. Follow the instructions on the screen.
3. In case [TABLE] have been selected press **CAL** to show the options.
4. Select [LOG#] and press **CAL**. Select the number of the first line to be shown on the display and press **CAL**.
5. Select [SHOW] and press **CAL**. Select [Log#] (line numbers) or [Date] to be shown on the display and press **CAL**.

Erase the stored values:

1. Select [ERASE] and press **CAL**. Follow the instructions on the screen.

1. Select the menu by pressing **MODE**.
2. Select [RS232] and press **CAL**.
3. Select [BAUDRATE] and press **CAL**.
4. Select the desired baudrate and press **CAL**.
5. Select [INTERVAL] and press **CAL**.
6. Select the desired interval between the transmitted data and press **CAL**. Pre-set to zero if no automatic transmitting is required.
7. Select [ALARM] and press **CAL**.
8. Select the desired interval between the transmitted data, when the alarm is activated, and press **CAL**. Pre-set to zero if no automatic transmitting is required.
9. Select [IDENTIF. No.] and press **CAL**.
10. Enter an identification number for the transmitted data and press **CAL**.

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DECLARATION OF CONFORMITY

We declare under our sole responsibility that the product

Thermometer

content of the type numbers

T871, T872

to which this declaration relates is in conformity

with the following standards

EN61010

LOW VOLTAGE DIRECTIVE 73/23/EEG

EN50081-1

EN50082-1

EN60555-2

EMC DIRECTIVE 89/336/EEG

Turnhout, October 17, 2005

A handwritten signature in black ink, appearing to be 'G. Ooms', is located below the date.

on behalf of **CONSORT nv**

WARRANTY CERTIFICATE

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. Accessories and breakable items such as electrodes are not warranted unless proven to be defective before shipment.

The original purchase order numbers, Consort invoice numbers and serial numbers of the products must be provided.

CONSORT will repair all defective equipment returned to it during the warranty period without charge (CIF Turnhout prepaid by sender), provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures described in the instruction manual and when not having been subject to accident, alteration, misuse or abuse.

If the products have been used with or have come into contact with fluids, an MSDS (material safety data sheet) must be supplied prior to issuing a return authorisation.

A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

Consort

www.consort.be

Parklaan 36 • B-2300 Turnhout • Belgium

Tel: (+32) (0)14 41 12 79

Fax: (+32) (0)14 42 91 79

Sales: sales@consort.be

Support: support@consort.be

Information: info@consort.be