CONSORT

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

Z921

November 2002

QUALITYMASTERS

CERTIFICATE OF REGISTRATION

This is to certify that:

Consort NV

Parklaan 36
B-2300 Turnhout

Has been assessed by QualityMasters in respect of its Quality Management System and found to comply with:

NEN-EN-ISO 9001 : 2000

Approval is hereby granted for registration providing the rules and conditions relating to certification are observed at all times.

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Lack of fulfillment of conditions as set forth in the certification regulations may render this certificate invalid. The use of the accreditation mark indicates accreditation in respect to the activities covered by the accreditation number EU9706010F.

Q

M

REGISTRATION

CERTIFICATE OF

Authorized signatory

A. LEUWATY
QualityMasters B.V.

This certificate remains the property of QualityMasters B.V.

Instruments manufactured by CONSORT

Electrophoresis power supplies

•

Ion analysers

•

pH meters

•

Conductometers

•

Dissolved oxygen meters

•

Thermometers

•

Industrial controllers

CONSORT nv
Parklaan 36
B2300 Turnhout
Belgium

Find out more at www.consort.be

Tel (+32)(14)41 12 79

Fax (+32)(14)42 91 79

E-mail: info@consort.be

DECLARATION OF CONFORMITY

We declare under our sole responsibility that the product

**Multi-channel analysers
content of the type numbers**

Z921

to which this declaration relates is in conformity

with the following standards

EN61010

LOW VOLTAGE DIRECTIVE 73/23/EEG

EN50081-1

EN50082-1

EN60555-2

EMC DIRECTIVE 89/336/EEG

Turnhout, May 12, 1998

A handwritten signature in black ink, appearing to be 'J. G. G.', is written over a horizontal line.

on behalf of **CONSORT nv**

Table of contents

Introduction	5
Manufacturer	5
Warranty	5
Servicing.....	5
Specifications.....	6
Keyboard	6
Batteries	6
AC adaptor	6
Inputs	6
Recorder output.....	7
Digital output.....	7
System	7
Date and Time.....	7
Recorder.....	7
Password	7
mg/l measurement	8
%O ₂ measurement.....	8
Air pressure measurement	8
Temperature measurement	8
Theory	9
Maintenance	9
GLP.....	10
Data-logging	11
RS232	11

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Forms

Technical assistance form	33
Your satisfaction	34
DIS data information system	35

This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 4° and 40° C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

Manufacturer

CONSORT nv

Parklaan 36

B2300 Turnhout

Belgium

Tel (+32)(14)41 12 79

Fax (+32)(14)42 91 79

E-mail: info@consort.be

Warranty

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

CONSORT is not liable for consequential damages arising out of the use or handling of its products.

Servicing

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Manual

Cables

Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

Ranges	O₂	0...60 mg/l
	Saturation	0...600 %
	Air pressure	800...1100 hPa
	°C	0...50 °C
Resolution	0.01 mg/l, 0.1 %, 0.1 °C	
Inputs	BNC + BANANA input for oxygen electrode	
Temp. comp.	automatic with Pt1000	
Pressure comp.	automatic with built-in sensor	
Digital output	programmable RS232, 300...9600 b/s, for bi-directional communication with a computer or printer	
Analogue output	0...2 V, ca 8 kΩ, 8 bit D/A	
Display	LCD display, 128 x 64 pixels	
Keys	8 tactile keys	
Ambient temp.	4...40 °C	
Rel. Humidity	0...90 % (non-condensing!)	
Power supply	210-250 V~ (105-125 V~), 50/60 Hz, max. 2 VA and 4 NiMH batteries, size AA	
Cabinet	IP65 cabinet	
Dimensions	252 x 121 x 50 mm	
Weight	600 g	

Keyboard

MODE	=	Selects all modes or escapes from error traps, calibration procedures, etc..by returning to the original mode.
▲/▼/◀/▶	=	Button for entering a value or for selecting a function.
CAL	=	Starts or proceeds a calibration or a function.
PRINT	=	Displayed value is printed through the RS232 output.
HELP	=	Built-in manual.
HOLD	=	Holds display when measuring.
ON/OFF	=	Switches the instrument on or off.

Batteries

The batteries can be replaced by opening the bottom compartment of the cabinet. Only use NiMH batteries!

AC adaptor

Remove the protective silicone stopper from the DC socket and connect the jack of a suitable AC adaptor (model **A4070** for 230 V~ or model **A4071** for 115 V~) to the DC socket for recharging the NiMH batteries. **Avoid to use an AC adaptor without NiMH batteries being inserted in the instrument!** For field work, remove the adaptor and re-insert the protective silicone stopper into the DC socket. Do not hold the adapter by wet hand!

Inputs

The measuring electrode should be connected to the coaxial INPUT connector. Its banana plugs should be inserted in the °C terminals.

Recorder output

A recorder can be connected to the red (+) and black (-) terminals. Use only laboratory recorders with a high input impedance!

Digital output

A standard RS232 output terminal (DP9) is provided for interfacing the instrument with a printer or computer. Data is sent in the ASCII code at a BAUD rate of 150...4800 bps (8 bit, no parity, 1 stopbit).

Serial port pinout specifications:

- pin 1 : connected to pin 4 and pin 6
- pin 2 : TxD, transmit data
- pin 3 : RxD, receive data
- pin 4 : connected to pin 1 and pin 6
- pin 5 : Gnd, signal ground
- pin 6 : connected to pin 1 and pin 4
- pin 7 : connected to pin 8
- pin 8 : connected to pin 7
- pin 9 : not connected

System

1. Select [MAIN MENU] by pressing **MODE**.
2. Select [SET-UP], [SYSTEM] and press **CAL**.
3. Follow the instructions on the screen to adjust language, contrast and automatic power-off timer.

Date and Time

1. Select [MAIN MENU] by pressing **MODE**.
2. Select [SET-UP], [DATE/TIME] and press **CAL**.
3. Adjust the date (respectively year, month and day) by pressing **CAL** each time.
4. Adjust the time (respectively hour, minutes and seconds) by pressing **CAL** each time.

Recorder

1. Select [MAIN MENU] by pressing **MODE**.
2. Select [SET-UP], [RECORDER] and press **CAL**.
3. Select the desired range and press **CAL**.
4. Select the desired minimum level corresponding to a 0 V recorder output and press **CAL**.
5. Select the desired maximum level corresponding to a 2 V recorder output and press **CAL**.

Password

1. Select [MAIN MENU] by pressing **MODE**.
2. Select [SET-UP], [PASSWORD] and press **CAL**.
3. A private code can be programmed to avoid undesired access to the instrument. Enter your secret sequence of 5 keys and press **CAL**.

1. Select [MAIN MENU], [mg/l] by pressing **MODE**. The display will immediately show the measured value according to the previous calibration. Should you want to recalibrate, press **CAL**.
2. The display shows the salinity correction. Select the proper values and press **CAL**. Leave salinity correction to zero unless you are going to measure in heavily salted solutions such as e.g. sea-water (35 g/l). Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
3. The electrode exposed to the air reaches an equilibrium corresponding to the partial pressure of oxygen and thus to saturation in water at the given temperature. The instrument shows the measured saturation, current, temperature, and will calibrate automatically when readings are stable.
4. After rinsing the electrodes with distilled water, immerse them in the samples and read the display. Stirring the solution during the measurements promotes the homogeneity and is **obligatory!** The advection rate must be at least 10 cm/s.
5. Rinse the electrode always after use and store it in distilled water.

Relative measurements:

1. While in the O₂ range, press \wedge . At that moment the instrument stores the actual value and automatically subtracts it from all next measurements. You can repeat this as many times as required. [REL] in the display appears.
 2. To cancel relative readings, press **MODE** and you will be in the normal measuring mode again.
- *A blinking decimal point warns you for unstable measurements. Wait to read the display!*
 - *Stirring the solution during the measurements promotes the homogeneity and is therefore always recommended.*
 - *Press ∇ to change the resolution from 0.1 to 0.01 mg/l or from 1 to 0.1 %.*

%O₂ measurement

1. Select [MAIN MENU], [%O₂] by pressing **MODE**. The display will immediately show the measured value according to the previous calibration. Should you want to recalibrate, press **CAL**. Proceed as for mg/l.

Air pressure measurement

1. Select [MAIN MENU], [hPa] by pressing **MODE**. Should you want to recalibrate, press **CAL**.
2. Adjust to the real air pressure and press **CAL**.

Temperature measurement

1. Select [MAIN MENU], [°C] by pressing **MODE**. Should you want to recalibrate, press **CAL**.
2. Immerse the oxygen electrode in a standard solution of known temperature. Calibrate to the proper value and press **CAL**.

Principle:

The oxygenmeter and its electrode function according to the Clark principle with silver as cathode and zinc as anode in an electrolyte cell. Oxygen gas present in the electrolyte is reduced to OH ions at the cathode. The resulting current is diffusion limited and therefore proportional to the oxygen concentration in the sample solution. This current is amplified, corrected, and displayed in mg/l, ppm or % dissolved oxygen.

Electrode:

The anode and cathode of the electrode are immersed in a special buffer solution inside the electrode, to give high stability readings and short polarisation times. The electrolyte compartment is separated by a PTFE film from the sample solution. Only gases, but no ions, can diffuse through this film from the sample solution into the internal electrolyte.

Temperature compensation:

Speed of electrode reactions depends on temperature. Therefore, for a given oxygen concentration, the electrode current becomes greater at higher temperature, and conversely. This temperature effect is compensated with the Pt1000 resistor incorporated in the electrode. For aqueous solutions the oxygen meter therefore always gives correct readings, independent of the actual temperature.

Interferences:

All substances which can diffuse through the membrane and for which 800 mV potential suffices for polarographic reduction, will be reduced in the electrode. This will give a corresponding current contribution, if they are present. Interference can be caused by ions entering the electrode through porous or mechanically damaged membranes and by diffusion of other reactive gases apart from oxygen, e.g. CO₂, Cl₂, SO₂, and H₂S. These substances react in undesired manner with the electrode. Acidic or basic gases change the pH value of the electrolyte solution and thus disturb the reading, particularly when measuring small oxygen concentrations. High salt concentrations in the sample solution can falsify readings too.

Maintenance

The electrode body has almost unlimited life and requires no maintenance. After prolonged use of the electrode, the silver anode may become deactivated by undesired deposits. An indication is that the electrode no longer responds correctly to calibration. In this case:

1. Carefully screw off the measuring head.
2. Cautiously remove the precipitated deposits with the aid of filter paper. **Do not use grinding paper or a glass fibre!**
3. Rinse several times with distilled water and shake carefully off the water drops.
4. Take the measuring head, fill it with new electrolyte solution and very slowly screw it onto the electrode while holding it vertically.
5. **Let the electrode rest for a few hours!** The electrode is now ready for use again.

Replace the membrane only when damaged, not when calibration is no longer possible!

1. Select [MAIN MENU] by pressing **MODE**.
2. Select [SET-UP], [GLP] and press **CAL**.
3. Select [SHOW REPORT] and press **CAL**. Browse with **▲**, **▼**, **◀** or **▶** to show a complete calibration report. Press **PRINT** to print the report.

CONSORT O₂ meter Z921

Date : 15/02/1998
Time : 11:32:04
Version : 1.0

SETTINGS

Identification No : 003
Password : OFF
Salinity (g/l) : 0.00

O₂ CALIBRATION

Date : 14/02/1998
Time : 03:47:22
Cal. points : 2

Temperature (°C) : 19.4
Air pressure (hP) : 1035
Sat. point (mg/l) : 7.48
Sat. current (µA) : 8.63
Response time (s) : 34

STATISTICS

Temperature (°C) : +0.7
Air pressure (hP) : -5
Sat. current (µA) : -1.7
Response time (s) : +12

Calibration reminder:

1. Select [INTERVAL] and press **CAL**.
2. Select the desired time interval between each automatic warning for a new calibration of the electrodes and press **CAL**.

1. Select [MAIN MENU] by pressing **MODE**.
2. Select [SET-UP], [DATALOGGER] and press **CAL**.

Start the data-logging:

1. Select [STORE] and press **CAL**.
2. Select [RANGE] and press **CAL**.
3. Select the desired range and press **CAL**.
4. Select [RESOLUTION] and press **CAL**.
5. Select the desired resolution and press **CAL** (high = 16 bit, low = 8 bit).
6. Select [INTERVAL] and press **CAL**.
7. Select the desired time interval between the data-logging and press **CAL**.
8. Select [NUMBER] and press **CAL**.
9. Select the desired number of values to be data-logged and press **CAL**.
10. Select [START DATA-LOGGING] and press **CAL**.
11. Starts the data-logging according to the previous settings while a blinking [LOG] appears. When manual data-logging has been selected, press **PRINT** to put a next measurement into memory. In the meantime the display shows the logging-number e.g. [#0027].

View the stored values on the display:

1. Select [RECALL] and press **CAL**.
2. Select the desired form to display or print the stored data and press **CAL** to continue. The scales of a graph can be modified by pressing **CAL**. Follow the instructions on the screen.

Erase the stored values:

1. Select [ERASE] and press **CAL**. Follow the instructions on the screen.

RS232

1. Select [MAIN MENU] by pressing **MODE**.
2. Select [SET-UP], [RS232] and press **CAL**.
3. Select [BAUDRATE] and press **CAL**.
4. Select the desired baudrate and press **CAL**.
5. Select [INTERVAL] and press **CAL**.
6. Select the desired interval between the transmitted data and press **CAL**. Pre-set to zero if no automatic transmitting is required.

Remote control by a computer:

Send a series of characters to the instrument and it will execute the corresponding command as follows:

Command	Instruction	Answer
">1"	press MODE	"<1"
">2"	press ▲	"<2"
">3"	press ▼	"<3"
">4"	press CAL	"<4"
">5"	press PRINT	"<5"
">6"	press HELP	"<6"
">7"	press HOLD	"<7"
">?"	send displayed measurement or selection	"<?"
">+"	switch keyboard on	"<+"
">-"	switch keyboard off	"<-"
">M"	go directly to Set-up	"<M"
">G"	go directly to mg/l measurement	"<G"
">%"	go directly to %O₂ measurement	"<%"
">T"	go directly to °C measurement	"<T"
">H"	go directly to hPa measurement	"<H"
">Nxy"	input of value "xy" (x = low byte / y = high byte)	"<Nxy"

Follow exactly the same measuring or calibration procedures, as described in this manual, to include any desired command in your own software.

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If you need assistance with a technical problem, make sure to fill the **complete** checklist (20 points!) and fax a copy to Consort: +32/14/429179.

1. purchase date:

2. purchased from:

2. indicate model:

4. serial no:

5. organisation:

6. name:

7. division:

8. address:

9. postal code:

10. city:

11. telephone:

12. telefax:

13. e-mail:

14. type of connected electrode:

Perform the following steps:

15. Carefully screw off the measuring head of the oxygen electrode. Cautiously remove the precipitated deposits with the aid of filter paper. Do not use grinding paper or a glass fibre! Rinse several times with distilled water and shake carefully off the water drops. Take the measuring head, fill it with new electrolyte solution and very slowly screw it onto the electrode while holding it vertically. Let the electrode rest for a few hours! yes?

16. reset the meter (switch on while holding MODE pressed). yes?

17. short-circuit the O₂ input, select the O₂ mode and read display. mg/l:

18. select the °C mode, measure air temperature and read the display. °C:

19. select the O₂ mode (do not calibrate!). yes?

20. rinse the O₂ electrode, hold in the air and read display. mg/l:

We would like to know your opinion, provide us with your remarks and fax a copy of this page to Consort: +32/14/429179.

1. purchase date: 2. purchased from:

2. indicate model: 4. serial no:

5. organisation:

6. name:

7. division:

8. address:

9. postal code: 10. city:

11. telephone:

12. telefax:

13. e-mail:

- | | ☺ | ☹ | ☹ |
|--|--------------------------|--------------------------|--------------------------|
| • Are you generally satisfied with the performance of the instrument ? | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| • Are you satisfied with the service of the company you purchased from ? | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| • What is the degree of satisfaction with: | | | |
| the design of the instrument ? | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| the operating instructions ? | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| the supplied accessories ? | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| the delivery ? | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| the price/quality ratio ? | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| • What are your suggestions: | | | |

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....


.....

Data acquisition software to transfer data to a personal computer:

- For Windows 95 or higher.
- Allows to log, view and print the incoming values in real time or to retrieve data-logged values from the memory of the instrument.
- All channels of multi-channel instruments can be processed simultaneously.
- Values are shown in separate windows, tables or graphs.
- A graph window also allows to set alarm limits and to view several statistic values.

To receive your free copy of this Data Information System, make sure to fill in all lines, sign and send or fax this page to Consort nv, Parklaan 36, B-2300 Turnhout, Belgium.

Fax.: +32/14/429179

Please PRINT in capitals!

1. purchase date:

2. purchased from:

2. indicate model:

4. serial no:

5. organisation:

6. name:

7. division:

8. address:

9. postal code:

10. city:

11. telephone:

12. telefax:

13. e-mail:

14. I understand that I may use this software only for evaluation and that any other use requires payment of the registration fee.

Signature: