

- On screen help in 4 languages
- Real Time Clock
- Various running modes:
 - Simple Mode
 - 9x9 Method Programming Mode
 - Voltage Ramp Mode
 - Timer Mode
- Constant voltage/Constant current/Constant power
- Automatic cross-over
- Automatic recovery after power failure
- Password protection
- Data-logging
- Remote control via open communication protocol
- Safety features:
 - Ground leakage detection
 - Overload protection
 - Smooth voltage rise
 - No load detection
 - Isolated USB communication
- Free software and firmware updates




EVS3100-SYS



EHCA1200-SYS



BLOOK



EHCA1200 Strips Sheets

The EV1450 is a small power supply suitable for most smaller tanks and applications. EV1450 has a firmware upgrade capability so future improvements and features will always be available.

Moreover EV1450 has a continuous logging combined with a real time clock so it's possible to get an overview of previous runs, including possible down-times in case of mains power failures.

Consort Power Supplies are the most robust, long lasting and durable electrophoresis power supplies in the market.

Power 50 W
 Voltage 400 V
 Current 500 mA
 Outputs 4
 Operating Modes 4
 Simple Mode
 9x9 Method Programming Mode
 Voltage Ramp mode
 Timer Mode (time or Vh)

USB interface
 Datalogging
 Real Time Clock
 Multiple safety features

Warranty 36 months
 Made in Belgium

Code	Description
EV0220	Power supply, 200 V, 200 mA, 20 W
EHS3100-SYS	Mini horizontal unit 7x8 cm
EHCA1200-SYS	Horizontal unit for cellulose acetate
EHCA1200-ST	Cellulose acetate strips, available in different sizes
EHCA1200-SH	Cellulose acetate sheets, available in different sizes
BLOOK	LED gel documentation table
→ Supplied with a european mains cord + USB cable	

Power	50 W
Voltage	400 V
Current	500 mA
Outputs	4
Operating Modes	4 Simple Mode 9x9 Method Programming Mode Voltage Ramp mode Timer Mode (time or Vh)
USB interface	
Datalogging	
Real Time Clock	
Multiple safety features	
Warranty	36 months
Made in Belgium	



● Description

The EV1450 is a small power supply suitable for most smaller tanks and applications. The front panel and graphical display are designed for ease of use. The display provides all useful information during runs and will show an on screen help to guide the user in setting up the power supply. In Simple Mode you just have to set your power supply to the desired parameters and press run. EV1450 has a firmware upgrade capability so future improvements and features will always be available. Moreover EV1450 has a continuous logging combined with a real time clock so it's possible to get an overview of previous runs, including possible down-times in case of mains power failures. Consort Power Supplies are the most robust, long lasting and durable electrophoresis power supplies in the market.

● Features

On screen help in 4 languages to assist setting up the power supply parameters and solve errors.

Firmware updates allows for upgrades to the latest version via the USB interface. Feature requests can also be implemented via the firmware system.

Real Time Clock date and time are kept in a battery backup system and is used logging an electrophoresis run.

Various running modes:

Simple Mode: just set voltage, current, power and time for a routine electrophoresis run.

9x9 Method Programming Mode: Up to 9 different programs, each with 9 steps, can be stored in the non-volatile memory. Each step is able to recall a next one, providing a flexible multiple step function for special techniques. Parameters of the running step can be changed temporarily without interrupting the run.

Voltage Ramp Mode: a linear voltage gradient for any step provided the limiting current or power is not attained.

Timer Mode: Timer or volt-hour controlled operation will automatically stop the run and sound an alarm.

Automatic cross-over Each model has constant voltage, constant current, constant power capabilities with automatic cross-over and shows which parameter is kept constant.

Automatic recovery after power failure

Password protection

Data-logging Data logging of about 100 hours of runs are automatically stored. Data includes data/time, voltage, current, power and date/time of following events: start, stop, pause, program number, step, changes, mains failure and auto restart.

Data Transfer Free data acquisition software for PC can be downloaded from our website. It allows to visualize and examine the stored run details.

Remote control EV1450 can be controlled by a computer using special commands. These commands can be found in the support section of our website.

Safety features:

Ground leakage detection: protection from potential shock hazard when a ground leakage path is detected.

Overload protection: full protection against any overload condition including accidental short circuit of the output.

Smooth voltage rise: high voltage cannot suddenly appear at the outputs but will increase smoothly up to the pre-set limits.

No load detection: prevents errors such as a bad or a dangling connection.

Isolated communication: Optically isolated USB input/output connection to prevent any high voltage on the communication lines.

● Specifications

VOLTAGE	0...400 V	PROGRAMS	9x9 set of parameters	DATA-LOGGING	3600 values
CURRENT	0...500 mA	OUTPUTS	4 in parallel, 4 mm sockets	INTERVAL	1...60 seconds
POWER	0...50 W	MINIMUM LOAD	30 Ω	REAL TIME CLOCK	✓
PARAMETER RANGE	1...100% of full scale	GROUND LEAKAGE DETECTION	✓	AMBIENT TEMPERATURE	0...40°C
TIMER	0...99:59 h	OVERLOAD DETECTION	✓	RELATIVE HUMIDITY	0...95%, non condensing
VOLT-HOURS	0...99.99 kWh	COMPUTER CONTROL	✓	POWER REQUIREMENTS	210-250 VAC, 50/60 Hz, 75 W 100-125 VAC, 50/60 Hz, 75 W
SETUP RESOLUTION	1 V, 1 mA, 1 W	USB INTERFACE	✓	DIMENSIONS (WxDxH)	24x20x13 cm
MEASUREMENT RES.		PASSWORD	✓	WEIGHT	3 kg
		DISPLAY	graphical		