LabArt C51

**MANUAL
HANDLEIDING
MODE D'EMPLOI
ANLEITUNG**
Rev. 03/2021

C51

Input CH1: Conductivity, TDS, SAL
Input CH2: pH, mV, Dissolved Oxygen
Input °C: Temperature (Pt1000)

Index

General information	1
Set-up	2
Measurement	2
Calibration	2
Good Laboratory Practice	3
Data Storage in the internal memory	3
Algemene informatie	4
Instellen	5
Metingen	5
Kalibreren	5
Goede Laboratorium Praktijk	6
Data-opslag in het intern geheugen	6
Information générales	7
Configure	8
Mesure	8
Calibration	8
Bonnes Pratiques de Laboratoire	9
Collection de données dans la mémoire interne	9
Allgemeine Information	10
Einstellen	11
Messung	11
Kalibrierung	11
Gute Laborpraxis	12
Datenspeicherung im internen Speicher	12
Specifications	13
NIST pH buffers	14
Standard pH buffers	14
Conductivity standards	14
Redox standards	14

General Information

Introduction

This instrument is manufactured with the latest technology and needs no particular maintenance. **Consort** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 0° and 40°C. Never use the instrument in a room with high humidity (>95 %) or at very low temperatures (condensation water!).

Manufacturer

Consort bvba
Parklaan 36 Tel (+32)(14)41 12 79
B2300 Turnhout Fax (+32)(14)42 91 79
Belgium E-mail: info@consort.be

Warranty

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **Consort** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **Consort** before returning any product for warranty repair on a freight prepaid basis!

Consort is not liable for consequential damages arising out of the use or handling of its products.

Servicing

In the event of this instrument being returned for servicing, the owner is requested **NOT** to send the following items unless they are suspect:

Manual
Cables
Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **Consort** dealer.

Batteries

The batteries can be replaced by opening the bottom compartment of the cabinet. Only use NiMH batteries!

AC adaptor

Connect the mains adapter to the DC socket.

Keyboard

F1, F2 = Variable function, selects the function displayed on the screen
OK = Validate
STORE = Store measurement
↑↓ = Button for entering a value or for selecting a function.
ON/OFF = Switches the instrument on or off.

Set-up

1. Select [SETTINGS] by pressing **MENU (F2)**.
2. Follow the instructions on the screen to adjust language, contrast and automatic power-off timer or to erase the memory. Press **OK or SELECT (F2)** to confirm or **CANCEL (F1)** to exit.

Measurement

1. Connect the electrodes to the corresponding BNC inputs. Automatic temperature compensation and temperature measurements are possible with a Pt1000 temperature probe.
2. Press MENU (F2), select the desired mode with **↑↓** and press **OK or SELECT (F2)**.

Never immerse a pH electrode and a dissolved oxygen electrode together in the same solution!

A blinking decimal point warns you for unstable measurements. Wait to read the display! Stirring the solution during the measurements promotes the homogeneity and is obligatory!

Calibration

Press **CALIBRATE (F1)** to go to the calibration menu

The instrument will refuse automatic calibration when the electrode is unstable. Insufficient stirring or a worn electrode may be the cause.

pH

1. The instrument permits to choose between the buffers in memory (1.68, 2.00, 4.00, 4.01, 6.87, 7.00, 9.18, 9.21, 9.95, 10.01, 12.00, 12.45). Select the proper values and press **OK**.
2. Rinse the electrodes with distilled water and immerse them in the first buffer solution. Select [CALIBRATE], press **OK or SELECT (F2)** and follow the instructions on the screen until the calibration is finished.
3. After rinsing the electrodes with distilled water, immerse them in the samples and read the display.
4. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.

mV

1. Immerse the electrodes in a standard solution of known potential. Select [CALIBRATE] and press **OK or SELECT (F2)**. Adjust to the proper value and press **OK or SELECT (F2)**.
 - Select [RESET] and press **OK or SELECT (F2)** to reset the calibration.

Temperature

1. Immerse the Pt1000 in a solution of known temperature. Adjust to the proper value and press **OK**.
 - Select [RESET] and press **OK** to reset the calibration.
 - Without Pt1000, press **OK**, adjust the manual temperature compensation and proceed by pressing **OK**.

Conductivity

1. The instrument permits to choose between the standards in memory (0.01, 0.1, 1 M KCl). Select the proper values and press **OK or SELECT (F2)**
2. After rinsing the electrode several times with the first standard solution, immerse it in that standard. The solution temperature is not so critical but should lie between 0°C and 40°C. When no Pt1000 is used, do not forget to compensate manually first! Select [CALIBRATE], press **OK or SELECT (F2)** and follow the instructions on the screen until the calibration is finished.
3. Rinse the electrode several times with the sample, immerse it in that solution and read the display.
4. Rinse the electrode always after use and store it in distilled water.

Dissolved oxygen

1. The display shows the air pressure correction and the salinity correction. Select the proper values and press **OK or SELECT (F2)**. Leave salinity correction to zero unless you are going to measure in heavily salted solutions such as e.g. sea-water (35 g/l). Select [CALIBRATE], press **OK or SELECT (F2)** and follow the instructions on the screen until the calibration is finished.

2. The electrode exposed to the air reaches an equilibrium corresponding to the partial pressure of oxygen and thus to saturation in water at the given temperature. The instrument shows the measured saturation, current, temperature, and will calibrate automatically when readings are stable.
3. After rinsing the electrodes with distilled water, immerse them in the samples and read the display. Stirring the solution during the measurements promotes the homogeneity and is **obligatory!** The advection rate must be at least 10 cm/s.
4. Rinse the electrode always after use and store it in distilled water.

Good Laboratory Practice

1. Press CALIBRATE (F1) to go to the calibration menu
2. Select GLP and press **OK or SELECT (F2)** to display a complete calibration report.

Data Storage in the internal memory

Storing data

1. Press **STORE** to store the displayed value in memory. The display shows a next storage address.
2. Select the desired address and press **OK or SELECT (F2)** to enter the measured value in memory.

Recalling data

1. Press **↑** to recall the stored data. The display shows a list of stored values.
2. Select the desired address to read the stored values.
3. Press **OK** to return to the normal measurements.

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Specifications		LabArt C51
pH	<i>Range</i>	0...14 pH
	<i>Resolution</i>	0.01 pH
	<i>Accuracy</i>	0.2% ± 1 digit
	<i>Calibration</i>	1...3 points
	<i>Buffers</i>	12 pre-programmed
	<i>Temperature Compensation</i>	0...100° C
	<i>ISO-pH</i>	6...8 pH
	<i>Slope</i>	80...120%
mV	<i>Range</i>	±2000 mV
	<i>Resolution</i>	1 mV
	<i>Accuracy</i>	0.2% ± 1 digit
	<i>Calibration</i>	1 point
CONDUCTIVITY	<i>Range</i>	0...1000 mS/cm
	<i>Resolution</i>	0.01 µS/cm
	<i>Accuracy</i>	1% f.s. of range
	<i>Calibration</i>	1 point
	<i>Standards</i>	0.01/0.1/1M KCl
	<i>Cell Constant</i>	0.1/1/10 cm ⁻¹ ± 30%
	<i>Temperature Compensation</i>	0...100° C
	<i>Reference Temperature</i>	20° or 25° C
	<i>Temperature Coefficient</i>	natural waters (EN27888)
SALINITY	<i>Range</i>	0...70 ppt
	<i>Resolution</i>	0.1 ppt
TDS	<i>Range</i>	0...100 µg/l
	<i>Resolution</i>	0.1 mg/l
DISSOLVED OXYGEN	<i>Range</i>	0...20 mg/l (0...200%)
	<i>Resolution</i>	0.05 mg/l (0.5%)
	<i>Accuracy</i>	1% ± 1 digit
	<i>Calibration</i>	1 point
	<i>Temperature Compensation</i>	0...50° C
	<i>Salinity Compensation</i>	0...40 ppt
	<i>Air Pressure Compensation</i>	800...1200 hPa
TEMPERATURE	<i>Range</i>	-5...100° C
	<i>Resolution</i>	0.1° C
	<i>Accuracy</i>	0.5° C
	<i>Calibration</i>	1 point
INPUTS	<i>pH/mV</i>	BNC, 10 ¹² Ω
	<i>Conductivity/Dissolved Oxygen</i>	BNC
	<i>Temperature</i>	2 banana, for Pt1000
STORAGE MEMORY	<i>Data sets</i>	50
DISPLAY	<i>LCD</i>	160x160 pixels
	<i>White Backlight</i>	✓
	<i>Languages</i>	EN, NL, FR, DE
POWER SUPPLY	<i>Low Voltage</i>	9...15 VDC
	<i>Batteries</i>	4x1.2 V, NiMH
DIMENSIONS	<i>WxDxH</i>	10x20x4 cm
WEIGHT	<i>Meter</i>	350 g

NIST pH buffers

Temp.	1.68	4.01	6.87	9.18	12.45
0°C	1.666	4.003	6.984	9.464	13.423
5°C	1.668	3.999	6.951	9.395	13.207
10°C	1.670	3.998	6.923	9.332	13.003
15°C	1.672	3.999	6.900	9.276	12.810
20°C	1.675	4.002	6.881	9.225	12.627
25°C	1.679	4.008	6.865	9.180	12.454
30°C	1.683	4.015	6.853	9.139	12.289
40°C	1.694	4.035	6.838	9.068	11.984
50°C	1.707	4.060	6.833	9.011	11.705
60°C	1.723	4.091	6.836	8.962	11.449
70°C	1.743	4.126	6.845	8.921	
80°C	1.766	4.164	6.859	8.885	

Standard pH buffers

Temp.	2.00	4.00	7.00	9.21	10.01	12.00
5°C	1.85	3.99	7.08	9.45	10.24	12.56
10°C	1.89	3.99	7.06	9.38	10.18	12.35
20°C	1.96	3.99	7.02	9.26	10.06	12.14
25°C	2.00	4.00	7.00	9.21	10.01	12.00
30°C	2.03	4.01	6.99	9.16	9.87	11.89
40°C	2.08	4.03	6.98	9.06	9.89	11.71
50°C	2.12	4.06	6.97	8.99	9.83	11.52
60°C	2.14	4.09	6.98	8.93	9.79	11.33

Conductivity standards

Temp.	0.01 M	0.1 M	1 M	(mS/cm)
0°C	0.776	7.15	65.4	
5°C	0.896	8.22	74.1	
10°C	1.020	9.33	83.2	
15°C	1.147	10.48	92.5	
20°C	1.278	11.67	102.1	
25°C	1.413	12.88	111.8	
30°C	1.552	14.12	121.7	
35°C	1.694	15.38	131.7	
40°C	1.839	16.67	142.1	

Redox standards

Temp.	124	358	(mV)
0°C	156	388	
5°C	148	381	
10°C	141	376	
15°C	135	370	
20°C	130	364	
25°C	124	358	
30°C	122	352	
35°C	116	346	
40°C	99	340	

DECLARATION OF CONFORMITY

We declare under our sole responsibility that the product

**Multi-channel analysers
content of the type numbers**

LabArt C51

to which this declaration relates is in conformity
with the following standards

**LOW VOLTAGE DIRECTIVE 2006/95/EC
EN61010-1**

**EMC DIRECTIVE 2004/108/EC
EN61326-1**

Consort bvba

Parklaan 36, B-2300 Turnhout, Belgium

Tel: (+32) (0)14 41 12 79

Fax: (+32) (0)14 42 91 79

Sales: sales@consort.be

The electro-magnetic susceptibility has been chosen at a level that gains proper operation in residential areas, on business and light industrial premises and on small-scale enterprises, inside as well as outside of buildings. All places of operation are characterised by their connection to the public low voltage power supply system.

WARRANTY CERTIFICATE

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory.

Accessories and breakable items such as electrodes are not warranted unless proven to be defective before shipment.

The original purchase order numbers, Consort invoice numbers and serial numbers of the products must be provided.

CONSORT will repair all defective equipment returned to it during the warranty period without charge (CIF Turnhout prepaid by sender), provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures described in the instruction manual and when not having been subject to accident, alteration, misuse or abuse.

If the products have been used with or have come into contact with fluids, an MSDS (material safety data sheet) must be supplied prior to issuing a return authorisation.

A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

Consort bv

Hertenstraat 56/9
B-2300 Turnhout • Belgium

Tel : (+32) (0)14 41 12 79
Fax : (+32) (0)14 42 91 79

Sales : sales@consort.be
Support : support@consort.be
Information : info@consort.be