# *Power Supplies*

## *EV2000/EV3000 series*

**MANUAL  
HANDLEIDING  
ANLEITUNG**  
Rev. 03/2021

***EV2230, EV2310, EV2320,  
EV2650, EV3020, EV3150,  
EV3330, EV3610, EV3620***

**Consort**

## Models – Modelen - Modelle

### *Maxima*

	Voltage (V)	Current (mA)	Power (W)
EV2230	1500	200	150
EV2310	300	1000	150
EV2320	3000	100	150
EV2650	600	500	150

	Voltage (V)	Current (mA)	Power (W)
EV3020	300	2000	300
EV3150	1200	500	300
EV3330	3000	200	300
EV3610	600	1000	300
EV3620	6000	100	300

# Contents – Inhoudsopgave – Inhaltsverzeichnis

Models – Modellen - Modelle.....	1
Maxima .....	1
Contents – Inhoudsopgave – Inhaltsverzeichnis .....	b
English.....	1
General Information.....	1
Keyboard.....	1
Very important .....	2
Use of the instrument.....	2
Procedure?.....	3
Setup .....	3
Method programming.....	7
Manual programming .....	9
Running condition .....	10
Alarm messages.....	12
Troubleshooting.....	12
RESET .....	12
USB Digital Port .....	13
Frequently Asked Questions.....	14
Nederlands.....	15
Algemene informatie.....	15
Toetsenbord.....	15
Zeer belangrijk .....	16
Gebruik van het instrument.....	16
Procedure? .....	17
Instellen.....	17
Methode programmatie .....	21
Manuele programmatie .....	23
Verloop van een procedure .....	24
Meldingen .....	26
TROUBLESHOOTING .....	26
RESET .....	26
USB Digitale Poort .....	27
Veelgestelde vragen.....	28
Deutsch.....	29
Allgemeine Informationen .....	29
Bedienfeld .....	29
Sehr wichtig .....	30
Verwendung des Geräts .....	30
Verfahren? .....	31
Einstellen .....	31
Methodenprogrammierung .....	35
Manuelles Programmieren .....	37
Laufbedingungen.....	38

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## General Information

### Introduction

This instrument is manufactured with the latest technology and needs no specific maintenance. **Consort** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 0° and 40°C.

Never use the instrument in a room with high humidity (>95 %) or at very low temperatures (condensation water!).

### Manufacturer

#### Consort bvba

Hertenstraat 56 / 9  
2300 Turnhout  
Belgium

Tel (+32)(0)14 41 12 79  
Fax (+32)(0)14 42 91 79  
E-mail [info@consort.be](mailto:info@consort.be)

### Warranty

This instrument (excluding all accessories) is warranted against defective material and workmanship for period of thirty-six (36) months from the date of shipment ex-factory. **Consort** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **Consort** before returning any product for warranty repair on a freight prepaid basis!

**Consort** is not liable for consequential damages arising out of the use or handling of its products.

### Servicing

In the event of this instrument being returned for servicing, the owner is requested **NOT** to send the following items unless they are suspect:

- *Manual*
- *Cables*
- *Accessories*

If serious malfunctioning occurs, stop using the unit immediately and consult your local **Consort** dealer.

### Keyboard

- During **SETUP** or programming a procedure:

- |  |   |
|--|---|
|  | = Buttons for entering a value or for selecting a function.   |
|  | = Button for programming and confirming the desired parameters.<br><b>Hold this button pressed to return to the main menu!</b>  |
|  | = Button to <ul style="list-style-type: none"><li>• start the procedure (MANUAL or METHOD),</li><li>• cancel while entering a parameter (SETUP),</li><li>• return to the previous menu (SETUP).</li></ul> |

- While running a **PROCEDURE**:

- |  |  |
|--|--|
|  | = Buttons for selecting a display configuration. |
|  | = Button to modify the parameters.               |
|  | = Button to stop or to pause the experiment.     |

### *Very important*

This electrophoresis power supply is a high technology instrument available in several versions. As it is capable of giving dangerous voltage levels by which high power is involved, we suggest that you take a few moments to read this manual thoroughly. Although this instrument is equipped with all necessary safety features against abuse and other accidental failures, caution should be exercised when working with high voltage equipment.

1. **Avoid to touch the outlets with any conducting object** and make sure there is a second person present for your safety in case of any severe electric shock.
2. **Never touch** any part of the assembly (power supply, leads or tank) **before having switched off**.
3. **Never manipulate with wet hands.**
4. **Do not ground any of the outputs** or the buffer in the tank.
5. Connect the outlets only to an insulated electrophoresis tank with safety cover.
6. Never make any other connections, such as e.g. putting several power supplies in series or in parallel.
7. In order to prevent electric shock, never open the back plate nor remove the cover.
8. Do not expose the unit to rain or any other liquid.
9. **Do not spill liquid or insert metal objects inside the unit.**
10. Take care so that the power supply is not dropped to avoid damaging the cabinet which defeats safeguards or injuring yourself. If the unit has been dropped or the cabinet has been damaged, unplug it and have it checked by an authorised service technician to restore the safeguards.
11. The fact that the unit operates satisfactorily does not imply that the unit is properly earthed or that it is completely safe. If in any doubt about the effective earthing of the unit, contact a qualified electrician.
12. **Never block the ventilation holes** or place the unit in any enclosure unless proper ventilation is provided.
13. Never place the unit near or over a radiator, heat register or stove.
14. Avoid locations where the instrument is exposed directly to the sunlight to prevent overheating.

### *Use of the instrument*


- This instrument is especially designed to create the electrical field for the electrophoresis process. Do not use it for other purposes than for which it is intended! The behaviour of the device might not be suitable for it.
- The instrument will show the following information during the start-up:
  - Model
  - Version number
  - Serial number

This information is required when contacting our service department ([support@consort.be](mailto:support@consort.be)).

- Make use of the instrument while studying this manual.

## *Procedure?*

This is the main menu that is shown after powering up the device and has 3 options:

### 1. MANUAL

To start a regular procedure with fixed parameter limits for Voltage, Current, Power and Timer. See [Manual Programming](#).

### 2. METHOD


To start a procedure with multiple steps to vary the parameter limits automatically while running. This method also allows Voltage Gradient settings. See [Method Programming](#).

### 3. SETUP

To modify the settings of the [system](#), [detections](#) and [data-logging](#).

## *Setup*

Use this menu to set the system settings of the power supply for language, display contrast, real time clock, password, detections, data recording.


## *System*

### *Language*

The power supply has several languages built-in for the communication with the user.

1. When in **SETUP**, follow these steps to select **LANGUAGE** and press ✓:


2. Select the desired language and press ✓.

## Contrast

According to the location where the device is installed, it might be required to adjust the display contrast.

1. When in **SETUP**, follow these steps to select **CONTRAST** and press ✓:

2. Adjust the contrast setting and press ✓ to confirm.


## Date/Time

This device has a precise real time clock built in. The real time clock is used for reports and data recording. The date and the time can be adjusted.

1. When in **SETUP**, follow these steps to select **DATE/TIME** and press ✓:


2. Set consecutively the date, month, year, hour, minute and seconds. Confirm each with ✓ or stop entering the other values by pressing ▶/■.

## Password


When desired, a password can be set which will be requested when trying to modify the setup or the parameters of a procedure. Without a password will the user only be able to choose, run or stop an experiment.

1. When in **SETUP**, follow these steps to select **PASSWORD** and press ✓:

2. A question to activate the password will appear. Choose **YES** and confirm to enable a password, set to **NO** to disable the password.

- a. When enabled, the device will request twice the same combination of five keystrokes. This combination will be required to modify any setting.

An error will be given if the second password doesn't correspond. The password setting will be disabled, retry to enable.


## Detection

Choose **DETECTION** to alter the actions when a low current or a power failure is detected.

### Low Current

When a very low current is detected as soon as the Voltage limit has been reached, the power will shut down the procedure. An [audible error message](#) will be given.

1. When in **SETUP**, follow these steps to select **LOW CURRENT** and press ✓:


2. Decide whether or not the output Voltage should shut down at a too low output current and press ✓.

### Power fail

When the detection of a power failure is enabled, the power supply will continue to proceed the procedure if the timer was still running at the moment that the power failure occurred. After powering it again, the power supply will automatically proceed with the experiment for the remaining time.

1. When in **SETUP**, follow these steps to select **POWER FAIL** and press ✓.

2. Decide whether or not the instrument should proceed with the experiment for the remaining time after a power failure and press ✓.

## Data


Starting from version 1.5 is this **DATA** menu available. Check our support website ([www.consort.be/wiki](http://www.consort.be/wiki)) for firmware updates of the instrument.

The instrument logs all experiments into its internal memory and will hold the last 250 experiments or less if the memory size would be exceeded. Older data will be lost. The available data can be read using the DIS3 [USB-version] software which is freely distributed on our website.

### Info

The menu in the instrument allows to show information about the memory status of the logged data.

1. When in **SETUP**, follow these steps to see the status of the data log memory:


2. When the info appears, press ✓ to return.

## Erase

The menu in the instrument allows to erase the data log memory.

1. When in **SETUP**, follow these steps when required to erase the data log memory:

2. A question to erase the data will appear on the display. Choose **YES** and press **✓** to confirm that the data should be erased.

## Info

Choose **INFO** to find specific information about the device, the parameter limits, the firmware version and more:

1. When in **SETUP**, follow these steps when required to erase the data log memory:

2. The information is shown over several lines. Use the **↑↓** - keys to scroll. Use **✓** or **▶/■** to return to the previous menu.



## Method programming

Method programming permits to store up to 9 different methods, each with 9 steps. Each step is able to recall a next one, providing a flexible multiple step function for special techniques.


The method mode also permits to program a linear voltage gradient for any step provided the limiting current or power is not attained.

Parameters can be changed temporarily without interrupting the run.

1. Verify if the instrument is switched off.
2. Connect the electrophoresis tank to the instrument.
3. Switch the instrument on.  
If the timer is still running after switching on, the power supply will automatically proceed with a procedure for the remaining time. A temporary menu is shown: Choose YES within 15 sec. to stop the procedure. Choose NO or wait 15 sec. to continue the procedure for the remaining timer value.
4. Select **METHOD** and press ✓:



5. Select the desired method and press ✓:

6. Choose YES to modify the name of the procedure, NO to proceed with step 8.

7. To modify the name, select each character from left to right, max. 12 characters.  
Select the space character to end a name with less than 12 characters and return to point 6.


Selectable characters are:  
'a' to 'z', 'A' to 'Z', '0' to '9', ' ', ',', '.', ','

8. When previously multiple steps have been programmed, select the desired step and press ✓:

9. Decide whether or not a linear voltage gradient has to be performed and press ✓:

10. Enter the desired Voltage limit and press ✓ :

METHOD 1 1
VOLTAGE
300 V

When performing a Voltage gradient, enter the desired end Voltage. The start Voltage is the end point of the previous step. This is 0 V when programming the first step as gradient step.

11. Enter the desired current limit and press ✓ :

METHOD 1 1
CURRENT
500 mA

This limit is automatically set to its maximum and cannot be changed in case of a Voltage gradient step.

12. Enter the desired power limit and press ✓ :

METHOD 1 1
POWER
150 W

This limit is automatically set to its maximum and cannot be changed in case of a Voltage gradient step.

13. Enter the desired timer unit (h or kWh) and press ✓ :

METHOD 1 1
TIMER UNIT
h

h

The timer will count in hours and minutes.

kWh

The timer will run in kiloVolthours.

14. Enter the desired timer value and press ✓ :

METHOD 1 1
TIMER
00:30 h

The value format is hh:mm when entering a value in hours (h). The kiloVolthours (kWh) are entered as a value with two decimal places.

15. Decide if the experiment should automatically proceed with a next step, press ✓ :

METHOD 1 1
END METHOD
NO

Choose NO to resume at point 8 to proceed with the next step, choose YES to stop after this step.



16. Press ▶/■ to start the experiment.

- At any time, you can return to the main menu by holding ✓ pressed.
- At any time, you can avoid programming unnecessary parameters. Simply press ▶/■ instead of ✓ after entering the last relevant value and the instrument will immediately start the experiment. The other parameters are kept to their previous programmed values.
- To ensure that the desired parameter is kept constant, we recommend to preset all other parameters to their maximum.
- The behaviour of the power supply during a procedure is described in the chapter [Running Condition](#).



## Manual programming

Manual programming permits to set voltage, current, power and time limits for a simple routine electrophoresis run. Parameters can be changed temporarily without interrupting the run.



1. Verify if the instrument is switched off.
2. Connect the electrophoresis tank to the instrument.
3. Switch the instrument on.
4. Select MANUAL and press ✓ :
5. Enter the desired voltage limit and press ✓ :

6. Enter the desired current limit and press ✓ :
7. Enter the desired power limit and press ✓ :

8. Enter the desired timer unit (h or kWh) and press ✓ :
9. Enter the desired timer value and press ✓ :

**h**  
The timer will count in hours and minutes.

**kWh**  
The timer will count in kilovolthours.

**h**  
The value format is **hh:mm**.

**kWh**  
The format has a value with 2 decimal places..


Set to 00:00 when required to continue until **▶/■** is pressed or the device is switched off.

10. Press **▶/■** to start the experiment.

- At any time, you can return to the main menu by holding ✓ pressed.
- At any time, you can avoid programming unnecessary parameters. To do so, simply press **▶/■** instead of ✓ after entering the last relevant value and the instrument will immediately start the experiment keeping the other parameters to their previous values.
- To ensure that the desired parameter is kept constant, we suggest to preset all other parameters to its maximum.
- The behaviour of the power supply during a procedure is described in the chapter [Running Condition](#).

## Running condition

### Display

When running an experiment, the display allows to see the actual parameters. It is divided into three fields:

- The **upper left** field shows:
  - a flashing **⚡** symbol to indicate there is a possible dangerous condition on the output sockets. **Avoid to contact the outlet sockets and cables!**
  - the highlighted measured parameter. Use the **↑↓** keys to choose another parameter in this field.
- The **upper right** field shows the other 3 parameters.
- The **bottom line** shows date and time at the right and the actual combination of method and step at the left when such a procedure is running, no indication when running a manual procedure. The name of a method is given between the upper left field and the bottom line.
- The **bottom line** is also used for system or error messages when required.

More details to be read:

- The unit of the parameter that is being kept **constant** is given with **flashing bold** characters.
- The timer field counts down to show the remaining time.
- The timer will count upwards when no timer value has been pre-set. The procedure will run continuously until **▶/■** is begin pressed or the device is switched off.
- When performing a method with multiple steps will the display show 2 different symbols in the combination of the method and the step:
  - **/** A step with Voltage gradient.
  - **|** A step with a pre-set fixed maximum Voltage.

### Buttons

The keys can be used as follows:


- **↑↓** = Buttons for selecting another parameter in the upper left corner of the display.
- **✓** = Button to modify the parameters.
- **▶/■** = Button to stop or to pause the procedure.

### Adjusting parameters

#### During the run

It is possible to modify the pre-set maximum parameters while running. These modifications are temporarily and are not stored in the non-volatile memory. It is not possible to modify the settings of a Voltage gradient step.

- Press **✓** to enter the mode to adjust:

- The display shows the values set for V, mA, and W. The timer remains counting down when a timer value has been set.
- Choose the limit to adjust with the **✓** button.
- Adjust the value with the **↑↓** buttons.

- The display will return automatically to the normal run mode when no button is pressed for 5 seconds.

### **When paused**

A procedure can be paused when a timer has been set. This allows also to check the migration during the electrophoresis process.

- Press ▶/■ to **pause** the running procedure.
- The device will show the descending Voltage at the output sockets until it is sufficiently lowered to a safe value before proceeding.
- Press ✓ to continue. To **stop** the procedure completely, hold ✓ pressed until the main menu appears.
- In case of a Method procedure is it possible to modify the settings of the next steps in the same manner as during the programming of the method. The modifications will be stored in the non-volatile memory.
- Press ▶/■ to continue with the procedure.

### **End of a procedure**

- The procedure will stop automatically when the timer has reached its end value. A signal will sound for 10 seconds. The sound signal can be ended by pressing a button.
- The procedure can also be stopped by pressing ▶/■. The power supply will then be in a paused condition. To stop the procedure completely, hold ✓ pressed until the main menu appears.
- The power supply will remain showing the display of the running condition until the voltage has lowered to a safety value.
- When the procedure has finished, a report will be shown on the display. Scroll the display lines with the ↑↓ buttons to see it completely. Press ✓ or ▶/■ to return to the main menu.

## *Alarm messages*

In case a problem or a hazardous situation is detected, the device will interrupt the experiment with one of the following errors.

A sound will also be produced that can be stopped by pressing any key.

The hazardous situations require that the device needs to be switched off before starting the experiment again.

See also the [frequently asked questions](#) at the end of this manual.

### **GROUND LEAKAGE**

A dangerous ground leakage is present.

- Your power supply is **NOT** defective!
- This is a safety feature to protect the user.
- Check thoroughly the complete assembly: device, tank, connections, cables, ... .

### **SHORT-CIRCUIT**

The connected tank is short-circuited or has a very low resistance.

- Your power supply is **NOT** defective!
- This feature prevents damage to the device or tank.
- Check thoroughly the complete assembly: device, tank, connections, cables, buffer, ... .
- When the fuses are blown, check the website on how to replace these.

### **OVERLOAD**

The connected tank has a too low resistance.

- Your power supply is **NOT** defective!
- Check thoroughly the complete assembly: device, tank, connections, cables, buffer, ... .
- Check also the best choice of power supply for the required application. e.g. It is not possible to perform Electroblotting requiring high currents with a high Voltage device that is designed for IEF.

### **LOW CURRENT**

The output current is very low or zero. Your power supply is **NOT** defective!

This can be caused by:

- an interruption of the external cables or tank,
- the cable connectors are not firmly or correctly placed into the sockets,
- erroneous procedure.

It is possible to disable this alarm in the [SETUP](#) menu. See [Frequently Asked Questions](#) at the end of this manual.

### **TECHNICAL PROBLEM**

The device has detected a possible failure in the electronics. Contact your supplier.

## *Troubleshooting*

See our website: [Consort Power Supply Troubleshooting](#).

## **RESET**

A general RESET of the device might help to solve some errors.

- Switch the device on while holding **✓** pressed.
- Confirm the question on the display with **YES** to restore the original factory settings of the device.
- **Attention!** All procedure parameters, data and settings will be reset.

## USB Digital Port

This power supply has a digital USB port allowing to interface with a computer. Programs can read and control the power supply completely. More information can be found on our support website: [Consort Support](#).


## Driver

The device has a USB port which requires that a driver is installed prior to connecting the USB cable between the computer and the device.

The USB driver can be downloaded using the link for the USB drivers on our website: [Consort Software](#).

## DIS-3

A freely distributed program (for Windows) to read the logged information of the power supply can be downloaded from our website: [Consort Software Support](#).

## Frequently Asked Questions

### What are the relations between Voltage, Current, Power and Resistance?

Power (W) = Voltage (V) x Current (A)

Resistance (Ω) = Voltage (V) / Current (A)

### How does a power supply react after pressing RUN?

The internal generator will start building up the high voltage at the output terminals while voltage and current are constantly measured and power is calculated. When one of the preset parameters is exceeded, the generator stops and will keep that parameter constant.

### How important is the resistance of an electrophoresis unit?

The resistance of an electrophoresis unit depends on its size, gel thickness, amount of buffer, buffer conductivity and temperature. This resistance will normally decrease in time due to a slowly increasing temperature. Electrophoresis units which have a resistance below the minimum load resistance of a power supply will trigger an alarm! Read the output voltage and current during a run to measure the resistance and use above formula to calculate the value.

### How to keep a constant voltage during a run?

Program the desired voltage and a higher current and power than the maximum expected values:

Voltage > Current x Resistance

Power > Voltage x Current

### How to keep a constant current during a run?

Program the desired current and a higher voltage and power than the maximum expected values:

Voltage > Current x Resistance

Power > Voltage x Current

### How to keep a constant power during a run?

Program the desired power and a higher voltage and current than the maximum expected values:

Voltage > Current x Resistance

Current > Voltage / Resistance

### Why are my output values different from those of a similar experiment?

Either your programmed parameters are not equal to those described or the resistance of your electrophoresis unit is different (see above). It cannot be due to e.g. another model of power supply as the relations between Voltage, Current, Power and Resistance are monitored in the same way by any instrument (the electrical laws cannot be disregarded!).

### What about connecting more than one unit to the same power supply?

The outlets being in parallel each electrophoresis unit will be supplied with exactly the same voltage. However, current and power may differ due to differences between them even when the same model, gel, buffers, etc... are used. Therefore, it is recommended to run several electrophoresis units only in the constant voltage mode on the same power supply.

Mind also that adding tanks will increase the current and the power. These parameters may be limited by the power supply and cause it to switch over to constant current or constant power.

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EU

### DECLARATION OF CONFORMITY

We declare under our sole responsibility that the products

**Electrophoresis power supplies**  
content of the type numbers

**EV0220, EV1450, EV2230,  
EV2310, EV2320, EV2650, EV3020,  
EV3150, EV3330, EV3610, EV3620**

to which this declaration relates is in conformity  
with the following standards of the European Directives

**LOW VOLTAGE DIRECTIVE 2014/35/EU**

**EN61010-1: 2010**

**EMC DIRECTIVE 2014/30/EU**

**EN61326-1: 2013**

**ROHS DIRECTIVE 2011/65/EU**

**EN50581: 2012**

**Name** Jan De Ceuster  
**Title** Manager  
**Date** 27 May 2016

**Signature**

A handwritten signature in blue ink, appearing to read 'Jan De Ceuster'.

## *Warranty*

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory.

Accessories and breakable items such as electrodes are not warranted unless proven to be defective before shipment.

The original purchase order numbers, Consort invoice numbers and serial numbers of the products must be provided.

**Consort** will repair all defective equipment returned to it during the warranty period without charge (CIF Turnhout in Belgium prepaid by sender), provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures described in the instruction manual and when not having been subject to accident, alteration, misuse or abuse.

If the products were used or came into contact with fluids, a MSDS (material safety data sheet) must be supplied prior to issuing a return authorisation.

A return authorisation must be obtained from **Consort** before returning any product for warranty repair on a freight prepaid basis!

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## EVHV\_ADAPTER

The power supplies with a maximum Voltage of 1500V or more, have special safety output sockets for High Voltage. These are deep recessed 4 mm banana sockets for which many plugs of tanks will be suitable. Because the output sockets have a maximum diameter that can be accepted, it will be required in some cases to make use of a pair of adapters. 1 such pair is supplied with this type of power supply.

Here are the instructions how to install these.

### Contents

- 1 red colored screw-clamping adapter
- 1 black colored screw-clamping adapter
- 1 Hex key tool 1.5 mm

### Installation

**For protection against electric shock, the parts must not be installed or removed into or from an active power supply. Make sure the power supply is switched off and the power cable is unplugged to prevent accidental powering!**

When required, install these 4mm adapters in a pair of the output sockets as follows:

- When it was required to switch off or to remove the power cable, wait at least 10 seconds afterwards before continuing with the next steps!
- Insert the similar colored adapter into the socket.
- Insert the 1.5 mm Hex key tool in the middle of the adapter and turn clockwise to expand the plug clamps until the adapter is fixed in place.
- To avoid a rotation, hold the socket at the same time.

## **Consort** bvba

Hertenstraat 56 / 9  
2300 Turnhout  
Belgium

**Tel** : (+32) [0]14 41 12 79  
**Fax** : (+32) [0]14 42 91 79

**Sales** : sales@consort.be  
**Support** : support@consort.be  
**Information** : info@consort.be