



These stand-alone, complete units for electro-blotting, mini, mini-wide and maxi formats, include vertical tank and lid, electro-blotting insert, cassettes and fibre pads.

Units are interchangeable with the vertical slab and tube gel inserts.

Hinged rigid cassettes allow rapid set up and convenience and ensure even contact between the gel and membrane is maintained.

Recommended power supply: EV3020

● Ordering codes

Code	Description
EVS1100-BLOT	Mini blotter: tank & lid + 4 cassettes 10x10 cm + 8 fibre pads + cooling pack
EVS1200-BLOT	Mini-wide Blotter: tank & lid + 3 cassettes 20x10 cm + 6 fibre pads + cooling pack
EVS1300-BLOT	Maxi blotter: tank & lid + 3 cassettes 20x20 cm + 6 fibre pads + cooling pack

Replacement parts & Accessories

Code	Description
EVS1100-BLOTINSERT	Mini blot insert + 4 cassettes 10x10 cm + 8 fibre pads
EVS1100-CASSETTE	Mini blot cassette, 10x10 cm
EVS1100-FIBREPAD	Mini fibre pads, 10x10 cm, pk/6
EVS1200-BLOTINSERT	Mini-wide blot insert + 3 cassettes 20x10 cm + 6 fibre pads
EVS1200-CASSETTE	Mini-wide blot cassette, 20x10 cm
EVS1200-FIBREPAD	Mini-wide fibre pads, 20x10 cm, pk/6
EVS1200-BLOTINSERT	Maxi blot insert + 3 cassettes 20x20 cm + 6 fibre pads
EVS1300-CASSETTE	Maxi blot cassette, 20x20 cm
EVS1300-FIBREPAD	Maxi fibre pads, 20x20 cm, pk/6