

- Quick & Efficient
- Convenient & Versatile
- Safe & Reliable

Quick transfer times make these systems ideal for the rapid and efficient transfer of nucleic acids and proteins from agarose or acrylamide gels.

Solid plate style electrodes assure even pressure and complete molecular transfer.

Galileo semi-dry electroblotters will accommodate most pre-cast gels, and for many larger format pre-cast gels, once the stacking area of the gel is removed the "working" area will also fit easily into the devices. Since only the transfer "sandwich" of gel, membrane and blotting papers must be kept wet, much less buffer is required than traditional tank transfer systems.

Interlocking safety lid prevents reverse orientation of electric filed. High quality plate electrodes (stainless steel cathode and platinum anode) and gold plated components assure long, trouble-free performance. System is supplied complete with power leads and a sample pack of our ultra-pure cotton fiber filter paper.



Code	Description
ESDB3100	Semi Dry Electroblotter, 11x11cm blotting area
ESDB3300	Semi Dry Electroblotter, 21x21cm blotting area
EV2310	Power supply, 300 V, 1000 mA, 150 W
EV3020	Power supply, 300 V, 2000 mA, 300 W

## Two sizes to fits your needs:

Unit	ESDB3100	ESDB3300
Transfer area (cm)	11X11	21X21
Buffer Volume (ml)	50	100
Unit Dimensions	19x19x6 cm	29x29x6 cm
Running condition	0.8-3mA/cm	0.8-3mA/cm
Recommended voltage	<100V	<100V
Warranty	36 months	36 months
Recommended PS	EV2650	EV3150



- Description

**Quick & Efficient**

Quick transfer times make these systems ideal for the rapid and efficient transfer of nucleic acids and proteins from agarose or acrylamide gels.

**Convenient & Versatile**

Solid plate style electrodes assure even pressure and complete molecular transfer.

Galileo semi-dry electroblotters will accommodate most pre-cast gels, and for many larger format pre-cast gels, once the stacking area of the gel is removed the "working" area will also fit easily into the devices. Since only the transfer "sandwich" of gel, membrane and blotting papers must be kept wet, much less buffer is required than traditional tank transfer systems.

**Safe & Reliable**

Interlocking safety lid prevents reverse orientation of electric field.

High quality plate electrodes (stainless steel cathode and platinum anode) and gold plated components assure long, trouble-free performance. System is supplied complete with power leads and a sample pack of our ultra-pure cotton fiber filter paper.

- Ordering codes

Code	Description
ESDB3100	Semi Dry Electroblotter, 11x11cm blotting area
ESDB3300	Semi Dry Electroblotter, 21x21cm blotting area

